


=> d que stat 113  
L11 STR



VAR G2=O/N  
NODE ATTRIBUTES:  
DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED


GRAPH ATTRIBUTES:  
RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 13

STEREO ATTRIBUTES: NONE  
L13 832 SEA FILE=REGISTRY SSS FUL L11

100.0% PROCESSED 56776 ITERATIONS  
SEARCH TIME: 00.00.03

832 ANSWERS


=> d que stat 118  
L11 STR



VAR G2=O/N  
NODE ATTRIBUTES:  
DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 13

STEREO ATTRIBUTES: NONE  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR



VAR G1=CH/14  
 VAR G2=O/NH  
 VAR G3=NH3/16/19/22  
 NODE ATTRIBUTES:  
 CONNECT IS E1 RC AT 10  
 DEFAULT MLEVEL IS ATOM  
 DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
 RING(S) ARE ISOLATED OR EMBEDDED  
 NUMBER OF NODES IS 24

STEREO ATTRIBUTES: NONE  
 L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16

100.0% PROCESSED 832 ITERATIONS  
 SEARCH TIME: 00.00.01

553 ANSWERS

=> d l64 1-17  
 L64 ANALYZE L18 1- LC : 17 TERMS

TERM #	# OCC	# DOC	% DOC	LC
1	535	535	96.75	CA
2	535	535	96.75	CAPLUS
3	255	255	46.11	TOXCENTER
4	115	115	20.80	USPATFULL
5	29	29	5.24	USPAT2
6	3	3	0.54	BIOSIS
7	3	3	0.54	CHEMLIST
8	3	3	0.54	MEDLINE
9	2	2	0.36	CASREACT
10	2	2	0.36	DIOGENES
11	2	2	0.36	IPA
12	2	2	0.36	TSCA
13	1	1	0.18	BIOBUSINESS
14	1	1	0.18	CANCERLIT
15	1	1	0.18	CHEMINFORMRX
16	1	1	0.18	PIRA
17	1	1	0.18	PROMT

\*\*\*\*\* END OF L64\*\*\*

=&gt; d que nos 133

```
L11          STR
L13          832 SEA FILE=REGISTRY SSS FUL L11
L16          STR
L18          553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16
L24          715 SEA FILE=HCAPLUS ABB=ON  PLU=ON  L18
L25          56930 SEA FILE=HCAPLUS ABB=ON  PLU=ON  IMMUNOASSAY+PFT,NT/CT
L26          47927 SEA FILE=HCAPLUS ABB=ON  PLU=ON  "IMMUNOCHEMICAL ANALYSIS (L)
              IMMUNOASSAY"+PFT,NT/CT
L27          5294 SEA FILE=HCAPLUS ABB=ON  PLU=ON  AGGLUTINATION+PFT,NT/CT
L28          17 SEA FILE=HCAPLUS ABB=ON  PLU=ON  L24 AND (L25 OR L26 OR L27)
L29          QUE ABB=ON  PLU=ON  ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E
              LISA OR RIA
L31          16 SEA FILE=HCAPLUS ABB=ON  PLU=ON  L24 (L) L29
L33          23 SEA FILE=HCAPLUS ABB=ON  PLU=ON  L28 OR L31
```

=&gt; d his 139

(FILE 'USPATFULL, USPAT2' ENTERED AT 14:41:23 ON 20 SEP 2005)

L39 3 S L35 AND L38

=&gt; d que nos 139

```
L11          STR
L13          832 SEA FILE=REGISTRY SSS FUL L11
L16          STR
L18          553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16
L35          78 SEA L18
L38          23444 SEA G01N033-5?/IPC
L39          3 SEA L35 AND L38
```

=&gt; d que nos 162

```
L40          2142 SEA FILE=WPIX ABB=ON  PLU=ON  (B415 (P) B701 (P) B713 (P) B815
              (P) B831 (P) H1 (P) J011)/M0,M1,M2,M3,M4,M5,M6
L50          408 SEA FILE=WPIX ABB=ON  PLU=ON  C08F030-02/IPC
L59          58 SEA FILE=WPIX ABB=ON  PLU=ON  L40 AND (G01N033-53?/ICM,ICS)
L62          2 SEA FILE=WPIX ABB=ON  PLU=ON  L59 AND L50
```

=&gt; d que nos 166

```
L11          STR
L13          832 SEA FILE=REGISTRY SSS FUL L11
L16          STR
L18          553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16
L63          QUE ABB=ON  PLU=ON  ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E
              LISA OR RIA OR ?COAG?
L65          260 SEA FILE=TOXCENTER ABB=ON  PLU=ON  L18
L66          31 SEA FILE=TOXCENTER ABB=ON  PLU=ON  L65 AND L63
```

=&gt; d que nos 172

```
L11          STR
L13          832 SEA FILE=REGISTRY SSS FUL L11
L16          STR
L18          553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16
L63          QUE ABB=ON  PLU=ON  ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E
```

LISA OR RIA OR ?COAG?

L67	55	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L18
L68	282436	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	IMMUNOASSAY+PFT,NT/CT
L69	6576	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	AGGLUTINATION+PFT,NT/CT
L70	3	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L67 AND (L68 OR L69)
L71	21	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L67 AND L63
L72	21	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	(L70 OR L71)

=> d que nos 173

L11		STR			
L13	832	SEA	FILE=REGISTRY	SSS	FUL L11
L16		STR			
L18	553	SEA	FILE=REGISTRY	SUB=L13	SSS FUL L16
L73	0	SEA	FILE=EMBASE	ABB=ON	PLU=ON L18

=> d his 175

(FILE 'BIOSIS, CANCERLIT' ENTERED AT 15:20:16 ON 20 SEP 2005)

L75	13	S	L74 AND L63
-----	----	---	-------------

=> d que nos 175

L11		STR			
L13	832	SEA	FILE=REGISTRY	SSS	FUL L11
L16		STR			
L18	553	SEA	FILE=REGISTRY	SUB=L13	SSS FUL L16
L63		QUE	ABB=ON	PLU=ON	?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E
			LISA OR RIA OR ?COAG?		
L74	57	SEA	L18		
L75	13	SEA	L74 AND L63		

=> d his 181

(FILE 'HCAPLUS, TOXCENTER, WPIX, MEDLINE, BIOSIS, CANCERLIT, EMBASE, PASCAL, JICST-EPLUS, DRUGU, BIOTECHNO, BIOTECHDS, SCISEARCH, CONF, CONFSCI, DISSABS' ENTERED AT 15:25:01 ON 20 SEP 2005)

L81	6	DUP	REM L80 (6 DUPLICATES REMOVED)
-----	---	-----	--------------------------------

=> d que nos 181

L63		QUE	ABB=ON	PLU=ON	?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E
			LISA OR RIA OR ?COAG?		
L76	1325	SEA	SUMIDA, K?/AU		
L77	17872	SEA	WADA, K?/AU		
L78	14175	SEA	ISHIHARA, K?/AU		
L79	3656	SEA	(L76 OR L77 OR L78) AND L63		
L80	12	SEA	L79 AND WAKO/CS,SO,PA		
L81	6	DUP	REM L80 (6 DUPLICATES REMOVED)		

=> d his ful

(FILE 'HOME' ENTERED AT 12:50:28 ON 20 SEP 2005)

FILE 'STNGUIDE' ENTERED AT 12:50:39 ON 20 SEP 2005

FILE 'HCAPLUS' ENTERED AT 12:51:02 ON 20 SEP 2005

L1 0 SEA ABB=ON PLU=ON US2003-626502/APPS  
FILE 'STNGUIDE' ENTERED AT 12:51:19 ON 20 SEP 2005  
FILE 'ZCAPLUS' ENTERED AT 12:52:40 ON 20 SEP 2005  
E JP2001-169051/APPS  
FILE 'HCAPLUS' ENTERED AT 12:52:48 ON 20 SEP 2005  
L2 1 SEA ABB=ON PLU=ON JP2001-169051/APPS  
SAVE TEMP L2 FET502HCAAPP/A  
D IBIB ED AB IND  
FILE 'STNGUIDE' ENTERED AT 12:53:16 ON 20 SEP 2005  
FILE 'HCAPLUS' ENTERED AT 12:54:13 ON 20 SEP 2005  
L3 1 SEA ABB=ON PLU=ON US2004-626502/APPS  
SAVE TEMP L3 FET502HCAAP2/A  
FILE 'STNGUIDE' ENTERED AT 12:54:48 ON 20 SEP 2005  
FILE 'HCAPLUS' ENTERED AT 12:55:14 ON 20 SEP 2005  
D IBIB ED AB IND  
FILE 'STNGUIDE' ENTERED AT 12:55:14 ON 20 SEP 2005  
FILE 'HCAPLUS' ENTERED AT 12:56:22 ON 20 SEP 2005  
L4 0 SEA ABB=ON PLU=ON L2 NOT L3  
FILE 'STNGUIDE' ENTERED AT 12:56:30 ON 20 SEP 2005  
FILE 'WPIX' ENTERED AT 12:56:35 ON 20 SEP 2005  
L5 1 SEA ABB=ON PLU=ON US2004-626502/APPS  
SAVE TEMP L5 FET502WPIAPP/A  
D IALL CMC  
FILE 'STNGUIDE' ENTERED AT 12:57:05 ON 20 SEP 2005  
FILE 'REGISTRY' ENTERED AT 12:57:35 ON 20 SEP 2005  
FILE 'HCAPLUS' ENTERED AT 12:57:42 ON 20 SEP 2005  
L6 TRA L3 1- RN : 9 TERMS  
FILE 'REGISTRY' ENTERED AT 12:57:45 ON 20 SEP 2005  
L7 9 SEA ABB=ON PLU=ON L6  
SAVE TEMP L7 FET502REGAPP/A  
D SCAN  
FILE 'STNGUIDE' ENTERED AT 12:58:14 ON 20 SEP 2005  
D SAVED  
FILE 'LREGISTRY' ENTERED AT 13:02:29 ON 20 SEP 2005  
L8 STRUCTURE UPLOADED  
L9 STR L8  
FILE 'REGISTRY' ENTERED AT 13:05:36 ON 20 SEP 2005  
L10 0 SEA SSS SAM L9  
D QUE STAT  
FILE 'STNGUIDE' ENTERED AT 13:05:56 ON 20 SEP 2005

L11 FILE 'LREGISTRY' ENTERED AT 13:06:35 ON 20 SEP 2005  
STR L9

L12 FILE 'REGISTRY' ENTERED AT 13:07:51 ON 20 SEP 2005  
25 SEA SSS SAM L11

FILE 'STNGUIDE' ENTERED AT 13:08:58 ON 20 SEP 2005  
D QUE STAT

L13 FILE 'REGISTRY' ENTERED AT 13:09:57 ON 20 SEP 2005  
832 SEA SSS FUL L11  
SAVE TEMP L13 FET502PSET1/A

L14 5 SEA ABB=ON PLU=ON L13 AND L7

L15 4 SEA ABB=ON PLU=ON L7 NOT L14  
D SCAN

FILE 'STNGUIDE' ENTERED AT 13:10:56 ON 20 SEP 2005  
D SAVED

L16 FILE 'LREGISTRY' ENTERED AT 13:11:22 ON 20 SEP 2005  
STR L11

L17 FILE 'REGISTRY' ENTERED AT 13:15:51 ON 20 SEP 2005  
26 SEA SUB=L13 SSS SAM L16

L18 553 SEA SUB=L13 SSS FUL L16  
SAVE TEMP L18 FET502RSET1/A

FILE 'STNGUIDE' ENTERED AT 13:17:43 ON 20 SEP 2005  
D SAVED

L19 FILE 'REGISTRY' ENTERED AT 13:18:01 ON 20 SEP 2005  
5 SEA ABB=ON PLU=ON L14 AND L18

FILE 'STNGUIDE' ENTERED AT 13:18:30 ON 20 SEP 2005

FILE 'REGISTRY' ENTERED AT 13:18:37 ON 20 SEP 2005  
D SCAN

FILE 'STNGUIDE' ENTERED AT 13:18:45 ON 20 SEP 2005

L20 FILE 'HCAPLUS' ENTERED AT 13:19:12 ON 20 SEP 2005  
715 SEA ABB=ON PLU=ON L18

FILE 'STNGUIDE' ENTERED AT 13:19:23 ON 20 SEP 2005

L21 FILE 'LREGISTRY' ENTERED AT 13:19:43 ON 20 SEP 2005  
STR L16

FILE 'STNGUIDE' ENTERED AT 13:21:43 ON 20 SEP 2005

L22 FILE 'REGISTRY' ENTERED AT 13:22:37 ON 20 SEP 2005  
495 SEA ABB=ON PLU=ON L18 AND PMS/CI

L23 FILE 'HCAPLUS' ENTERED AT 13:23:06 ON 20 SEP 2005  
595 SEA ABB=ON PLU=ON L22

FILE 'STNGUIDE' ENTERED AT 13:23:11 ON 20 SEP 2005

FILE 'ZCAPLUS' ENTERED AT 14:31:19 ON 20 SEP 2005  
E IMMUNOASSAY/CT

E E15+ALL  
E AGGLUTINATION/CT  
E E102+ALL

FILE 'HCAPLUS' ENTERED AT 14:32:44 ON 20 SEP 2005  
L24 715 SEA ABB=ON PLU=ON L18  
L25 56930 SEA ABB=ON PLU=ON IMMUNOASSAY+PFT,NT/CT  
L26 47927 SEA ABB=ON PLU=ON "IMMUNOCHEMICAL ANALYSIS (L) IMMUNOASSAY"+P  
FT,NT/CT  
L27 5294 SEA ABB=ON PLU=ON AGGLUTINATION+PFT,NT/CT  
L28 17 SEA ABB=ON PLU=ON L24 AND (L25 OR L26 OR L27)

FILE 'STNGUIDE' ENTERED AT 14:33:30 ON 20 SEP 2005

FILE 'HCAPLUS' ENTERED AT 14:34:15 ON 20 SEP 2005  
D SCAN

FILE 'STNGUIDE' ENTERED AT 14:34:37 ON 20 SEP 2005

FILE 'HCAPLUS' ENTERED AT 14:37:13 ON 20 SEP 2005  
L29 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR ELISA  
OR RIA  
L30 62 SEA ABB=ON PLU=ON L24 AND L29  
L31 16 SEA ABB=ON PLU=ON L24 (L) L29  
L32 6 SEA ABB=ON PLU=ON L31 NOT L28  
D SCAN  
L33 23 SEA ABB=ON PLU=ON L28 OR L31

FILE 'STNGUIDE' ENTERED AT 14:39:15 ON 20 SEP 2005

FILE 'HCAPLUS' ENTERED AT 14:39:27 ON 20 SEP 2005  
SAVE TEMP L33 FET502HCA1/A

FILE 'STNGUIDE' ENTERED AT 14:39:44 ON 20 SEP 2005  
D QUE STAT  
D SAVED

FILE 'HCAPLUS' ENTERED AT 14:40:58 ON 20 SEP 2005  
L34 1 SEA ABB=ON PLU=ON L3 AND L33

FILE 'STNGUIDE' ENTERED AT 14:41:07 ON 20 SEP 2005

FILE 'USPATFULL' ENTERED AT 14:41:16 ON 20 SEP 2005

FILE 'USPATFULL, USPAT2' ENTERED AT 14:41:23 ON 20 SEP 2005  
L35 78 SEA ABB=ON PLU=ON L18  
L36 13193 SEA ABB=ON PLU=ON G01N033-53?/IPC  
L37 1 SEA ABB=ON PLU=ON L35 AND L36  
L38 23444 SEA ABB=ON PLU=ON G01N033-5?/IPC  
L39 3 SEA ABB=ON PLU=ON L35 AND L38

FILE 'STNGUIDE' ENTERED AT 14:42:35 ON 20 SEP 2005

FILE 'USPATFULL, USPAT2' ENTERED AT 14:43:29 ON 20 SEP 2005  
D SCAN L39  
SAVE TEMP L39 FET502USP1/A

FILE 'STNGUIDE' ENTERED AT 14:44:05 ON 20 SEP 2005  
D SAVED

FILE 'WPIX' ENTERED AT 14:53:50 ON 20 SEP 2005  
L40 2142 SEA ABB=ON PLU=ON (B415 (P) B701 (P) B713 (P) B815 (P) B831  
(P) H1 (P) J011)/M0,M1,M2,M3,M4,M5,M6  
L41 1 SEA ABB=ON PLU=ON L40 AND L5  
L42 63776 SEA ABB=ON PLU=ON G01N033-5?/IPC  
L43 17781 SEA ABB=ON PLU=ON (B11-C08E OR C11-C08E OR E11-C08E)/MC  
L44 31675 SEA ABB=ON PLU=ON S03-E14H/MC  
L45 127 SEA ABB=ON PLU=ON L40 AND (L42 OR L43 OR L44)  
L46 119 SEA ABB=ON PLU=ON L45 AND (L42 OR L44)  
L47 29713 SEA ABB=ON PLU=ON G01N033-53/IPC  
L48 87 SEA ABB=ON PLU=ON L46 AND (L47 OR L44)  
L49 62 SEA ABB=ON PLU=ON L48 AND L47  
L50 408 SEA ABB=ON PLU=ON C08F030-02/IPC  
L51 1 SEA ABB=ON PLU=ON L49 AND L50  
L52 36 SEA ABB=ON PLU=ON L49 AND (?ASSAY?/BIX OR ?IMMUNO?/BIX OR  
?AGGLUT?/BIX OR ELISA/BIX OR RIA/BIX)  
D TRI 1-3

FILE 'STNGUIDE' ENTERED AT 14:58:14 ON 20 SEP 2005

FILE 'WPIX' ENTERED AT 14:58:44 ON 20 SEP 2005

FILE 'STNGUIDE' ENTERED AT 14:59:18 ON 20 SEP 2005

FILE 'WPIX' ENTERED AT 15:00:55 ON 20 SEP 2005  
L53 34679 SEA ABB=ON PLU=ON (A04-A OR B04-C03B OR C04-C03B)/MC  
L54 2 SEA ABB=ON PLU=ON L52 AND (L50 OR L53)  
D TRI 1-2  
L55 0 SEA ABB=ON PLU=ON L52 AND L5  
L56 0 SEA ABB=ON PLU=ON L49 AND L5  
L57 32454 SEA ABB=ON PLU=ON G01N033-53?/IPC  
L58 67 SEA ABB=ON PLU=ON L40 AND L57  
L\*\*\* DEL 0 S K40 AND (G01N033-53?/ICM,ICS)  
L59 58 SEA ABB=ON PLU=ON L40 AND (G01N033-53?/ICM,ICS)  
D QUE  
L60 35 SEA ABB=ON PLU=ON L59 AND ((?ASSAY?/BIX OR ?IMMUNO?/BIX OR  
?AGGLUT?/BIX OR ELISA/BIX OR RIA/BIX) OR ?COAG?/BIX)  
L61 2 SEA ABB=ON PLU=ON L60 AND L50  
L62 2 SEA ABB=ON PLU=ON L59 AND L50  
D TRI 1-2  
SAVE TEMP L62 FET502WPI1/A

FILE 'STNGUIDE' ENTERED AT 15:07:37 ON 20 SEP 2005

D SAVED  
D QUE L29

FILE 'HCAPLUS' ENTERED AT 15:08:56 ON 20 SEP 2005  
L63 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR ELISA  
OR RIA OR ?COAG?  
SAVE TEMP L63 FET502QUE1/Q

FILE 'STNGUIDE' ENTERED AT 15:09:24 ON 20 SEP 2005

D SAVED

FILE 'REGISTRY' ENTERED AT 15:10:06 ON 20 SEP 2005  
L64 ANALYZE PLU=ON L18 1- LC : 17 TERMS  
D

FILE 'TOXCENTER' ENTERED AT 15:12:10 ON 20 SEP 2005  
L65 260 SEA ABB=ON PLU=ON L18



L66           31 SEA ABB=ON   PLU=ON   L65 AND L63  
              SAVE TEMP L66 FET502TOX1/A

FILE 'STNGUIDE' ENTERED AT 15:12:53 ON 20 SEP 2005  
              D SAVED

FILE 'MEDLINE' ENTERED AT 15:13:17 ON 20 SEP 2005  
              E IMMUNOASSAY/CT  
              E E121+ALL  
              E AGGLUTINATION/CT  
              E E167+ALL

L67           55 SEA ABB=ON   PLU=ON   L18  
L68       282436 SEA ABB=ON   PLU=ON   IMMUNOASSAY+PFT,NT/CT  
L69       6576 SEA ABB=ON   PLU=ON   AGGLUTINATION+PFT,NT/CT  
L70           3 SEA ABB=ON   PLU=ON   L67 AND (L68 OR L69)  
L71       21 SEA ABB=ON   PLU=ON   L67 AND L63  
L72       21 SEA ABB=ON   PLU=ON   (L70 OR L71)  
              D TRI 1-21

FILE 'STNGUIDE' ENTERED AT 15:14:56 ON 20 SEP 2005

FILE 'MEDLINE' ENTERED AT 15:16:05 ON 20 SEP 2005  
              SAVE TEMP L72 FET502MED1/A

FILE 'STNGUIDE' ENTERED AT 15:16:22 ON 20 SEP 2005  
              D SAVED

FILE 'EMBASE' ENTERED AT 15:16:43 ON 20 SEP 2005  
L73           0 SEA ABB=ON   PLU=ON   L18

FILE 'STNGUIDE' ENTERED AT 15:16:59 ON 20 SEP 2005

FILE 'REGISTRY' ENTERED AT 15:17:59 ON 20 SEP 2005  
              D ALL L64  
              D L64 1-17

FILE 'EMBASE' ENTERED AT 15:19:37 ON 20 SEP 2005  
              SAVE TEMP L73 FET502EMB1/A

FILE 'STNGUIDE' ENTERED AT 15:20:00 ON 20 SEP 2005

FILE 'BIOSIS, CANCERLIT' ENTERED AT 15:20:16 ON 20 SEP 2005  
L74           57 SEA ABB=ON   PLU=ON   L18  
L75       13 SEA ABB=ON   PLU=ON   L74 AND L63  
              SAVE TEMP L75 FET502MUL1/A  
              D SCAN

FILE 'STNGUIDE' ENTERED AT 15:21:38 ON 20 SEP 2005  
              D SAVED

FILE 'HCAPLUS, TOXCENTER, WPIX, MEDLINE, BIOSIS, CANCERLIT, EMBASE,  
PASCAL, JICST-EPLUS, DRUGU, BIOTECHNO, BIOTECHDS, SCISEARCH, CONF,  
CONFSCI, DISSABS' ENTERED AT 15:25:01 ON 20 SEP 2005  
L76       1325 SEA ABB=ON   PLU=ON   SUMIDA, K?/AU  
L77       17872 SEA ABB=ON   PLU=ON   WADA, K?/AU  
L78       14175 SEA ABB=ON   PLU=ON   ISHIHARA, K?/AU  
L79       3656 SEA ABB=ON   PLU=ON   (L76 OR L77 OR L78) AND L63  
L80       12 SEA ABB=ON   PLU=ON   L79 AND WAKO/CS,SO,PA  
L81       6 DUP REM L80 (6 DUPLICATES REMOVED)  
              ANSWERS '1-3' FROM FILE HCAPLUS

ANSWERS '4-5' FROM FILE BIOSIS  
ANSWER '6' FROM FILE PASCAL  
D SCAN  
SAVE TEMP L81 FET502MULINV/A  
D SAVED

FILE 'STNGUIDE' ENTERED AT 15:29:29 ON 20 SEP 2005  
FILE 'LREGISTRY' ENTERED AT 15:30:14 ON 20 SEP 2005  
FILE 'REGISTRY' ENTERED AT 15:30:17 ON 20 SEP 2005  
FILE 'ZCAPLUS' ENTERED AT 15:30:20 ON 20 SEP 2005  
FILE 'USPATFULL' ENTERED AT 15:30:25 ON 20 SEP 2005  
FILE 'USPAT2' ENTERED AT 15:30:29 ON 20 SEP 2005  
FILE 'HCAPLUS' ENTERED AT 15:30:32 ON 20 SEP 2005  
FILE 'TOXCENTER' ENTERED AT 15:30:36 ON 20 SEP 2005  
FILE 'WPIX' ENTERED AT 15:30:39 ON 20 SEP 2005  
FILE 'MEDLINE' ENTERED AT 15:30:45 ON 20 SEP 2005  
FILE 'BIOSIS' ENTERED AT 15:30:48 ON 20 SEP 2005  
FILE 'CANCERLIT' ENTERED AT 15:30:52 ON 20 SEP 2005  
FILE 'EMBASE' ENTERED AT 15:30:54 ON 20 SEP 2005  
FILE 'PASCAL' ENTERED AT 15:30:58 ON 20 SEP 2005  
FILE 'JICST-EPLUS' ENTERED AT 15:31:01 ON 20 SEP 2005  
FILE 'DRUGU' ENTERED AT 15:31:04 ON 20 SEP 2005  
FILE 'BIOTECHNO' ENTERED AT 15:31:09 ON 20 SEP 2005  
FILE 'BIOTECHDS' ENTERED AT 15:31:13 ON 20 SEP 2005  
FILE 'SCISEARCH' ENTERED AT 15:31:22 ON 20 SEP 2005  
FILE 'CONF' ENTERED AT 15:31:24 ON 20 SEP 2005  
FILE 'CONFSCI' ENTERED AT 15:31:29 ON 20 SEP 2005  
FILE 'DISSABS' ENTERED AT 15:31:36 ON 20 SEP 2005  
FILE 'STNGUIDE' ENTERED AT 15:31:38 ON 20 SEP 2005  
D QUE STAT L13  
D QUE STAT L18  
D L64 1-17  
D QUE NOS L33  
D QUE NOS L39  
D QUE NOS L62  
D QUE NOS L66  
D QUE NOS L72  
D QUE NOS L73

D QUE NOS L75

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS, CANCERLIT'  
ENTERED AT 15:33:39 ON 20 SEP 2005

L82           70 DUP REM L33 L39 L62 L66 L72 L73 L75 (23 DUPLICATES REMOVED)  
              ANSWERS '1-23' FROM FILE HCAPLUS  
              ANSWERS '24-26' FROM FILE USPATFULL  
              ANSWER '27' FROM FILE WPIX  
              ANSWERS '28-57' FROM FILE TOXCENTER  
              ANSWERS '58-67' FROM FILE MEDLINE  
              ANSWERS '68-70' FROM FILE BIOSIS

FILE 'STNGUIDE' ENTERED AT 15:34:17 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:34:35 ON 20 SEP 2005

D IBIB ED AB HITIND HITSTR RETABLE

FILE 'STNGUIDE' ENTERED AT 15:34:36 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:34:58 ON 20 SEP 2005

D IBIB ED AB HITIND HITSTR RETABLE 2-23

FILE 'STNGUIDE' ENTERED AT 15:35:09 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:35:41 ON 20 SEP 2005

D IBIB AB HITSTR 24-26

FILE 'STNGUIDE' ENTERED AT 15:35:44 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:36:20 ON 20 SEP 2005

D IALL ABEQ TECH ABEX 27

FILE 'STNGUIDE' ENTERED AT 15:36:23 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:36:55 ON 20 SEP 2005

D IBIB ED AB HITIND 28

FILE 'STNGUIDE' ENTERED AT 15:36:56 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:37:26 ON 20 SEP 2005

FILE 'STNGUIDE' ENTERED AT 15:37:30 ON 20 SEP 2005

FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' ENTERED AT  
15:37:46 ON 20 SEP 2005

D IBIB ED AB HITIND 29-

FILE 'STNGUIDE' ENTERED AT 15:37:53 ON 20 SEP 2005

D QUE L81

FILE 'HCAPLUS, BIOSIS, PASCAL' ENTERED AT 15:39:22 ON 20 SEP 2005

D IBIB ED AB L81 1-6

FILE 'STNGUIDE' ENTERED AT 15:39:24 ON 20 SEP 2005

FILE 'STNGUIDE' ENTERED AT 15:39:41 ON 20 SEP 2005

D QUE STAT L13  
D QUE STAT L18  
D QUE NOS L33  
D QUE NOS L39  
D QUE NOS L62  
D QUE NOS L66  
D QUE NOS L72  
D QUE NOS L73  
D QUE NOS L75  
D QUE NOS L81  
D L64 1-17

FILE HOME

FILE STNGUIDE

FILE CONTAINS CURRENT INFORMATION.

LAST RELOADED: Sep 16, 2005 (20050916/UP).

FILE HCAPLUS

Copyright of the articles to which records in this database refer is held by the publishers listed in the PUBLISHER (PB) field (available for records published or updated in Chemical Abstracts after December 26, 1996), unless otherwise indicated in the original publications. The CA Lexicon is the copyrighted intellectual property of the American Chemical Society and is provided to assist you in searching databases on STN. Any dissemination, distribution, copying, or storing of this information, without the prior written consent of CAS, is strictly prohibited.

FILE COVERS 1907 - 20 Sep 2005 VOL 143 ISS 13

FILE LAST UPDATED: 19 Sep 2005 (20050919/ED)

New CAS Information Use Policies, enter HELP USAGETERMS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

FILE ZCAPLUS

Copyright of the articles to which records in this database refer is held by the publishers listed in the PUBLISHER (PB) field (available for records published or updated in Chemical Abstracts after December 26, 1996), unless otherwise indicated in the original publications. The CA Lexicon is the copyrighted intellectual property of the American Chemical Society and is provided to assist you in searching databases on STN. Any dissemination, distribution, copying, or storing of this information, without the prior written consent of CAS is strictly prohibited.

FILE COVERS 1907 - 20 Sep 2005 VOL 143 ISS 13

FILE LAST UPDATED: 19 Sep 2005 (20050919/ED)

New CAS Information Use Policies, enter HELP USAGETERMS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

## FILE WPIX

FILE LAST UPDATED: 15 SEP 2005 <20050915/UP>  
MOST RECENT DERWENT UPDATE: 200559 <200559/DW>  
DERWENT WORLD PATENTS INDEX SUBSCRIBER FILE, COVERS 1963 TO DATE

>>> FOR A COPY OF THE DERWENT WORLD PATENTS INDEX STN USER GUIDE,  
PLEASE VISIT:  
[http://www.stn-international.de/training\\_center/patents/stn\\_guide.pdf](http://www.stn-international.de/training_center/patents/stn_guide.pdf) <<<

>>> FOR DETAILS OF THE PATENTS COVERED IN CURRENT UPDATES, SEE  
<http://thomsonderwent.com/coverage/latestupdates/> <<<

>>> FOR INFORMATION ON ALL DERWENT WORLD PATENTS INDEX USER  
GUIDES, PLEASE VISIT:  
<http://thomsonderwent.com/support/userguides/> <<<

>>> NEW! FAST-ALERTING ACCESS TO NEWLY-PUBLISHED PATENT  
DOCUMENTATION NOW AVAILABLE IN DERWENT WORLD PATENTS INDEX  
FIRST VIEW - FILE WPIFV.  
FOR FURTHER DETAILS: <http://www.thomsonderwent.com/dwpifv> <<<

>>> THE CPI AND EPI MANUAL CODES HAVE BEEN REVISED FROM UPDATE 200501.  
PLEASE CHECK:  
<http://thomsonderwent.com/support/dwpioref/reftools/classification/code-rev>.  
FOR DETAILS. <<<

## FILE REGISTRY

Property values tagged with IC are from the ZIC/VINITI data file  
provided by InfoChem.

STRUCTURE FILE UPDATES: 19 SEP 2005 HIGHEST RN 863478-08-4  
DICTIONARY FILE UPDATES: 19 SEP 2005 HIGHEST RN 863478-08-4

New CAS Information Use Policies, enter HELP USAGETERMS for details.

TSCA INFORMATION NOW CURRENT THROUGH JULY 14, 2005

Please note that search-term pricing does apply when  
conducting SmartSELECT searches.

\*\*\*\*\*  
\*  
\* The CA roles and document type information have been removed from \*  
\* the IDE default display format and the ED field has been added, \*  
\* effective March 20, 2005. A new display format, IDERL, is now \*  
\* available and contains the CA role and document type information. \*  
\*  
\*\*\*\*\*

Structure search iteration limits have been increased. See HELP SLIMITS  
for details.

Experimental and calculated property data are now available. For more  
information enter HELP PROP at an arrow prompt in the file or refer  
to the file summary sheet on the web at:  
<http://www.cas.org/ONLINE/DBSS/registryss.html>

## FILE LREGISTRY

LREGISTRY IS A STATIC LEARNING FILE

NEW CAS INFORMATION USE POLICIES, ENTER HELP USAGETERMS FOR DETAILS.

FILE USPATFULL

FILE COVERS 1971 TO PATENT PUBLICATION DATE: 15 Sep 2005 (20050915/PD)

FILE LAST UPDATED: 15 Sep 2005 (20050915/ED)

HIGHEST GRANTED PATENT NUMBER: US6944881

HIGHEST APPLICATION PUBLICATION NUMBER: US2005204445

CA INDEXING IS CURRENT THROUGH 15 Sep 2005 (20050915/UPCA)

ISSUE CLASS FIELDS (/INCL) CURRENT THROUGH: 15 Sep 2005 (20050915/PD)

REVISED CLASS FIELDS (/NCL) LAST RELOADED: Aug 2005

USPTO MANUAL OF CLASSIFICATIONS THESAURUS ISSUE DATE: Aug 2005

>>> USPAT2 is now available. USPATFULL contains full text of the <<<  
>>> original, i.e., the earliest published granted patents or <<<  
>>> applications. USPAT2 contains full text of the latest US <<<  
>>> publications, starting in 2001, for the inventions covered in <<<  
>>> USPATFULL. A USPATFULL record contains not only the original <<<  
>>> published document but also a list of any subsequent <<<  
>>> publications. The publication number, patent kind code, and <<<  
>>> publication date for all the US publications for an invention <<<  
>>> are displayed in the PI (Patent Information) field of USPATFULL <<<  
>>> records and may be searched in standard search fields, e.g., /PN, <<<  
>>> /PK, etc. <<<

>>> USPATFULL and USPAT2 can be accessed and searched together <<<  
>>> through the new cluster USPATALL. Type FILE USPATALL to <<<  
>>> enter this cluster. <<<  
>>> <<<  
>>> Use USPATALL when searching terms such as patent assignees, <<<  
>>> classifications, or claims, that may potentially change from <<<  
>>> the earliest to the latest publication. <<<

This file contains CAS Registry Numbers for easy and accurate  
substance identification.

FILE USPAT2

FILE COVERS 2001 TO PUBLICATION DATE: 15 Sep 2005 (20050915/PD)

FILE LAST UPDATED: 15 Sep 2005 (20050915/ED)

HIGHEST GRANTED PATENT NUMBER: US2005193552

HIGHEST APPLICATION PUBLICATION NUMBER: US2005204275

CA INDEXING IS CURRENT THROUGH 15 Sep 2005 (20050915/UPCA)

ISSUE CLASS FIELDS (/INCL) CURRENT THROUGH: 15 Sep 2005 (20050915/PD)

REVISED CLASS FIELDS (/NCL) LAST RELOADED: Aug 2005

USPTO MANUAL OF CLASSIFICATIONS THESAURUS ISSUE DATE: Aug 2005

USPAT2 is a companion file to USPATFULL. USPAT2 contains full text  
of the latest US publications, starting in 2001, for the inventions  
covered in USPATFULL. USPATFULL contains full text of the original  
published US patents from 1971 to date and the original applications  
from 2001. In addition, a USPATFULL record for an invention contains  
a complete list of publications that may be searched in standard  
search fields, e.g., /PN, /PK, etc.

USPATFULL and USPAT2 can be accessed and searched together through  
the new cluster USPATALL. Type FILE USPATALL to enter this cluster.

Use USPATALL when searching terms such as patent assignees,  
classifications, or claims, that may potentially change from the

earliest to the latest publication.

#### FILE TOXCENTER

FILE COVERS 1907 TO 20 Sep 2005 (20050920/ED)

This file contains CAS Registry Numbers for easy and accurate substance identification.

New CAS Information Use Policies, enter HELP USAGETERMS for details.

TOXCENTER has been enhanced with new files segments and search fields. See HELP CONTENT for more information.

TOXCENTER thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2005 vocabulary. See <http://www.nlm.nih.gov/mesh/> and [http://www.nlm.nih.gov/pubs/techbull/nd04/nd04\\_mesh.html](http://www.nlm.nih.gov/pubs/techbull/nd04/nd04_mesh.html) for a description of changes.

#### FILE MEDLINE

FILE LAST UPDATED: 17 SEP 2005 (20050917/UP). FILE COVERS 1950 TO DATE.

On December 19, 2004, the 2005 MeSH terms were loaded.

The MEDLINE reload for 2005 is now available. For details enter HELP RLOAD at an arrow prompt (=>). See also:

<http://www.nlm.nih.gov/mesh/>  
[http://www.nlm.nih.gov/pubs/techbull/nd04/nd04\\_mesh.html](http://www.nlm.nih.gov/pubs/techbull/nd04/nd04_mesh.html)

OLDMEDLINE now back to 1950.

MEDLINE thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2005 vocabulary.

This file contains CAS Registry Numbers for easy and accurate substance identification.

#### FILE EMBASE

FILE COVERS 1974 TO 15 Sep 2005 (20050915/ED)

EMBASE has been reloaded. Enter HELP RLOAD for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

#### FILE BIOSIS

FILE COVERS 1969 TO DATE.

CAS REGISTRY NUMBERS AND CHEMICAL NAMES (CNs) PRESENT FROM JANUARY 1969 TO DATE.

RECORDS LAST ADDED: 14 September 2005 (20050914/ED)

FILE RELOADED: 19 October 2003.

#### FILE CANCERLIT

FILE COVERS 1963 TO 15 Nov 2002 (20021115/ED)

On July 28, 2002, CANCERLIT was reloaded. See HELP RLOAD for details.

CANCERLIT thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2002 vocabulary. Enter HELP THESAURUS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

## FILE PASCAL

FILE LAST UPDATED: 19 SEP 2005 &lt;20050919/UP&gt;

FILE COVERS 1977 TO DATE.

>>> SIMULTANEOUS LEFT AND RIGHT TRUNCATION IS AVAILABLE  
IN THE BASIC INDEX (/BI) FIELD <<<

## FILE JICST-EPLUS

FILE COVERS 1985 TO 19 SEP 2005 (20050919/ED)

THE JICST-EPLUS FILE HAS BEEN RELOADED TO REFLECT THE 1999 CONTROLLED  
TERM (/CT) THESAURUS RELOAD.

## FILE DRUGU

FILE LAST UPDATED: 20 SEP 2005 &lt;20050920/UP&gt;

&gt;&gt;&gt; DERWENT DRUG FILE (SUBSCRIBER) &lt;&lt;&lt;

>>> FILE COVERS 1983 TO DATE <<<  
>>> THESAURUS AVAILABLE IN /CT <<<

## FILE BIOTECHNO

FILE LAST UPDATED: 7 JAN 2004 &lt;20040107/UP&gt;

FILE COVERS 1980 TO 2003.

&gt;&gt;&gt; BIOTECHNO IS NO LONGER BEING UPDATED AS OF 2004 &lt;&lt;&lt;

>>> SIMULTANEOUS LEFT AND RIGHT TRUNCATION AVAILABLE IN  
/CT AND BASIC INDEX <<<

## FILE BIOTECHDS

FILE LAST UPDATED: 14 SEP 2005 &lt;20050914/UP&gt;

&gt;&gt;&gt; USE OF THIS FILE IS LIMITED TO BIOTECH SUBSCRIBERS &lt;&lt;&lt;

&gt;&gt;&gt; NEW CLASSIFICATION SYSTEM FROM 2002 ONWARDS - SEE HELP CLA &lt;&lt;&lt;

>>> NEW DISPLAY FIELDS LS AND LS2 (LEGAL STATUS DATA FROM  
THE INPADOC DATABASE) AVAILABLE - SEE NEWS <<<

## FILE SCISEARCH

FILE COVERS 1974 TO 15 Sep 2005 (20050915/ED)

SCISEARCH has been reloaded, see HELP RLOAD for details.

## FILE CONF

FILE LAST UPDATED: 16 SEP 2005 &lt;20050916/UP&gt;

FILE COVERS 1976 TO DATE.

## FILE CONFSCI

FILE COVERS 1973 TO 25 May 2005 (20050525/ED)

## FILE DISSABS



FILE COVERS 1861 TO 26 AUG 2005 (20050826/ED)

Only fair use as provided by the United States copyright law is permitted. PROQUEST INFORMATION AND LEARNING COMPANY MAKES NO WARRANTY REGARDING THE ACCURACY, COMPLETENESS OR TIMELINESS OF THE LICENSED MATERIALS OR ANY WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND SHALL NOT BE LIABLE FOR DAMAGES OF ANY KIND OR LOST PROFITS OR OTHER CLAIMS RELATED TO THE LICENSED MATERIALS OR THEIR USE.

=>

This Page Blank (uspto)

=> fil lreg  
FILE 'LREGISTRY' ENTERED AT 15:30:14 ON 20 SEP 2005  
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 1985 AMERICAN CHEMICAL SOCIETY (ACS)

LREGISTRY IS A STATIC LEARNING FILE

NEW CAS INFORMATION USE POLICIES, ENTER HELP USAGETERMS FOR DETAILS.

=> fil reg  
FILE 'REGISTRY' ENTERED AT 15:30:17 ON 20 SEP 2005  
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2005 American Chemical Society (ACS)

Property values tagged with IC are from the ZIC/VINITI data file  
provided by InfoChem.

STRUCTURE FILE UPDATES: 19 SEP 2005 HIGHEST RN 863478-08-4  
DICTIONARY FILE UPDATES: 19 SEP 2005 HIGHEST RN 863478-08-4

New CAS Information Use Policies, enter HELP USAGETERMS for details.

TSCA INFORMATION NOW CURRENT THROUGH JULY 14, 2005

Please note that search-term pricing does apply when  
conducting SmartSELECT searches.

\*\*\*\*\*  
\*  
\* The CA roles and document type information have been removed from \*  
\* the IDE default display format and the ED field has been added, \*  
\* effective March 20, 2005. A new display format, IDERL, is now \*  
\* available and contains the CA role and document type information. \*  
\*  
\*\*\*\*\*

Structure search iteration limits have been increased. See HELP SLIMITS  
for details.

Experimental and calculated property data are now available. For more  
information enter HELP PROP at an arrow prompt in the file or refer  
to the file summary sheet on the web at:  
<http://www.cas.org/ONLINE/DBSS/registryss.html>

=> fil zcap  
FILE 'ZCAPLUS' ENTERED AT 15:30:20 ON 20 SEP 2005  
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

Copyright of the articles to which records in this database refer is  
held by the publishers listed in the PUBLISHER (PB) field (available  
for records published or updated in Chemical Abstracts after December  
26, 1996), unless otherwise indicated in the original publications.  
The CA Lexicon is the copyrighted intellectual property of the  
American Chemical Society and is provided to assist you in searching  
databases on STN. Any dissemination, distribution, copying, or storing

of this information, without the prior written consent of CAS is strictly prohibited.

FILE COVERS 1907 - 20 Sep 2005 VOL 143 ISS 13  
FILE LAST UPDATED: 19 Sep 2005 (20050919/ED)

New CAS Information Use Policies, enter HELP USAGETERMS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> fil uspatfull

FILE 'USPATFULL' ENTERED AT 15:30:25 ON 20 SEP 2005  
CA INDEXING COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

FILE COVERS 1971 TO PATENT PUBLICATION DATE: 15 Sep 2005 (20050915/PD)  
FILE LAST UPDATED: 15 Sep 2005 (20050915/ED)  
HIGHEST GRANTED PATENT NUMBER: US6944881  
HIGHEST APPLICATION PUBLICATION NUMBER: US2005204445  
CA INDEXING IS CURRENT THROUGH 15 Sep 2005 (20050915/UPCA)  
ISSUE CLASS FIELDS (/INCL) CURRENT THROUGH: 15 Sep 2005 (20050915/PD)  
REVISED CLASS FIELDS (/NCL) LAST RELOADED: Aug 2005  
USPTO MANUAL OF CLASSIFICATIONS THESAURUS ISSUE DATE: Aug 2005

>>> USPAT2 is now available. USPATFULL contains full text of the <<<  
>>> original, i.e., the earliest published granted patents or <<<  
>>> applications. USPAT2 contains full text of the latest US <<<  
>>> publications, starting in 2001, for the inventions covered in <<<  
>>> USPATFULL. A USPATFULL record contains not only the original <<<  
>>> published document but also a list of any subsequent <<<  
>>> publications. The publication number, patent kind code, and <<<  
>>> publication date for all the US publications for an invention <<<  
>>> are displayed in the PI (Patent Information) field of USPATFULL <<<  
>>> records and may be searched in standard search fields, e.g., /PN, <<<  
>>> /PK, etc. <<<

>>> USPATFULL and USPAT2 can be accessed and searched together <<<  
>>> through the new cluster USPATALL. Type FILE USPATALL to <<<  
>>> enter this cluster. <<<  
>>> <<<  
>>> Use USPATALL when searching terms such as patent assignees, <<<  
>>> classifications, or claims, that may potentially change from <<<  
>>> the earliest to the latest publication. <<<

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> fil uspat2

FILE 'USPAT2' ENTERED AT 15:30:29 ON 20 SEP 2005  
CA INDEXING COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

FILE COVERS 2001 TO PUBLICATION DATE: 20 Sep 2005 (20050920/PD)  
FILE LAST UPDATED: 20 Sep 2005 (20050920/ED)  
HIGHEST GRANTED PATENT NUMBER: US2005159081  
HIGHEST APPLICATION PUBLICATION NUMBER: US2005204275  
CA INDEXING IS CURRENT THROUGH 20 Sep 2005 (20050920/UPCA)  
ISSUE CLASS FIELDS (/INCL) CURRENT THROUGH: 20 Sep 2005 (20050920/PD)  
REVISED CLASS FIELDS (/NCL) LAST RELOADED: Aug 2005  
USPTO MANUAL OF CLASSIFICATIONS THESAURUS ISSUE DATE: Aug 2005

USPAT2 is a companion file to USPATFULL. USPAT2 contains full text of the latest US publications, starting in 2001, for the inventions covered in USPATFULL. USPATFULL contains full text of the original published US patents from 1971 to date and the original applications from 2001. In addition, a USPATFULL record for an invention contains a complete list of publications that may be searched in standard search fields, e.g., /PN, /PK, etc.

USPATFULL and USPAT2 can be accessed and searched together through the new cluster USPATALL. Type FILE USPATALL to enter this cluster.

Use USPATALL when searching terms such as patent assignees, classifications, or claims, that may potentially change from the earliest to the latest publication.

=> fil hcap

FILE 'HCAPLUS' ENTERED AT 15:30:32 ON 20 SEP 2005

USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.

PLEASE SEE "HELP USAGETERMS" FOR DETAILS.

COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

Copyright of the articles to which records in this database refer is held by the publishers listed in the PUBLISHER (PB) field (available for records published or updated in Chemical Abstracts after December 26, 1996), unless otherwise indicated in the original publications. The CA Lexicon is the copyrighted intellectual property of the the American Chemical Society and is provided to assist you in searching databases on STN. Any dissemination, distribution, copying, or storing of this information, without the prior written consent of CAS, is strictly prohibited.

FILE COVERS 1907 - 20 Sep 2005 VOL 143 ISS 13

FILE LAST UPDATED: 19 Sep 2005 (20050919/ED)

New CAS Information Use Policies, enter HELP USAGETERMS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> fil toxcenter

FILE 'TOXCENTER' ENTERED AT 15:30:36 ON 20 SEP 2005

COPYRIGHT (C) 2005 ACS

FILE COVERS 1907 TO 20 Sep 2005 (20050920/ED)

This file contains CAS Registry Numbers for easy and accurate substance identification.

New CAS Information Use Policies, enter HELP USAGETERMS for details.

TOXCENTER has been enhanced with new files segments and search fields. See HELP CONTENT for more information.

TOXCENTER thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2005 vocabulary. See <http://www.nlm.nih.gov/mesh/> and [http://www.nlm.nih.gov/pubs/techbull/nd04/nd04\\_mesh.html](http://www.nlm.nih.gov/pubs/techbull/nd04/nd04_mesh.html) for a description of changes.

=> fil wpix

FILE 'WPIX' ENTERED AT 15:30:39 ON 20 SEP 2005  
COPYRIGHT (C) 2005 THE THOMSON CORPORATION

FILE LAST UPDATED: 15 SEP 2005 <20050915/UP>  
MOST RECENT DERWENT UPDATE: 200559 <200559/DW>  
DERWENT WORLD PATENTS INDEX SUBSCRIBER FILE, COVERS 1963 TO DATE

>>> FOR A COPY OF THE DERWENT WORLD PATENTS INDEX STN USER GUIDE,  
PLEASE VISIT:  
[http://www.stn-international.de/training\\_center/patents/stn\\_guide.pdf](http://www.stn-international.de/training_center/patents/stn_guide.pdf) <<<

>>> FOR DETAILS OF THE PATENTS COVERED IN CURRENT UPDATES, SEE  
<http://thomsonderwent.com/coverage/latestupdates/> <<<

>>> FOR INFORMATION ON ALL DERWENT WORLD PATENTS INDEX USER  
GUIDES, PLEASE VISIT:  
<http://thomsonderwent.com/support/userguides/> <<<

>>> NEW! FAST-ALERTING ACCESS TO NEWLY-PUBLISHED PATENT  
DOCUMENTATION NOW AVAILABLE IN DERWENT WORLD PATENTS INDEX  
FIRST VIEW - FILE WPIFV.  
FOR FURTHER DETAILS: <http://www.thomsonderwent.com/dwpifv> <<<

>>> THE CPI AND EPI MANUAL CODES HAVE BEEN REVISED FROM UPDATE 200501.  
PLEASE CHECK:  
<http://thomsonderwent.com/support/dwpiref/reftools/classification/code-revision/>  
FOR DETAILS. <<<

=> fil medline

FILE 'MEDLINE' ENTERED AT 15:30:45 ON 20 SEP 2005

FILE LAST UPDATED: 17 SEP 2005 (20050917/UP). FILE COVERS 1950 TO DATE.

On December 19, 2004, the 2005 MeSH terms were loaded.

The MEDLINE reload for 2005 is now available. For details enter HELP  
RLOAD at an arrow prompt (=>). See also:

<http://www.nlm.nih.gov/mesh/>  
[http://www.nlm.nih.gov/pubs/techbull/nd04/nd04\\_mesh.html](http://www.nlm.nih.gov/pubs/techbull/nd04/nd04_mesh.html)

OLDMEDLINE now back to 1950.

MEDLINE thesauri in the /CN, /CT, and /MN fields incorporate the  
MeSH 2005 vocabulary.

This file contains CAS Registry Numbers for easy and accurate  
substance identification.

=> fil biosis

FILE 'BIOSIS' ENTERED AT 15:30:48 ON 20 SEP 2005  
Copyright (c) 2005 The Thomson Corporation

FILE COVERS 1969 TO DATE.  
CAS REGISTRY NUMBERS AND CHEMICAL NAMES (CNs) PRESENT  
FROM JANUARY 1969 TO DATE.

RECORDS LAST ADDED: 14 September 2005 (20050914/ED)

FILE RELOADED: 19 October 2003.

=> fil cancerlit

FILE 'CANCERLIT' ENTERED AT 15:30:52 ON 20 SEP 2005

FILE COVERS 1963 TO 15 Nov 2002 (20021115/ED)

On July 28, 2002, CANCERLIT was reloaded. See HELP RLOAD for details.

CANCERLIT thesauri in the /CN, /CT, and /MN fields incorporate the MeSH 2002 vocabulary. Enter HELP THESAURUS for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> fil embase

FILE 'EMBASE' ENTERED AT 15:30:54 ON 20 SEP 2005

COPYRIGHT (C) 2005 Elsevier Inc. All rights reserved.

FILE COVERS 1974 TO 15 Sep 2005 (20050915/ED)

EMBASE has been reloaded. Enter HELP RLOAD for details.

This file contains CAS Registry Numbers for easy and accurate substance identification.

=> fil pascal

FILE 'PASCAL' ENTERED AT 15:30:58 ON 20 SEP 2005

Any reproduction or dissemination in part or in full,  
by means of any process and on any support whatsoever  
is prohibited without the prior written agreement of INIST-CNRS.  
COPYRIGHT (C) 2005 INIST-CNRS. All rights reserved.

FILE LAST UPDATED: 19 SEP 2005 <20050919/UP>  
FILE COVERS 1977 TO DATE.

>>> SIMULTANEOUS LEFT AND RIGHT TRUNCATION IS AVAILABLE  
IN THE BASIC INDEX (/BI) FIELD <<<

=> fil jicst

FILE 'JICST-EPLUS' ENTERED AT 15:31:01 ON 20 SEP 2005

COPYRIGHT (C) 2005 Japan Science and Technology Agency (JST)

FILE COVERS 1985 TO 19 SEP 2005 (20050919/ED)

THE JICST-EPLUS FILE HAS BEEN RELOADED TO REFLECT THE 1999 CONTROLLED  
TERM (/CT) THESAURUS RELOAD.

=> fil drugu

FILE 'DRUGU' ENTERED AT 15:31:04 ON 20 SEP 2005

COPYRIGHT (C) 2005 THE THOMSON CORPORATION

FILE LAST UPDATED: 20 SEP 2005 <20050920/UP>  
>>> DERWENT DRUG FILE (SUBSCRIBER) <<<

>>> FILE COVERS 1983 TO DATE <<<  
>>> THESAURUS AVAILABLE IN /CT <<<

=> fil biotechno

FILE 'BIOTECHNO' ENTERED AT 15:31:09 ON 20 SEP 2005

COPYRIGHT (C) 2005 Elsevier Science B.V., Amsterdam. All rights reserved.

FILE LAST UPDATED: 7 JAN 2004

<20040107/UP>

FILE COVERS 1980 TO 2003.

>>> BIOTECHNO IS NO LONGER BEING UPDATED AS OF 2004 <<<

>>> SIMULTANEOUS LEFT AND RIGHT TRUNCATION AVAILABLE IN  
/CT AND BASIC INDEX <<<

=> fil biotechds

FILE 'BIOTECHDS' ENTERED AT 15:31:13 ON 20 SEP 2005

COPYRIGHT (C) 2005 THE THOMSON CORPORATION

FILE LAST UPDATED: 14 SEP 2005

<20050914/UP>

>>> USE OF THIS FILE IS LIMITED TO BIOTECH SUBSCRIBERS <<<

>>> NEW CLASSIFICATION SYSTEM FROM 2002 ONWARDS - SEE HELP CLA <<<

>>> NEW DISPLAY FIELDS LS AND LS2 (LEGAL STATUS DATA FROM  
THE INPADOC DATABASE) AVAILABLE - SEE NEWS <<<

=> fil scisearch

FILE 'SCISEARCH' ENTERED AT 15:31:22 ON 20 SEP 2005

Copyright (c) 2005 The Thomson Corporation

FILE COVERS 1974 TO 15 Sep 2005 (20050915/ED)

SCISEARCH has been reloaded, see HELP RLOAD for details.

=> fil conf

FILE 'CONF' ENTERED AT 15:31:24 ON 20 SEP 2005

COPYRIGHT (c) 2005 FIZ Karlsruhe

FILE LAST UPDATED: 16 SEP 2005

<20050916/UP>

FILE COVERS 1976 TO DATE.

=> fil confsci

FILE 'CONFSCI' ENTERED AT 15:31:29 ON 20 SEP 2005

COPYRIGHT (C) 2005 Cambridge Scientific Abstracts (CSA)

FILE COVERS 1973 TO 25 May 2005 (20050525/ED)

=> fil dissabs

FILE 'DISSABS' ENTERED AT 15:31:36 ON 20 SEP 2005

COPYRIGHT (C) 2005 ProQuest Information and Learning Company; All Rights Reserved.

FILE COVERS 1861 TO 26 AUG 2005 (20050826/ED)

Only fair use as provided by the United States copyright law is permitted. PROQUEST INFORMATION AND LEARNING COMPANY MAKES NO WARRANTY REGARDING THE ACCURACY, COMPLETENESS OR TIMELINESS OF THE LICENSED MATERIALS OR ANY WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND SHALL NOT BE LIABLE FOR DAMAGES OF ANY KIND OR LOST PROFITS OR OTHER CLAIMS RELATED TO THE LICENSED



MATERIALS OR THEIR USE.

=> file stnguide

FILE 'STNGUIDE' ENTERED AT 15:31:38 ON 20 SEP 2005

USE IS SUBJECT TO THE TERMS OF YOUR CUSTOMER AGREEMENT


COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY, JAPAN SCIENCE

AND TECHNOLOGY CORPORATION, AND FACHINFORMATIONSZENTRUM KARLSRUHE

FILE CONTAINS CURRENT INFORMATION.

LAST RELOADED: Sep 16, 2005 (20050916/UP).

=> d que stat l13  
L11 STR



VAR G2=O/N  
NODE ATTRIBUTES:  
DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED


GRAPH ATTRIBUTES:  
RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 13

STEREO ATTRIBUTES: NONE  
L13 832 SEA FILE=REGISTRY SSS FUL L11

100.0% PROCESSED 56776 ITERATIONS  
SEARCH TIME: 00.00.03

832 ANSWERS


=> d que stat l18  
L11 STR



VAR G2=O/N  
NODE ATTRIBUTES:  
DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 13

STEREO ATTRIBUTES: NONE  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR



VAR G1=CH/14  
 VAR G2=O/NH  
 VAR G3=NH3/16/19/22  
 NODE ATTRIBUTES:  
 CONNECT IS E1 RC AT 10  
 DEFAULT MLEVEL IS ATOM  
 DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
 RING(S) ARE ISOLATED OR EMBEDDED  
 NUMBER OF NODES IS 24

STEREO ATTRIBUTES: NONE  
 L18      553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16

100.0% PROCESSED      832 ITERATIONS  
 SEARCH TIME: 00.00.01

553 ANSWERS

=> d 164 1-17  
 L64      ANALYZE L18 1- LC :      17 TERMS

TERM #	# OCC	# DOC	% DOC LC
1	535	535	96.75 CA
2	535	535	96.75 CAPLUS
3	255	255	46.11 TOXCENTER
4	115	115	20.80 USPATFULL
5	29	29	5.24 USPAT2
6	3	3	0.54 BIOSIS
7	3	3	0.54 CHEMLIST
8	3	3	0.54 MEDLINE
9	2	2	0.36 CASREACT
10	2	2	0.36 DIOGENES
11	2	2	0.36 IPA
12	2	2	0.36 TSCA
13	1	1	0.18 BIOBUSINESS
14	1	1	0.18 CANCERLIT
15	1	1	0.18 CHEMINFORMRX

16 1 1 0.18 PIRA  
17 1 1 0.18 PROMT  
\*\*\*\*\* END OF L64\*\*\*

=> d que nos l33

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR  
L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L24 715 SEA FILE=HCAPLUS ABB=ON PLU=ON L18  
L25 56930 SEA FILE=HCAPLUS ABB=ON PLU=ON IMMUNOASSAY+PFT,NT/CT  
L26 47927 SEA FILE=HCAPLUS ABB=ON PLU=ON "IMMUNOCHEMICAL ANALYSIS (L)  
IMMUNOASSAY"+PFT,NT/CT  
L27 5294 SEA FILE=HCAPLUS ABB=ON PLU=ON AGGLUTINATION+PFT,NT/CT  
L28 17 SEA FILE=HCAPLUS ABB=ON PLU=ON L24 AND (L25 OR L26 OR L27)  
L29 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E  
LISA OR RIA  
L31 16 SEA FILE=HCAPLUS ABB=ON PLU=ON L24 (L) L29  
L33 23 SEA FILE=HCAPLUS ABB=ON PLU=ON L28 OR L31

=> d his l39

(FILE 'USPATFULL, USPAT2' ENTERED AT 14:41:23 ON 20 SEP 2005)  
L39 3 S L35 AND L38

=> d que nos l39

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR  
L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L35 78 SEA L18  
L38 23444 SEA G01N033-5?/IPC  
L39 3 SEA L35 AND L38

=> d que nos l62

L40 2142 SEA FILE=WPIX ABB=ON PLU=ON (B415 (P) B701 (P) B713 (P) B815  
(P) B831 (P) H1 (P) J011)/M0,M1,M2,M3,M4,M5,M6  
L50 408 SEA FILE=WPIX ABB=ON PLU=ON C08F030-02/IPC  
L59 58 SEA FILE=WPIX ABB=ON PLU=ON L40 AND (G01N033-53?/ICM,ICS)  
L62 2 SEA FILE=WPIX ABB=ON PLU=ON L59 AND L50

=> d que nos l66

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR  
L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L63 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E  
LISA OR RIA OR ?COAG?  
L65 260 SEA FILE=TOXCENTER ABB=ON PLU=ON L18  
L66 31 SEA FILE=TOXCENTER ABB=ON PLU=ON L65 AND L63

=> d que nos l72

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR

L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L63 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E  
LISA OR RIA OR ?COAG?  
L67 55 SEA FILE=MEDLINE ABB=ON PLU=ON L18  
L68 282436 SEA FILE=MEDLINE ABB=ON PLU=ON IMMUNOASSAY+PFT,NT/CT  
L69 6576 SEA FILE=MEDLINE ABB=ON PLU=ON AGGLUTINATION+PFT,NT/CT  
L70 3 SEA FILE=MEDLINE ABB=ON PLU=ON L67 AND (L68 OR L69)  
L71 21 SEA FILE=MEDLINE ABB=ON PLU=ON L67 AND L63  
L72 21 SEA FILE=MEDLINE ABB=ON PLU=ON (L70 OR L71)

=> d que nos 173

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR  
L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L73 0 SEA FILE=EMBASE ABB=ON PLU=ON L18

=> d his 175

(FILE 'BIOSIS, CANCERLIT' ENTERED AT 15:20:16 ON 20 SEP 2005)  
L75 13 S L74 AND L63

=> d que nos 175

L11 STR  
L13 832 SEA FILE=REGISTRY SSS FUL L11  
L16 STR  
L18 553 SEA FILE=REGISTRY SUB=L13 SSS FUL L16  
L63 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E  
LISA OR RIA OR ?COAG?  
L74 57 SEA L18  
L75 13 SEA L74 AND L63

=> dup rem 133 139 162 166 172 173 175

L73 HAS NO ANSWERS

FILE 'HCAPLUS' ENTERED AT 15:33:39 ON 20 SEP 2005

USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.

PLEASE SEE "HELP USAGETERMS" FOR DETAILS.

COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

FILE 'USPATFULL' ENTERED AT 15:33:39 ON 20 SEP 2005

CA INDEXING COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY (ACS)

FILE 'WPIX' ENTERED AT 15:33:39 ON 20 SEP 2005

COPYRIGHT (C) 2005 THE THOMSON CORPORATION

FILE 'TOXCENTER' ENTERED AT 15:33:39 ON 20 SEP 2005

COPYRIGHT (C) 2005 ACS

FILE 'MEDLINE' ENTERED AT 15:33:39 ON 20 SEP 2005

FILE 'BIOSIS' ENTERED AT 15:33:39 ON 20 SEP 2005

Copyright (c) 2005 The Thomson Corporation

FILE 'CANCERLIT' ENTERED AT 15:33:39 ON 20 SEP 2005

PROCESSING COMPLETED FOR L33

PROCESSING COMPLETED FOR L39

PROCESSING COMPLETED FOR L62

PROCESSING COMPLETED FOR L66  
PROCESSING COMPLETED FOR L72  
PROCESSING COMPLETED FOR L73  
PROCESSING COMPLETED FOR L75

L82           70 DUP REM L33 L39 L62 L66 L72 L73 L75 (23 DUPLICATES REMOVED)  
              ANSWERS '1-23' FROM FILE HCAPLUS  
              ANSWERS '24-26' FROM FILE USPATFULL  
              ANSWER '27' FROM FILE WPIX  
              ANSWERS '28-57' FROM FILE TOXCENTER  
              ANSWERS '58-67' FROM FILE MEDLINE  
              ANSWERS '68-70' FROM FILE BIOSIS

=> file stnguide

FILE 'STNGUIDE' ENTERED AT 15:34:17 ON 20 SEP 2005  
USE IS SUBJECT TO THE TERMS OF YOUR CUSTOMER AGREEMENT  
COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY, JAPAN SCIENCE  
AND TECHNOLOGY CORPORATION, AND FACHINFORMATIONSZENTRUM KARLSRUHE

FILE CONTAINS CURRENT INFORMATION.  
LAST RELOADED: Sep 16, 2005 (20050916/UP).

=> d ibib ed ab hitind hitstr retable

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

L82 ANSWER 1 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 2

ACCESSION NUMBER: 2004:264055 HCAPLUS

DOCUMENT NUMBER: 141:3747

TITLE: Evaluation of 2-methacryloyloxyethyl phosphorylcholine polymeric nanoparticle for immunoassay of C-reactive protein detection

AUTHOR(S): Park, Jongwon; Kurosawa, Shigeru; Watanabe, Junji; Ishihara, Kazuhiko

CORPORATE SOURCE: Department of Materials Engineering, School of Engineering, University of Tokyo, Tokyo, Japan

SOURCE: Analytical Chemistry (2004), 76(9), 2649-2655  
CODEN: ANCHAM; ISSN: 0003-2700

PUBLISHER: American Chemical Society

DOCUMENT TYPE: Journal

LANGUAGE: English

ED Entered STN: 01 Apr 2004

AB To prepare novel 2-methacryloyloxyethyl phosphorylcholine (MPC)-polymeric nanoparticle (MPC-PNP), water-soluble amphiphilic phospholipid polymer, poly [MPC-co-Bu methacrylate (BMA)-co-p-nitrophenyloxycarbonyl poly(ethylene glycol) methacrylate (MEONP) (PMBN)], which has active ester groups for bioconjugation on the side chains, was synthesized. MPC-PNP was prepared by a solvent evaporation technique where the poly(L-lactic acid) was used as core and PMBN was applied as an emulsifier and a surface modifier under systematical design of well-arranged phospholipids polar groups in its surface. Characteristics for MPC-PNP were thoroughly investigated with dynamic light scattering, electrophoresis light scattering, XPS, and field emission SEM measurements. Through a protein adsorption test, the phosphorylcholine group on the surface of MPC-PNPs, which had their active ester groups substituted by glycine, were shown to suppress the nonspecific adsorption of bovine serum albumin. These particles were used for C-reactive protein (CRP) detection, where anti-CRP monoclonal antibodies were immobilized on the MPC-PNP using the active ester group, while the remaining active ester groups were thoroughly reacted with glycine. The detection limit about serum-free CRP in the calibration curve was shown to extend from 0.01 to 10 mg/dL when anti-CRP antibody immobilized MPC-PNP was used for serum-free CRP detection. This compares favorably with measurement using polystyrene nanoparticles that were shown to detect from 0.1 to 10 mg/dL by an immunoagglutination technique. Also, for the detection of CRP in serum, MPC-PNP was shown to give the same calibration curve explained by the efficient suppression of nonspecific binding. Furthermore, denaturation of immobilizing anti-CRP antibody on the MPC-PNP hardly occurred despite increasing the temperature. It is concluded that MPC-PNP is unique due to the design of its interfacial properties, also it will perform well in a diagnostic immunoassay because of its optimized material properties.

CC 9-16 (Biochemical Methods)

Section cross-reference(s): 13, 15

IT **Immunoassay**

Scanning electron microscopy

(2-methacryloyloxyethyl phosphorylcholine polymeric nanoparticle for immunoassay of C-reactive protein detection)

IT **685901-25-1P**

RL: ARU (Analytical role, unclassified); DEV (Device component use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES

(Uses)

(2-methacryloyloxyethyl phosphorylcholine polymeric nanoparticle for  
**immunoassay** of C-reactive protein detection)

IT 685901-25-1P

RL: ARU (Analytical role, unclassified); DEV (Device component use); SPN  
(Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES  
(Uses)(2-methacryloyloxyethyl phosphorylcholine polymeric nanoparticle for  
**immunoassay** of C-reactive protein detection)

RN 685901-25-1 HCAPLUS


CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl  
2-methyl-2-propenoate and  $\alpha$ -(2-methyl-1-oxo-2-propenyl)- $\omega$ -[(4-  
nitrophenoxy)carbonyl]oxy]poly(oxy-1,2-ethanediyl) (9CI) (CA INDEX NAME)

CM 1

CRN 666711-01-9

CMF (C2 H4 O)n C11 H9 N O6


CCI PMS



CM 2

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 3

CRN 97-88-1

CMF C8 H14 O2



RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
----------------------------	---------------	--------------	-------------	--------------------------	--------------------



Aceti, A	1991	22	135	J Infect	MEDLINE
Akhtar, S	1995	107	209	J Membr Sci	HCAPLUS
Babba, H	1994	50	64	Am J Trop Med	MEDLINE
Bangs, L	1996	69	1873	Pure Appl Chem	
Biasucci, L	1999	99	855	Circulation	MEDLINE
Bundy, J	1999	71	1460	Anal Chem	HCAPLUS
Chen, J	2001	17	369	Biotechnol Prog	HCAPLUS
de Winter, R	1999	42	240	Cardiovasc Res	HCAPLUS
Diamandis, E	1990	194	19	Clin Chim Acta	HCAPLUS
Elvassore, N	2001	40	795	Ind Eng Chem Res	HCAPLUS
Ghourchian, H	1997	41	401	Talanta	
Haverkate, F	1997	349	462	Lancet	MEDLINE
Hayward, J	1984	5	135	Biomaterials	HCAPLUS
Hirsch, L	2003	75	2377	Anal Chem	HCAPLUS
Holownia, P	2001	73	3426	Anal Chem	HCAPLUS
Ishihara, K	1999	10	1047	J Biomater Sci, Poly	HCAPLUS
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Ishihara, K	1992	26	1543	J Biomed Mater Res	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Ishihara, K	1990	22	355	Polym J	HCAPLUS
Kitano, H	2002	104	10425	J Phys Chem B	
Kitano, H	2003	19	10260	Langmuir	HCAPLUS
Koenig, W	1999	99	237	Circulation	MEDLINE
Kurosawa, S	1990	38	1117	Chem Pharm Bull	HCAPLUS
Lu, D	1991	3	127	J Biomater Sci, Poly	HCAPLUS
Molina-Bolivar, J	2001	17	2514	Langmuir	HCAPLUS
Pentikanen, M	2000	247	359	J Intern Med	
Perez-Amodio, S	2001	73	3417	Anal Chem	HCAPLUS
Ridker, P	1998	97	425	Circulation	MEDLINE
Ruiz, L	1998	19	987	Biomaterials	HCAPLUS
Sakai-Kato, K	2002	74	2943	Anal Chem	HCAPLUS
Sakaki, S	1999	47	523	J Biomed Mater Res	HCAPLUS
Singer, S	1972	175	720	Science	HCAPLUS
Velev, O	1999	15	3693	Langmuir	HCAPLUS
Wang, D	2002	2	857	Nano Lett	HCAPLUS
Wang, J	2001	17	5739	Langmuir	HCAPLUS
Whicher, J	1999	37	495	Clin Chem Lab Med	HCAPLUS

=> d ibib ed ab hitind hitstr retable 2-23

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

L82 ANSWER 2 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 3

ACCESSION NUMBER: 2004:1004654 HCAPLUS

DOCUMENT NUMBER: 142:370012

TITLE: Evaluation of a high-affinity QCM immunosensor using antibody fragmentation and 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer

AUTHOR(S): Kurosawa, Shigeru; Nakamura, Miki; Park, Jong-Won; Aizawa, Hidenobu; Yamada, Kazunori; Hirata, Mitsuo

CORPORATE SOURCE: National Institute of Advanced Industrial Science and Technology (AIST), Tsukuba, Ibaraki, 305-8565, Japan

SOURCE: Biosensors & Bioelectronics (2004), 20(6), 1134-1139  
CODEN: BBIOE4; ISSN: 0956-5663

PUBLISHER: Elsevier B.V.

DOCUMENT TYPE: Journal

LANGUAGE: English

ED Entered STN: 22 Nov 2004

AB This study evaluated construction of a high affinity quartz crystal microbalance (QCM) immunosensor using anti-C-reactive protein (CRP) antibody and its fragments for CRP detection. Three types of antibody were immobilized on the surface of a QCM via covalent-bonding. Then affinity was evaluated through antigen-antibody binding between CRP and its antibody. Affinity between antigen-antibody was shown to be highest when anti-CRP F(ab')<sub>2</sub>-IgG antibody (70 µg/mL) was immobilized on the QCM. In case of anti-CRP F(ab')<sub>2</sub>-IgG antibody, affinity which was attributable to antigen-antibody binding was almost twice that of anti-CRP IgG antibody, which is used conventionally for QCM immunosensors. In addition, when it was treated with 2-methacryloyloxyethyl phosphorylcholine-co-Bu methacrylate, so-called MPC polymer, highly affinitive and selective immunosensing for CRP was achieved without non-specific binding from plasma proteins in human serum. When anti-CRP F(ab')<sub>2</sub>-IgG antibody was immobilized on the QCM, the detection limit and the linearity of CRP calibration curve were achieved at concns. from 0.001 to 100 µg/dL even during investigation in serum samples. Exptl. results verified the successful construction of a highly affinitive and selective QCM-immunosensor which was modified with anti-CRP F(ab')<sub>2</sub>-IgG antibody and MPC polymer.

CC 9-1 (Biochemical Methods)

IT 67881-98-5, 2-Methacryloyloxyethyl phosphorylcholine

RL: TEM (Technical or engineered material use); USES (Uses)

(evaluation of a high-affinity QCM **immunosensor** using antibody fragmentation and 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer)


IT 67881-98-5, 2-Methacryloyloxyethyl phosphorylcholine

RL: TEM (Technical or engineered material use); USES (Uses)

(evaluation of a high-affinity QCM **immunosensor** using antibody fragmentation and 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Abad, J	1998	70	2848	Anal Chem	HCAPLUS
Aizawa, H	2001	76	173	Sens Actuat B	
Borgue, L	2000	46	1839	Clin Chem	
Bunde, R	1998	46	1223	Talanta	HCAPLUS
Chou, S	2002	48	913	Clin Chem	HCAPLUS
Haverkate, F	1997	349	462	Lancet	MEDLINE
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Ishihara, K	2000	1	131	Sci Technol Adv Mat	HCAPLUS
Kanai, I	1998			CRP (C-reactive prot	
Koenig, W	1999	99	237	Circulation	MEDLINE

Kurosawa, S	1990	38	1117	Chem Pharm Bull	HCAPLUS
Kurosawa, S	2000	47	1256	IEEE Trans UFFC	
Kurosawa, S	2003	14	1882	Meas Sci Technol	HCAPLUS
Ledue, T	1998	35	745	Annu Clin Biochem	HCAPLUS
Morgan, C	1996	42	193	Clin Chem	HCAPLUS
Muratsugu, M	1993	65	2933	Anal Chem	HCAPLUS
Muratsugu, M	1997	1	99	Colloid Interface Sc	HCAPLUS
Park, J	2003	50	193	IEEE Trans UFFC	
Park, J	2003	91	158	Sens Actuat B	
Pentikanen, M	2000	247	359	J Intern Med	
Ridker, P	1998	97	425	Circulation	MEDLINE
Ridker, P	1997	336	973	N Engl J Med	HCAPLUS
Rifai, N	1999	45	2136	Clin Chem	HCAPLUS
Sakai, G	1995	24-25	134	Sens Actuat B	
Storri, S	1998	13	347	Biosens Bioelectron	HCAPLUS
Thompson, M	1986	58	1206	Anal Chem	HCAPLUS
Uttenthaler, E	2001	16	735	Biosen Bioelectron	HCAPLUS
Whicher, J	1999	37	495	Clin Chem Lab Med	HCAPLUS
Wood, W	2000	46	131	Clin Lab	HCAPLUS
Xia, C	1996	336	185	Anal Chim Acta	

L82 ANSWER 3 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 4

ACCESSION NUMBER: 2005:4348 HCAPLUS

DOCUMENT NUMBER: 142:296317

TITLE: Cytokine Adsorptive Property of Various Adsorbents in Immunoabsorption Columns and a Newly Developed Adsorbent: An in vitro Study

AUTHOR(S): Oda, Shigeto; Hirasawa, Hiroyuki; Shiga, Hidetoshi; Nakanishi, Kazuya; Matsuda, Ken-ichi; Nakamura, Masataka; Ikeda, Hiroyuki; Sakai, Masamune

CORPORATE SOURCE: Department of Emergency and Critical Care Medicine, Graduate School of Medicine, Chiba University, Chiba, Japan

SOURCE: Blood Purification (2004), 22(6), 530-536

CODEN: BLPUDO; ISSN: 0253-5068

PUBLISHER: S. Karger AG

DOCUMENT TYPE: Journal

LANGUAGE: English

ED Entered STN: 04 Jan 2005

AB Background/Aims: Cytokines play important roles in the pathophysiol. of systemic inflammatory response syndrome (SIRS) and sepsis. Therefore, some effective measures to remove cytokines from the bloodstream could be effective in the treatment of SIRS and sepsis. The aim of this study was to evaluate the cytokine adsorptive property of various adsorbents for the purpose of the development of new selective cytokine adsorption columns. Methods: The cytokine adsorptive property of adsorbent in a CF-X column, which consists of cellulose beads cross-linked with hexamethylene-di-isocyanate, was compared with those of various adsorbents in currently available immunoabsorption columns, such as Immusorba TR, Immusorba PH, Selesorb, and Lixelle, in vitro batchwise test using patients' plasma. A newly developed adsorbent, MPCF-X, which was modified by coating the surface of the adsorbent in CF-X with 2-methacryloyloxyethyl phosphorylcholine (MPC), was also tested for its cytokine adsorptive property. Results: The adsorbent in CF-X showed a significantly higher adsorption rate for TNF- $\alpha$ , interleukin (IL)-6 and IL-10 compared with other adsorbents. Adsorbent in Lixelle showed good affinity to TNF- $\alpha$  and IL-8. Especially, the adsorbent in CF-X almost completely removed TNF- $\alpha$ , whereas it also had considerable affinity to normal IgG. MPCF-X showed decreased affinity to IgG with considerable adsorptive properties to cytokines. Conclusion: Selective cytokine adsorption

columns could be developed with improvement of currently available adsorbents. Such a new selective cytokine adsorption column could be clin. applied for the treatment of SIRS/sepsis.

CC 15-1 (Immunochemistry)

Section cross-reference(s): 14

IT 67881-98-5, 2-Methacryloyloxyethyl phosphorylcholine


RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(cytokine adsorption by various adsorbents in **immuno**adsorption columns in treatment of sepsis)

IT 67881-98-5, 2-Methacryloyloxyethyl phosphorylcholine

RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(cytokine adsorption by various adsorbents in **immuno**adsorption columns in treatment of sepsis)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



# RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Bellomo, R	1993	21	522	Crit Care Med	MEDLINE
Cain, B	1998	186	337	J Am Coll Surg	MEDLINE
Cole, L	2001	27	978	Intensive Care Med	MEDLINE
De Vriese, A	1999	10	846	J Am Soc Nephrol	MEDLINE
Dinareello, C	1993	269	1829	JAMA	MEDLINE
Heering, P	2003	26	128	Kidney Blood Press R	MEDLINE
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Kutsuki, H	1998	2	18	Ther Apher	HCAPLUS
Matsuda, K	2001	5	306	Ther Apher	HCAPLUS
McMaster, P	2003	4	2	Pediatr Crit Care Me	
Members of American Col	1992	20	864	Crit Care Med	
Nakaji, S	2001	5	301	Ther Apher	HCAPLUS
Oda, S	2002	6	193	Ther Apher	
O'Reilly, M	1999	12	411	Shock	MEDLINE
Ronco, C	2002	30	1250	Crit Care Med	
Ronco, C	2000	76	148	Kidney Int	
Shapiro, L	1993	1	13	New Horiz	MEDLINE
Shetz, M	1995	21	169	Intensive Care Med	
Stegmayr, B	2000	18	149	Blood Purif	MEDLINE
Suzuki, K	2003	7	104	Ther Apher	
Tomyo, M	1996	15	118	Jpn J Apheresis	
Tsuchida, K	2002	10	485	Int J Mol Med	HCAPLUS
Wheeler, A	1999	340	207	N Engl J Med	MEDLINE
Winchester, J	2003	21	79	Blood Purif	HCAPLUS
Winchester, J	2002	137	170	Contrib Nephrol	
Yagihashi, A	2002	6	358	Ther Apher	HCAPLUS
Yoshida, M	1998	2	185	Ther Apher	HCAPLUS

L82 ANSWER 4 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 6

ACCESSION NUMBER: 2003:945791 HCAPLUS

DOCUMENT NUMBER: 140:14529

TITLE: Developing solvent, measuring method, and kit for immunochromatography  
 INVENTOR(S): Mochizuki, Takeshi; Komatsu, Mariko; Sakaki, Shujiro  
 PATENT ASSIGNEE(S): Taunzu K. K., Japan; NOF Corporation  
 SOURCE: Jpn. Kokai Tokkyo Koho, 10 pp.  
 CODEN: JKXXAF  
 DOCUMENT TYPE: Patent  
 LANGUAGE: Japanese  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2003344406	A2	20031203	JP 2002-150996	20020524
PRIORITY APPLN. INFO.:			JP 2002-150996	20020524

ED Entered STN: 04 Dec 2003

AB An improved developing solvent for an immunochromatog. is provided, with which non-specific aggregation and non-specific reaction upon measurements are prevented, and the measurements are performed with high accuracy. The developing solvent for an immunochromatog. is characterized in that it comprises a buffer containing a polymer possessing phosphorylcholine groups. It is preferable that the polymer is contained in the concentration of 0.005-0.3w/v%, and its number average mol. weight is higher than 40,000. The polymer preferably contains 2-methacryloyloxyethylphosphorylcholine as the constituting monomer, and it can be either a homopolymer or a copolymer.

IC ICM G01N033-543  
ICS G01N033-531

CC 9-10 (Biochemical Methods)


IT 9005-64-5, Tween 20 26628-22-8, Sodium azide **67881-98-5D**, 2-Methacryloyloxyethylphosphorylcholine, copolymer with methoxypolyethyleneglycolmonomethacrylate, copolymer with methacrylate **67881-98-5D**, 2-Methacryloyloxyethylphosphorylcholine, homopolymer  
 RL: ARU (Analytical role, unclassified); ANST (Analytical study)  
 (improved developing solvent, measuring method, and kit for **immunochromatog.**)

IT 79-10-7, Acrylic acid, reactions 79-41-4, Methacrylic acid, reactions 25249-16-5, Polyethyleneglycolmonomethacrylate 26403-58-7, Polyethyleneglycolmonoacrylate 26915-72-0, Methoxypolyethyleneglycolmono methacrylate 32171-39-4 **67881-98-5**, 2-Methacryloyloxyethylphosphorylcholine **150120-15-3**  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (improved developing solvent, measuring method, and kit for **immunochromatog.**)


IT **67881-98-5D**, 2-Methacryloyloxyethylphosphorylcholine, copolymer with methoxypolyethyleneglycolmonomethacrylate, copolymer with methacrylate  
 RL: ARU (Analytical role, unclassified); ANST (Analytical study)  
 (improved developing solvent, measuring method, and kit for **immunochromatog.**)

RN 67881-98-5 HCAPLUS


CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



IT 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine  
 150120-15-3  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (improved developing solvent, measuring method, and kit for  
 immunochromatog.)  
 RN 67881-98-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
 tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RN 150120-15-3 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N-trimethyl-9-  
 oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



L82 ANSWER 5 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 8  
 ACCESSION NUMBER: 2002:958988 HCAPLUS  
 DOCUMENT NUMBER: 138:21783  
 TITLE: Agglutination-promoting agent for antigen or antibody  
 immunoassay  
 INVENTOR(S): Kakuta, Kyoichi; Wada, Hiroshi; Ishihara, Kazuhiko  
 PATENT ASSIGNEE(S): Wako Pure Chemical Industries, Ltd., Japan  
 SOURCE: Jpn. Kokai Tokkyo Koho, 13 pp.  
 CODEN: JKXXAF  
 DOCUMENT TYPE: Patent  
 LANGUAGE: Japanese  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2002365296	A2	20021218	JP 2001-169051	20010605
US 2004157276	A1	20040812	US 2004-626502	20040304
			JP 2001-169051	A 20010605


PRIORITY APPLN. INFO.:

ED Entered STN: 18 Dec 2002


AB Provided are aggregation-promoting compds. for use in agglutination immunoassay. These compds are branched polymers or copolymers having basic (monomer) structure of OPO2-O-R4-N(R1R2R3) where R1-3 are independently H, OH or alkyl group; and R4 is an alkyl group or alkylene group. The agglutination immunoassay reagent comprises carrier- or latex-immobilized antibody or antigen. The agglutination immunoassay is useful for determination of antigen or antibody, e.g. C-reactive protein, rheumatic factor and prostate-specific antigen.

IC ICM G01N033-531

ICS C08F030-02; G01N033-543  
 CC 9-10 (Biochemical Methods)  
 Section cross-reference(s): 15  
 IT **Immunoassay**  
 (agglutination test; agglutination-promoting agent for antigen or antibody immunoassay)  
 IT **Agglutination**  
 (promoters; agglutination-promoting agent for antigen or antibody immunoassay)  
 IT **67881-99-6P**, Poly-2-methacryloyloxyethylphosphorylcholine  
**125275-25-4P**, n-Butyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer **144514-08-9P**, Stearyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer **313216-64-7P**  
**478015-82-6P**  
 RL: ARU (Analytical role, unclassified); DGN (Diagnostic use); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)  
 (agglutination-promoting agent for antigen or antibody immunoassay)  
 IT **67881-99-6P**, Poly-2-methacryloyloxyethylphosphorylcholine  
**125275-25-4P**, n-Butyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer **144514-08-9P**, Stearyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer **313216-64-7P**  
**478015-82-6P**  
 RL: ARU (Analytical role, unclassified); DGN (Diagnostic use); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)  
 (agglutination-promoting agent for antigen or antibody immunoassay)  
 RN 67881-99-6 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, homopolymer (9CI) (CA INDEX NAME)  
 CM 1  
 CRN 67881-98-5  
 CMF C11 H22 N O6 P




RN 125275-25-4 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)  
 CM 1  
 CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2

CRN 97-88-1

CMF C8 H14 O2




RN 144514-08-9 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with octadecyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 32360-05-7

CMF C22 H42 O2



RN 313216-64-7 HCAPLUS


CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with phenylmethyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P






CM 2

CRN 2495-37-6

CMF C11 H12 O2




RN 478015-82-6 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with 2-hydroxy-N,N,N-trimethyl-3-[(2-methyl-1-oxo-2-propenyl)oxy]-1-propanaminium chloride (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 13052-11-4

CMF C10 H20 N O3 . Cl

● Cl<sup>-</sup>

L82 ANSWER 6 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 15

ACCESSION NUMBER: 1999:680723 HCAPLUS

DOCUMENT NUMBER: 132:32839

TITLE: Stabilization of an antibody conjugated with enzyme by 2-methacryloyloxyethyl phosphorylcholine copolymer in enzyme-linked immunosorbent assay

AUTHOR(S): Sakaki, Shujiro; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo Medical and Dental University, Tokyo, 101-0062, Japan

SOURCE: Journal of Biomedical Materials Research (1999), 47(4), 523-528  
CODEN: JBMRBG; ISSN: 0021-9304

PUBLISHER: John Wiley & Sons, Inc.

DOCUMENT TYPE: Journal

LANGUAGE: English

ED Entered STN: 27 Oct 1999

AB The purpose of this study was to develop a novel synthetic stabilizer of enzyme-linked antibody in the ELISA. The water-soluble amphiphilic phospholipid polymer, poly[2-methacryloyloxyethyl phosphorylcholine (MPC)-co-styrene (St)] was synthesized, and its stabilizing functions for the antibody were compared with conventional stabilizers of the antibody conjugated with enzyme (enzyme-antibody conjugate), such as bovine serum albumin (BSA) and casein. In the absence of the stabilizer, the remaining immunol. activity decreased to about 10% of its initial value after 37 days. The same tendency was observed even when the enzyme-antibody conjugate in 1.0 wt % BSA solution was used as a stabilizer. In 1.0 wt % casein solution, the immunol. activity decreased to 29% of the initial value after 37 days. On the other hand, in 0.1 wt % and 1.0 wt % poly(MPC-co-St) solution, the activity remained 74% and 92% of the initial value, resp. The effects of poly(MPC-co-St) on the stabilization of the enzyme-antibody conjugate depended on the concentration of poly(MPC-co-St). During the ELISA procedure, not only did poly(MPC-co-St) have no effect on the reaction between the antigen and the antibody, but it also had no effect on the reaction between the enzyme and the substrate. These results indicate that poly(MPC-co-St) has the ability to suppress the denaturation of protein, enzyme, and antibody. We concluded that water-soluble poly(MPC-co-St) is an effective synthetic stabilizer in the ELISA.

CC 9-10 (Biochemical Methods)

IT **Immunoassay**  
(enzyme-linked immunosorbent assay; stabilization of an antibody conjugated with enzyme by 2-methacryloyloxyethyl phosphorylcholine copolymer in ELISA)

IT **134483-35-5**  
RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use); ANST (Analytical study); USES (Uses)  
(stabilization of an antibody conjugated with enzyme by 2-methacryloyloxyethyl phosphorylcholine copolymer in **ELISA**)

IT **134483-35-5**  
RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use); ANST (Analytical study); USES (Uses)  
(stabilization of an antibody conjugated with enzyme by 2-methacryloyloxyethyl phosphorylcholine copolymer in **ELISA**)


RN 134483-35-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI)  
(CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P



CM 2

CRN 100-42-5

CMF C8 H8

 $\text{H}_2\text{C}=\text{CH}-\text{Ph}$ 

## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Anderton, B	1980	2	122	Immunol Today	
Donovan, J	1975	250	1966	J Biol Chem	HCAPLUS
Farr, A	1981	47	129	J Immunol Meth	HCAPLUS
Ghosh, S	1974	337	395	Biochem Biophys Acta	HCAPLUS
Ikemi, M	1982	15	281	Macromolecules	HCAPLUS
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Ishihara, K	1994	32	859	J Polym Sci Polym Ch	HCAPLUS
Ishihara, K	1990	22	355	Polym J	HCAPLUS
Ishikawa, E	1987	1	238	J Clin Lab Anal	HCAPLUS
Kojima, M	1991	12	121	Biomaterials	HCAPLUS
Naik, D	1975	47	267	Anal Biochem	HCAPLUS
Orci, L	1985	28	528	Diabetologia	HCAPLUS
Osborn, M	1978	77	R27	J Cell Biol	MEDLINE
Sakaki, S	1997		167	Advances in polymeri	
Shaw, D	1980			Introduction to coll	
Wisdom, G	1976	22	1243	Clin Chem	HCAPLUS

L82 ANSWER 7 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2005:300694 HCAPLUS

DOCUMENT NUMBER: 142:351673

TITLE: Automated analytical method and apparatus

INVENTOR(S): Kurosawa, Shigeru; Aizawa, Hidenobu

PATENT ASSIGNEE(S): National Institute of Advanced Industrial Science and Technology, Japan

SOURCE: PCT Int. Appl., 108 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2005031316	A1	20050407	WO 2004-JP14664	20040929

W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD,

GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC,  
 LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NA, NI,  
 NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY,  
 TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW  
 RW: BW, GH, GM, KE, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZW, AM,  
 AZ, BY, KG, KZ, MD, RU, TJ, TM, AT, BE, BG, CH, CY, CZ, DE, DK,  
 EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PL, PT, RO, SE,  
 SI, SK, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE,  
 SN, TD, TG

PRIORITY APPLN. INFO.:

JP 2003-338604

A 20030929

ED Entered STN: 07 Apr 2005

AB An automated anal. method is provided, which comprises making a sample in a container absorbed to a probe, dispensing the sample from the probe to a sensor part equipped with a piezoelec. element for converting the mass change on the sensor into an elec. change such as a basic resonant frequency and quantitating it, performing at least once making the dispensed sample reabsorbed to the probe and re-dispensing the reabsorbed sample to the sensor part, and thereby, promoting the progress of a chemical reaction or else generated on the sensor. Also provided is an automated anal. apparatus used for this method. Diagrams describing the apparatus assembly are given.

IC ICM G01N005-02

ICS G01N001-00; G01N033-543; G01N035-02

CC 9-1 (Biochemical Methods)

IT **Immunoassay**

(sandwich reaction; automated anal. method using apparatus equipped with piezoelec. element)

IT 56-81-5, Glycerol, analysis **67881-98-5D**, 2-Methacryloyloxyethylphosphorylcholine, polymer 658045-21-7, BlockAce 658051-44-6, StabilGuard

RL: ARU (Analytical role, unclassified); ANST (Analytical study)

(automated anal. method using apparatus equipped with piezoelec. element)


IT **67881-98-5D**, 2-Methacryloyloxyethylphosphorylcholine, polymer

RL: ARU (Analytical role, unclassified); ANST (Analytical study)

(automated anal. method using apparatus equipped with piezoelec. element)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Aizawa, K	2003		202	Dai 42 Kai Japan Oil	
Beckman Instruments Inc	1999			JP 11-502937 A	
Beckman Instruments Inc	1999			US 2002182117 A1	HCAPLUS
Beckman Instruments Inc	1999			EP 819256 A1	HCAPLUS
Beckman Instruments Inc	1999			AU 9715304 A	HCAPLUS
Beckman Instruments Inc	1999			WO 9726536 A1	HCAPLUS
Koyama, N	2000			JP 2000283905 A	HCAPLUS
Nakamura, K	2002	18	46	Chemical Sensors, Th	

Olympus Optical Co Ltd	1983			JP 58-196461 A	
Toshiba Corp	1993			JP 05-285000 A	HCAPLUS
Totsuka, K	2002	18	73	Chemical Sensors, Th	
Yase, T	2003		58.1	Sogo Kenkyu Project	

L82 ANSWER 8 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2005:591492 HCAPLUS

DOCUMENT NUMBER: 143:93600

TITLE: Blood-group antibody measurement/identification instrument tool, and method

INVENTOR(S): Ito, Yoshihiro; Yamauchi, Tetsuya; Ishikawa, Yoshihide; Uchikawa, Makoto

PATENT ASSIGNEE(S): Kanagawa Academy of Science and Technology, Japan; Japanese Red Cross Society

SOURCE: Jpn. Kokai Tokkyo Koho, 18 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2005181154	A2	20050707	JP 2003-423771	20031219
PRIORITY APPLN. INFO.:			JP 2003-423771	20031219

ED Entered STN: 08 Jul 2005

AB A blood-group antibody measurement/identification instrument tool is provided, with which the quant. measurement of a blood-group antibody is feasible with the less required quantity of panel blood cells than the traditional blood-group plate without requiring a mech. operation such as centrifugation or else. Also provided is a blood-group antibody measurement/identification method using this instrument tool. The instrument tool is constituted by resp. immobilizing multiple kinds of panel blood cells to a different region on a base body through an immobilization agent consisting of a zwitterionic water-soluble polymer or nonionic water-soluble polymer possessing at least two photo-reactive groups in a mol.

IC ICM G01N033-53

ICS C08F220-36; C08F220-60; C12M001-34; G01N021-76; G01N033-543; G01N033-547

CC 9-10 (Biochemical Methods)

IT **Immunoassay**

(blood-group antibody measurement tool using panel blood cells immobilized through photo-reactive water-soluble polymer immobilization agent)

IT **Immunoassay**

(chemiluminescence, enhanced; blood-group antibody measurement tool using panel blood cells immobilized through photo-reactive water-soluble polymer immobilization agent)

IT 6066-82-6, N-Hydroxysuccinimide 6427-66-3, 4-Azido-benzoic acid 14860-64-1, 4-Azido-aniline 18358-13-9, Methacrylate, reactions 25322-68-3, Polyethyleneglycol 34901-14-9 39927-08-7, Poly(ethyleneglycol)bis(carboxymethyl)ether 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine 67881-98-5D, 2-Methacryloyloxyethylphosphorylcholine, copolymer with photo-reactive methacrylamide derivative 88539-10-0 150120-15-3 158760-93-1 163801-79-4 179637-06-0 179953-15-2 857051-80-0

RL: RCT (Reactant); RACT (Reactant or reagent)

(blood-group antibody measurement tool using panel blood cells


immobilized through photo-reactive water-soluble polymer immobilization agent)

IT 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine  
 67881-98-5D, 2-Methacryloyloxyethylphosphorylcholine, copolymer  
 with photo-reactive methacrylamide derivative 88539-10-0  
 150120-15-3 158760-93-1 179637-06-0  
 179953-15-2

RL: RCT (Reactant); RACT (Reactant or reagent)  
 (blood-group antibody measurement tool using panel blood cells  
 immobilized through photo-reactive water-soluble polymer immobilization agent)


RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 88539-10-0 HCAPLUS

CN 3,5,12-Trioxa-4-phosphapentadec-14-en-1-aminium, 4-hydroxy-N,N,N,14-tetramethyl-13-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 150120-15-3 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N-trimethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 158760-93-1 HCAPLUS

CN 3,5,10-Trioxa-4-phosphatridec-12-en-1-aminium, 4-hydroxy-N,N,N,12-tetramethyl-11-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 179637-06-0 HCAPLUS

CN 3,5,16-Trioxa-4-phosphanonadec-18-en-1-aminium, 4-hydroxy-N,N,N,18-tetramethyl-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RN 179953-15-2 HCAPLUS

CN 3,5-Dioxa-8-aza-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



L82 ANSWER 9 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2005:299600 HCAPLUS

DOCUMENT NUMBER: 142:353877

TITLE: Antigen or antibody determination in body fluid with rapid, simple and highly sensitive immunoassay

INVENTOR(S): Ishihara, Kazuhiko; Watanabe, Junji; Kurosawa, Shigeru

PATENT ASSIGNEE(S): National Institute of Advanced Industrial Science and Technology, Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 12 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2005091236	A2	20050407	JP 2003-326829	20030918
PRIORITY APPLN. INFO.:			JP 2003-326829	20030918

ED Entered STN: 07 Apr 2005

AB Disclosed is a rapid, simple and highly sensitive immunoassay method for detecting antigen and antibody in body fluid. The immunoassay uses polymeric particles or nanoparticles comprising alkyl, hydroxy and aromatic side chain for antigen or antibody immobilization to preserve binding

activity and to reduce nonspecific binding by contaminant proteins. In example, the polymeric nanoparticles were prepared with 2-methacryloyloxyethylphosphorylcholine and Bu methacrylate. The polymeric nanoparticle-immobilized monoclonal antibodies were used for agglutination immunoassay of human C-reactive protein in blood serum.

IC ICM G01N033-543

CC 15-2 (Immunochemistry)

Section cross-reference(s): 9

IT **Immunoassay**

(agglutination test; antigen determination in body fluid with immunoassay comprising polymeric nanoparticle-immobilized antibody)

IT **125275-25-4P**, 2-Methacryloyloxyethylphosphorylcholine-butyl methacrylate copolymer

RL: ARU (Analytical role, unclassified); BSU (Biological study, unclassified); DGN (Diagnostic use); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(antigen determination in body fluid with **immunoassay** comprising polymeric nanoparticle-immobilized antibody)

IT **125275-25-4P**, 2-Methacryloyloxyethylphosphorylcholine-butyl methacrylate copolymer

RL: ARU (Analytical role, unclassified); BSU (Biological study, unclassified); DGN (Diagnostic use); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(antigen determination in body fluid with **immunoassay** comprising polymeric nanoparticle-immobilized antibody)


RN 125275-25-4 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 97-88-1

CMF C8 H14 O2



L82 ANSWER 10 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN  
ACCESSION NUMBER: 2004:857773 HCAPLUS



DOCUMENT NUMBER: 141:346148  
 TITLE: Substance immobilization agent, and substance immobilization method/substrate using agent  
 INVENTOR(S): Yamauchi, Tetsuya; Ito, Yoshihiro; Hasuda, Hirokazu; Konno, Tomohiro; Ishihara, Kazuhiko  
 PATENT ASSIGNEE(S): Kanagawa Academy of Science and Technology, Japan  
 SOURCE: PCT Int. Appl., 32 pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: Japanese  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2004088319	A1	20041014	WO 2004-JP4510	20040330
W:				
AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NA, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW				
RW:				
BW, GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PL, PT, RO, SE, SI, SK, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG				

PRIORITY APPLN. INFO.: JP 2003-93834 A 20030331  
 JP 2003-346560 A 20031006

ED Entered STN: 18 Oct 2004

AB A substance immobilization agent is disclosed, which is capable of immobilizing various substances to be immobilized to a substrate by covalent bond, and is highly effective in preventing non-specific adsorption. The agent comprises a polymer having phosphorylcholine groups and multiple photoreactive groups (e.g., azido groups) in a single mol. The polymer binds through the photoreactive groups with a substrate and a substance to be immobilized, whereby the substance to be immobilized is bonded to the substrate through the polymer by a covalent bond. Furthermore, non-specific adsorption is effectively prevented by the phosphorylcholine groups. Also provided is an immobilization agent, which is used for immobilizing a desired substance to a substrate and comprises a nonionic water-soluble polymer having at least two photoreactive groups in a single mol.

IC ICM G01N033-547

ICS G01N033-53; C12M001-40; C12M001-34

CC 9-16 (Biochemical Methods)

IT **Immunoassay**

(chemiluminescence; substance immobilization agent consisting of polymer with photoreactive groups, and use in substance immobilization method/substrate)

IT **Immunoassay**

(fluorescence; substance immobilization agent consisting of polymer with photoreactive groups, and use in substance immobilization method/substrate)

IT 106-91-2, Glycidylmethacrylate 814-68-6, Acrylic acid chloride  
 6066-82-6, N-Hydroxysuccinimide 6427-66-3, 4-Azidobenzoic acid  
 14860-64-1, 4-Azidoaniline 18358-13-9, Methacrylate, reactions  
 25322-68-3, Polyethyleneglycol 38862-24-7 67881-98-5,  
 2-Methacryloyloxyethylphosphorylcholine 88539-10-0  
 150120-15-3 158760-93-1 163674-39-3 163801-79-4

179637-06-0 179953-15-2

RL: RCT (Reactant); RACT (Reactant or reagent)  
 (substance immobilization agent consisting of polymer with  
 photoreactive groups, and use in substance immobilization  
 method/substrate)

IT 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine


88539-10-0 150120-15-3 158760-93-1

179637-06-0 179953-15-2

RL: RCT (Reactant); RACT (Reactant or reagent)  
 (substance immobilization agent consisting of polymer with  
 photoreactive groups, and use in substance immobilization  
 method/substrate)


RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
 tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 88539-10-0 HCAPLUS

CN 3,5,12-Trioxa-4-phosphapentadec-14-en-1-aminium, 4-hydroxy-N,N,N,14-  
 tetramethyl-13-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 150120-15-3 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N-trimethyl-9-  
 oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 158760-93-1 HCAPLUS


CN 3,5,10-Trioxa-4-phosphatridec-12-en-1-aminium, 4-hydroxy-N,N,N,12-  
 tetramethyl-11-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RN 179637-06-0 HCAPLUS  
 CN 3,5,16-Trioxa-4-phosphanonadec-18-en-1-aminium, 4-hydroxy-N,N,N,18-tetramethyl-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RN 179953-15-2 HCAPLUS  
 CN 3,5-Dioxa-8-aza-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Ito	2003	24	3021	Biomaterials	
Motorola Inc	2001			EP 1190254 A	HCAPLUS
Motorola Inc	2001			WO 2001001143 A	HCAPLUS
Motorola Inc	2001			JP 2003524150 A	
Motorola Inc	2001			US 6372813 A	HCAPLUS
Surmodics Inc	2002			EP 1141385 A	HCAPLUS
Surmodics Inc	2002			WO 2000040593 A	HCAPLUS
Surmodics Inc	2002			JP 2002534663 A	
Surmodics Inc	2002			US 6465178 A	HCAPLUS
The Kanagawa Academy Of	2004			JP 2004125781 A	HCAPLUS

L82 ANSWER 11 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:1038468 HCAPLUS

DOCUMENT NUMBER: 142:22289

TITLE: Anti-Fc antibodies and fragments or antibody-binding molecules for clinical diagnosis and biotechnological research

INVENTOR(S): Kawamura, Kenji

PATENT ASSIGNEE(S): Sumitomo Bakelite Co., Ltd., Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 6 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2004340818	A2	20041202	JP 2003-139139	20030516
PRIORITY APPLN. INFO.:			JP 2003-139139	20030516

ED Entered STN: 03 Dec 2004

AB Disclosed is an immunoassay container comprising ≥1

antibody-binding mol. such as protein A, protein G, anti-Fc antibodies or Fab fragments. The antibody-binding mol.-coated container is treated with hydrophilic polymer or grafted hydrophilic polymer. The immunoassay container is useful for ELISA without the requirement of blocking procedure and is especially useful for clin. diagnosis and biotechnol. purpose.

IC ICM G01N033-543

CC 15-1 (Immunochemistry)

Section cross-reference(s): 9

IT Biotechnology

Diagnosis

Human

**Immunoassay**

(anti-Fc antibodies and fragments or antibody-binding mols. for clin. diagnosis and biotechnol. research)

IT **Immunoassay**

(apparatus; anti-Fc antibodies and fragments or antibody-binding mols. for clin. diagnosis and biotechnol. research)

IT **Immunoassay**

(enzyme-linked immunosorbent assay; anti-Fc antibodies and fragments or antibody-binding mols. for clin. diagnosis and biotechnol. research)

IT 25249-16-5, 2-Polyhydroxyethylmethacrylate **125275-25-4**,

2-Methacryloyloxyethylphosphorylcholine-butylmethacrylate copolymer

RL: ARU (Analytical role, unclassified); BSU (Biological study,

unclassified); DGN (Diagnostic use); ANST (Analytical study); BIOL

(Biological study); USES (Uses)

(anti-Fc antibodies and fragments or antibody-binding mols. for clin. diagnosis and biotechnol. research)

IT **125275-25-4**, 2-Methacryloyloxyethylphosphorylcholine-butylmethacrylate copolymer

RL: ARU (Analytical role, unclassified); BSU (Biological study,

unclassified); DGN (Diagnostic use); ANST (Analytical study); BIOL

(Biological study); USES (Uses)

(anti-Fc antibodies and fragments or antibody-binding mols. for clin. diagnosis and biotechnol. research)


RN 125275-25-4 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 97-88-1

CMF C8 H14 O2



L82 /ANSWER 12 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:348566 HCAPLUS

DOCUMENT NUMBER: 141:153361

TITLE: Evaluation of stabilizing effect for several monoclonal antibody immobilized quartz crystal microbalance by stabilizer reagents

AUTHOR(S): Park, Jong-Won; Kurosawa, Shigeru; Aizawa, Hidenobu; Naganawa, Ryuichi; Yamada, Satoshi; Ishihara, Kazuhiko  
CORPORATE SOURCE: National Institute of Advanced Industrial Science and Technology (AIST), Tsukuba, 305-8565, Japan

SOURCE: Proceedings of the IEEE International Frequency Control Symposium & PDA Exhibition jointly with the 17th European Frequency and Time Forum, Tampa, FL, United States, May 4-8, 2003/(2003), 978-980. Institute of Electrical and Electronics Engineers: New York, N. Y.

CODEN: 69FHUZ; ISBN: 0-7803-7688-9

DOCUMENT TYPE: Conference

LANGUAGE: English

ED Entered STN: 29 Apr 2004

AB We tested five stabilizers for remaining immunol. activity of anti-dinitrophenol (DNP), anti-C-reactive protein (CRP), and anti-2,3,7,8 tetrachlorodibenzo-p-dioxin (TCDD) antibody immobilized QCM under several storage conditions. Investigated nine stabilizers were as following; 0.25 % BlockAce and StabilGuard as com. available reagents, 0.2 % glycerin and 0.2 % Bovine Serum Albumin (BSA) as conventional stabilizers, 0.2 % MPC copolymer stabilizer as developed reagent, and PBS solution as blank reagent. According to the exptl. results, we found MPC copolymer (NOF Corp., Japan) coated QCM showed highly immunol. activity through specific antigens and their antibodies after the heat acceleration test and long-term storage.

CC 9-10 (Biochemical Methods)

IT **Immunoassay**

Stability

Stabilizing agents

(evaluation of stabilizing effect for several monoclonal antibody immobilized quartz crystal microbalance by stabilizer reagents)

IT 56-81-5, Glycerin, analysis 125275-25-4 658045-21-7, BlockAce 658051-44-6, StabilGuard

RL: ARU (Analytical role, unclassified); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study)

(evaluation of stabilizing effect for several monoclonal antibody immobilized quartz crystal microbalance by stabilizer reagents)

IT 125275-25-4

RL: ARU (Analytical role, unclassified); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study)


(evaluation of stabilizing effect for several monoclonal antibody immobilized quartz crystal microbalance by stabilizer reagents)

RN 125275-25-4 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)


CM 1

CRN 67881-98-5  
CMF C11 H22 N O6 P



CM 2

CRN 97-88-1  
CMF C8 H14 O2



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Bunde, R	1998	46	1223	Talanta	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Kurosawa, S	2000	47	1256	IEEE Trans Ultrason,	
Luppa, P	2001	314	1	Clin Chim Acta	HCAPLUS
Park, J	2003	50	193	IEEE Trans Ultrason,	
Sakaki, S	1999	47	523	J Biomed Mater Res	HCAPLUS
Sakaki, S	2000	32	637	Polym J	HCAPLUS
Steegborn, C	1997	12	19	Biosens Bioelect	HCAPLUS
Towery, R	2001	16	1	Biosen Bioelect	HCAPLUS

L82 ANSWER 13 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:338278 HCAPLUS

DOCUMENT NUMBER: 141:273834

TITLE: Reproducibility evaluation of immunosensor with antibody-immobilized beads column (II)

AUTHOR(S): Fuchiwaki, Yusuke; Rikitake, Kotaro; Futagami, Norimichi; Yasuzawa, Mikito

CORPORATE SOURCE: Department of Chemical Science and Technology, Faculty of Engineering, The University of Tokushima, Tokushima, 770-8506, Japan

SOURCE: Chemical Sensors (2003), 19(Suppl. B), 64-66  
CODEN: KAGSEU

PUBLISHER: Denki Kagakkai Kagaku Sensa Kenkyukai

DOCUMENT TYPE: Journal

LANGUAGE: Japanese

ED Entered STN: 26 Apr 2004

AB The reproducibility of the antibody-immobilized beads by the acidic antigen-antibody bonds cleavage treatment was performed using an optical procedure. The fluorescein derivative was conjugated to the antigen in order to determine the affinity of the immobilized antibody. The attempt to reduce the nonspecific adsorption of fluorescein-labeled antigen to immobilized antibody was also performed by adding polyethylene glycol and

poly(2-methacryloyloxyethyl phosphorylcholine) in the reaction solns. Although, both polymers were effective to reduce the nonspecific adsorption, either could eliminate the nonspecific adsorption. Fair reproducibility was observed for approx. first five repetition, while the affinity reduction of the immobilized antibody was inevitable by the acidic treatment of pH 3.

CC 9-10 (Biochemical Methods)

Section cross-reference(s): 2, 15

IT 25322-68-3, Polyethylene glycol 67881-99-6, Poly(2-methacryloyloxyethyl phosphorylcholine)

RL: ARU (Analytical role, unclassified); ANST (Analytical study) (effect on prevention of non-specific adsorption on beads of; reproducibility evaluation of **immunosensor** with antibody-immobilized beads column)

IT 67881-99-6, Poly(2-methacryloyloxyethyl phosphorylcholine)

RL: ARU (Analytical role, unclassified); ANST (Analytical study) (effect on prevention of non-specific adsorption on beads of; reproducibility evaluation of **immunosensor** with antibody-immobilized beads column)


RN 67881-99-6 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, homopolymer (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N 06 P



L82 ANSWER 14 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2002:172237 HCAPLUS

DOCUMENT NUMBER: 136:213193

TITLE: Highly reproducible agglutination immunoassay method and reagent

INVENTOR(S): Shigenobu, Kayoko; Shuto, Kenshiro; Sakaki, Shujiro

PATENT ASSIGNEE(S): Kyowa Medex Co., Ltd, Japan; Nof Corporation

SOURCE: PCT Int. Appl., 35 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent


LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002018953	A1	20020307	WO 2001-JP7385	20010828
W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG,				

US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM  
 RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY,  
 DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF,  
 BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG  
 CA 2420770 AA 20020307 CA 2001-2420770 20010828  
 AU 2001080210 A5 20020313 AU 2001-80210 20010828  
 EP 1314982 A1 20030528 EP 2001-958575 20010828  
 R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,  
 IE, SI, LT, LV, FI, RO, MK, CY, AL, TR  
 US 2003166302 A1 20030904 US 2003-363038 20030228  
 PRIORITY APPLN. INFO.: JP 2000-259964 A 20000829  
 WO 2001-JP7385 W 20010828  
 ED Entered STN: 08 Mar 2002  
 AB A highly reproducible agglutination immunoassay method is provided, in  
 which the agglutination of insol. carrier particles (e.g., latex) takes  
 place in a stable and homogeneous way. An immunoassay reagent used for  
 this method is also provided. In this agglutination immunoassay method,  
 an antigenic substance in a test sample is bound to the insol. carrier  
 particles substantially not carrying any bound-antigen or -antibody, and  
 then, an antibody or an antibody complex capable of specifically reacting  
 with the antigenic substance is bound to the particles to selectively give  
 rise to the agglutination. The reagent contains a polymer which is prepared  
 either by homogeneously polymerizing a monomer possessing a phosphorylcholine  
 group and a vinyl group (e.g., 2-methacryloyloxyethylphosphorylcholine), or  
 co-polymerizing the monomer possessing a phosphorylcholine group and a vinyl  
 group, and another monomer possessing a vinyl group (e.g.,  
 n-butylmethacrylate). An improved reproducibility was obtained when the  
 HbA1c concentration in blood samples were determined with this reagent using  
 anti-HbA1c monoclonal antibody in comparison to the conventional reagents.  
 IC ICM G01N033-543  
 CC 9-10 (Biochemical Methods)  
 IT **Immunoassay**  
 (agglutination test; highly reproducible agglutination immunoassay  
 method and reagent)  
 IT 97-88-1, n-Butylmethacrylate **67881-98-5**, 2-  
 Methacryloyloxyethylphosphorylcholine  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (highly reproducible **agglutination immunoassay**  
 method and reagent)  
 IT **67881-98-5**, 2-Methacryloyloxyethylphosphorylcholine  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (highly reproducible **agglutination immunoassay**  
 method and reagent)  
 RN 67881-98-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
 tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Hitachi Chemical Co Ltd	1987			JP 62218865 A	HCAPLUS



Nitto Electric Co Ltd	1986		JP 61274261 A	HCAPLUS
Oriental Yeast Co Ltd	1996		JP 08101196 A	HCAPLUS

L82 ANSWER 15 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2002:692742 HCAPLUS

DOCUMENT NUMBER: 138:316997

TITLE: Water-soluble phospholipid polymers as a novel synthetic blocking reagent in an immunoassay system

AUTHOR(S): Sakaki, Shujiro; Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Tokyo Medical and Dental University, Tokyo, Japan

SOURCE: Biomedical Diagnostic Science and Technology (2002), 353-366. Editor(s): Law, Wai Tak; Akmal, Naim; Usmani, Arthur M. Marcel Dekker, Inc.: New York, N. Y.

CODEN: 69DBEP; ISBN: 0-8247-0725-7

DOCUMENT TYPE: Conference

LANGUAGE: English

ED Entered STN: 13 Sep 2002

AB Water-soluble amphiphilic PMSt, a random copolymer of 2-methacryloyloxyethyl phosphorylcholine and styrene, was found to be an effective blocking reagent in the ELISA method. Measurements of the nonspecific adsorption and immunol. activities of enzyme-antibody conjugate showed that PMSt has comparable or even better performance than bovine serum albumin and casein as a blocking reagent. No characteristic changes of PMSt were observed from its being frozen and melted repeatedly.

CC 9-10 (Biochemical Methods)

IT **Immunoassay**

(enzyme-linked immunosorbent assay; water-soluble phospholipid polymers as novel synthetic blocking reagent in immunoassay system)

IT **Immunoassay**

(water-soluble phospholipid polymers as novel synthetic blocking reagent in immunoassay system)

IT **134483-35-5**

RL: ARU (Analytical role, unclassified); ANST (Analytical study)

(water-soluble phospholipid polymers as novel synthetic blocking reagent in **immunoassay** system)IT **134483-35-5**

RL: ARU (Analytical role, unclassified); ANST (Analytical study)

(water-soluble phospholipid polymers as novel synthetic blocking reagent in **immunoassay** system)


RN 134483-35-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI)  
(CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P



CM 2

CRN 100-42-5

CMF C8 H8

 $\text{H}_2\text{C}=\text{CH}-\text{Ph}$ 

## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Anderton, B	1980	2	122	Immunol Today	
Donovan, J	1975	250	1966	J Biol Chem	HCAPLUS
Farr, A	1981	47	129	J Immunol Method	HCAPLUS
Ghosh, S	1974	337	395	Biochem Biophys Acta	HCAPLUS
Ikemi, M	1982	15	281	Macromolecules	HCAPLUS
Imagawa, M	1982	4	41	J Appl Biochem	HCAPLUS
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Ishihara, K	1994	32	859	J Polym Sci A, Polym	HCAPLUS
Ishihara, K	1990	22	355	Polym J	HCAPLUS
Ishikawa, E	1987	1	238	J Clin Lab Anal	HCAPLUS
Ishikawa, E	1983	4	209	J Immunoassay	HCAPLUS
Isikawa, E	1983	18	219	Develop Immunol	
Orci, L	1985	28	528	Diabetologia	HCAPLUS
Osborn, M	1978	77	R27	J Cell Biol	MEDLINE
Wisdom, G	1976	22	1243	Clin Chem	HCAPLUS

L82 ANSWER 16 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2001:617197 HCAPLUS

DOCUMENT NUMBER: 135:192510

TITLE: Microparticle dispersion agent for clinical test,  
reagent for clinical test, its manufacturing method,  
clinical test method and application

INVENTOR(S): Shudo, Kenshiro; Sakaki, Shujiro; Yamada, Satoru;  
Sakamoto, Nobuyuki; Suzuki, Tadashi

PATENT ASSIGNEE(S): Nof Corporation, Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 12 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
-----	----	-----	-----	-----
JP 2001228149	A2	20010824	JP 2000-34931	20000214
PRIORITY APPLN. INFO.:			JP 2000-34931	20000214

ED Entered STN: 24 Aug 2001

AB A microparticle dispersion agent for a clin. test is provided, which improves a dispersion stability of the microparticle-containing reagent and a redispersion ability of the microparticles for clin. test coagulated during the process of reagent preparation or measurement without decreasing the activity of the bound antigen or antibody. The microparticle dispersion agent possessing high reproducibility and high sensitivity is processed by a simple method suited for an automated analyzer. The agent contains as an effective component a polymer prepared by polymerizing the monomer composition

consisting of phosphorylcholin-analogous group-containing monomer (e.g., 2-(meth)acryloyloxyethyl-2'-(trimethylammonio)ethylphosphate).

IC ICM G01N033-531

ICS G01N033-543

CC 9-10 (Biochemical Methods)

IT **Immunoassay**

(agglutination test; microparticle dispersion agent for clin. test, reagent for clin. test, manufacturing method, clin. test method and application)

IT **Immunoassay**

(apparatus, automated; microparticle dispersion agent for clin. test, reagent for clin. test, manufacturing method, clin. test method and application)

IT 97-88-1, n-Butylmethacrylate 18358-13-9, Methacrylate, reactions

**67881-98-5 150120-15-3**

RL: RCT (Reactant); RACT (Reactant or reagent)

(microparticle dispersion agent for clin. test, reagent for clin. test, manufacturing method, clin. test method and application)


IT **67881-98-5 150120-15-3**

RL: RCT (Reactant); RACT (Reactant or reagent)

(microparticle dispersion agent for clin. test, reagent for clin. test, manufacturing method, clin. test method and application)


RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RN 150120-15-3 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N-trimethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



L82 ANSWER 17 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2001:759015 HCAPLUS

DOCUMENT NUMBER: 137:68069

TITLE: A comparison of the use of an ATP-based bioluminescent assay and image analysis for the assessment of bacterial adhesion to standard HEMA and biomimetic soft contact lenses

AUTHOR(S): Andrews, C. S.; Denyer, S. P.; Hall, B.; Hanlon, G. W.; Lloyd, A. W.

CORPORATE SOURCE: School of Pharmacy and Biomolecular Sciences, Biomedical Materials Research Group, University of Brighton, Moulsecoomb, Brighton, BN2 4GJ, UK

SOURCE: Biomaterials (2001), 22(24), 3225-3233

PUBLISHER: CODEN: BIMADU; ISSN: 0142-9612  
 Elsevier Science Ltd.  
 DOCUMENT TYPE: Journal  
 LANGUAGE: English

ED Entered STN: 18 Oct 2001

AB The aim of this study was to investigate in vitro adhesion of clin. relevant bacteria to standard HEMA and novel biomimetic soft contact lenses (SCL) using bioluminescent ATP assay and image anal. Unworn SCL were incubated with *Pseudomonas aeruginosa*, *Staphylococcus epidermidis* or *Serratia marcescens* suspended in sterile phosphate buffered saline (PBS). The level of bacterial adhesion after 1, 2, 4, 6 and 18 h, was assessed using both image anal. and a bioluminescent ATP assay. Species differences in the overall level of adhesion to the different types of lens were observed using both measurement techniques. Generally bacterial adhesion was shown to peak at 4-6 h, then decline to a much lower level by 18 h. After 4 h, adhesion of all species of bacteria to the biomimetic SCL (omafilcon A) was found to be significantly lower than to the standard HEMA SCL (polymacon) ( $p < 0.05$ , Student's t-test,  $n=4$ ). Both these techniques demonstrated that novel biomimetic SCL materials exhibit significantly lower bacterial adhesion in vitro compared to standard HEMA SCL materials. SCL manufactured with these novel biomimetic materials may reduce the risk of infection.

CC 63-7 (Pharmaceuticals)

IT 25053-81-0, Polymacon **144056-32-6**, Omaficon A

RL: DEV (Device component use); THU (Therapeutic use); BIOL (Biological study); USES (Uses)

(comparison of use of ATP-based bioluminescent **assay** and image anal. for assessment of bacterial adhesion to standard HEMA and biomimetic soft contact lenses)

IT **144056-32-6**, Omaficon A

RL: DEV (Device component use); THU (Therapeutic use); BIOL (Biological study); USES (Uses)

(comparison of use of ATP-based bioluminescent **assay** and image anal. for assessment of bacterial adhesion to standard HEMA and biomimetic soft contact lenses)


RN 144056-32-6 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with 1,2-ethanediyl bis(2-methyl-2-propenoate) and 2-hydroxyethyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 868-77-9


CMF C6 H10 O3



CM 3

CRN 97-90-5

CMF C10 H14 O4



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Cowell, B	1998	84	950	J Appl Microbiol	HCAPLUS
Denyer, S	1989		189	ATP luminescence rap	HCAPLUS
Elder, M	1993	273	509	The practitioner	
Gopinathan, U	1997	82	653	J Appl Microbiol	MEDLINE
Gristina, A	1987	237	1588	Science	HCAPLUS
Hall, B	1989	10	219	Biomaterials	HCAPLUS
Hayward, J	1984	5	135	Biomaterials	HCAPLUS
Holden, B	1996	22	47	CLAO J	MEDLINE
Lan, J	1999	27	218	Aust NZ J Ophthalmol	MEDLINE
Liesegang, T	1997	16	125	Cornea	MEDLINE
Lloyd, A	2000	23	119	Contact Lens Anterio	
Lloyd, A	1997	38	884	Invest Ophthalmol Vi	
Lundin, A	1989		11	ATP luminescence rap	HCAPLUS
Schultz, C	1995	15	243	J Ind Microbiol	HCAPLUS
Schultz, C	2000	24	113	J Ind Microbiol Biot	HCAPLUS
Schultz, C	2000	25	17	J Ind Microbiol Biot	HCAPLUS
Slusher, M	1987	105	110	Arch Ophthalmol	MEDLINE
Stern, G	1990	9	S36	Cornea	
Taylor, R	1998	75	23	Optometry Vision Sci	MEDLINE
Thakur, A	1999	27	224	Aust NZ J Ophthalmol	MEDLINE
Williams, T	1997	25	S30	Aust NZ J Ophthalmol	
Williams, T	1998	75	266	Optometry Vision Sci	MEDLINE

L82 / ANSWER 18 OF 70 HCAPLUS. COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2000:457297 HCAPLUS

DOCUMENT NUMBER: 133:86476

TITLE: Immunoassay container free from non-specific adsorption

INVENTOR(S): Tanaka, Hayao

PATENT ASSIGNEE(S): Sumitomo Bakelite Co., Ltd., Japan

SOURCE: PCT Int. Appl., 21 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
------------	------	------	-----------------	------

```

-----
WO 2000039582      A1      20000706      WO 1999-JP5979      19991028
W: AU, CA, CN, JP, KR, NO, NZ, RU, SG, US
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL,
    PT, SE
CA 2356857          AA      20000706      CA 1999-2356857      19991028
AU 9963667          A1      20000731      AU 1999-63667        19991028
EP 1152242          A1      20011107      EP 1999-951118       19991028
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,
    IE, FI
JP 3681983          B2      20050810      JP 2000-591430       19991028
JP 2005099040       A2      20050414      JP 2004-341443       20041126
PRIORITY APPLN. INFO.:
JP 1998-367404      A      19981224
JP 1999-56253       A      19990303
JP 1999-212096      A      19990727
JP 2000-591430      A3     19991028
WO 1999-JP5979      W      19991028

```

ED Entered STN: 07 Jul 2000

AB An immunoassay container is designed so that the saturation adsorption of mols. used in the assay is smaller than  $1 \times 10^{-1}$  pmol/cm<sup>2</sup>. It is actually free from non-specific adsorption causative of reagent loss, sensitivity decrease, and precision decrease. The inner surface of the container is formed or coated with highly hydrophilic polymer or highly hydrophobic polymer, preferably coated with highly hydrophilic polymer, or more preferably coated with extremely highly hydrophilic polymer to prevent the adsorption of assay reagent, antigen or antibody. The extremely highly hydrophilic polymer is selected from polyhydroxyalkylmethacrylate, polyoxyalkylene (C2-C4) group-containing methacrylate polymer or copolymer, polyvinylpyrrolidone, phospholipid-polymer complex, or 2-methacryloyloxyethylphosphorylcholine polymer or copolymer.

IC ICM G01N033-53

CC 9-10 (Biochemical Methods)

Section cross-reference(s): 15

IT **Immunoassay**

(enzyme-linked immunosorbent assay; immunoassay container free from non-specific adsorption)

IT Blood serum

Containers

Hydroxyl group

**Immunoassay**

(immunoassay container free from non-specific adsorption)

IT 79-41-4D, Methacrylic acid, polyoxyalkylene (C2-C4) group-containing polymer, and copolymer 97-88-1D, Butyl methacrylate, copolymer with 2-methacryloyloxyethyl phosphorylcholine 9002-84-0, Polytetrafluoroethylene 25087-26-7D, Polymethacrylic acid, hydroxyalkyl derivs. 25249-16-5 **67881-98-5D**, 2-Methacryloyloxyethyl phosphorylcholine, polymer and copolymer with butylmethacrylate  
 RL: NUU (Other use, unclassified); USES (Uses)

(immunoassay container free from non-specific adsorption)

IT **67881-98-5D**, 2-Methacryloyloxyethyl phosphorylcholine, polymer and copolymer with butylmethacrylate

RL: NUU (Other use, unclassified); USES (Uses)

(immunoassay container free from non-specific adsorption)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)

RL: SPN (Synthetic preparation); THU (Therapeutic use); BIOL (Biological study); PREP (Preparation); USES (Uses)  
(polymeric solid support-immobilized antigen or antibody and its use)


IT 67881-98-5DP, polymers and copolymers 67882-00-2P  
125275-25-4P 134483-35-5P 148569-41-9P

RL: SPN (Synthetic preparation); THU (Therapeutic use); BIOL (Biological study); PREP (Preparation); USES (Uses)

(polymeric solid support-immobilized antigen or antibody and its use)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 67882-00-2 HCAPLUS

CN Ethanaminium, 2-[[hydroxy[2-[(2-methyl-1-oxo-2-propenyl)oxy]ethoxy]phosphinyl]oxy]-N,N,N-trimethyl-, inner salt, polymer with methyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 80-62-6

CMF C5 H8 O2




RN 125275-25-4 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 97-88-1

CMF C8 H14 O2




RN 134483-35-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI)  
(CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 100-42-5

CMF C8 H8



RN 148569-41-9 HCAPLUS


CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with 2-hydroxyethyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P






CM 2

CRN 868-77-9

CMF C6 H10 O3



L82, ANSWER 22 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1998:112708 HCAPLUS

DOCUMENT NUMBER: 128:215281

TITLE: Stabilization of proteins with phosphorylcholine group-containing polymers and stabilized compositions  
 INVENTOR(S): Sakaki, Shujiro; Sudo, Kenshiro; Yamada, Akira; Ando, Ryota; Matsuyama, Kazuo; Koinuma, Yasuhiro; Nakabayashi, Norio; Ishihara, Kazuhiko

PATENT ASSIGNEE(S): Nippon Oil and Fats Co., Ltd., Japan; Nakabayashi, Norio; Ishihara, Kazuhiko; Research Development Corp. of Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 11 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 10045794	A2	19980217	JP 1996-202367	19960731
PRIORITY APPLN. INFO.:			JP 1996-202367	19960731

ED Entered STN: 25 Feb 1998

AB Proteins, e.g. albumins, blood coagulation factors, Igs, enzymes for cleaning contact lenses, etc., are stabilized in the presence of phosphorylcholine group-containing polymers. Also claimed are stabilized comps. containing (A) polymers prepared from monomers containing  $\text{CH}_2:\text{CXCO}_2\text{CH}_2\text{CH}_2\text{OP}(\text{O})(\text{O}-)\text{OCH}_2\text{CH}_2\text{NMe}_3$  (X = H, Me) 1.0 + 10<sup>-4</sup>-80, (B) proteins for plasma preps. or labeled immunoreactive substances 1.0 + 10<sup>-14</sup>-20, and (C) buffer solns. 0-99.9 weight%. A solution of 5.0 + 10<sup>-4</sup> weight% human IgG and 2.0 weight% poly(2-methacryloyloxyethylphosphorylcholine) (preparation given) in a Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub> buffer was incubated at 40° for 4 wk to show reactivity with anti-human IgG1 antibody-immobilized plate 104.8%, vs. 5.5% for a control solution containing BSA instead of the polymer.

ICM C07K001-00

ICS C07K014-00; C07K016-00; C08L089-00; C08L101-00; C12N009-96;  
 C12Q001-25; A61K038-00; A61K038-16; A61K038-43; C08F030-02;  
 C08L043-02

CC 9-16 (Biochemical Methods)  
 Section cross-reference(s): 7, 15, 38, 63


IT **67881-99-6P 125275-25-4P 134483-35-5P**  
 RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use);  
 PNU (Preparation, unclassified); THU (Therapeutic use); ANST (Analytical  
 study); BIOL (Biological study); PREP (Preparation); USES (Uses)  
 (stabilization of proteins such as albumins and **Igs** and  
 enzymes and their labeled products with phosphorylcholine group-containing  
 polymers)

IT **67881-99-6P 125275-25-4P 134483-35-5P**  
 RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use);  
 PNU (Preparation, unclassified); THU (Therapeutic use); ANST (Analytical  
 study); BIOL (Biological study); PREP (Preparation); USES (Uses)  
 (stabilization of proteins such as albumins and **Igs** and  
 enzymes and their labeled products with phosphorylcholine group-containing  
 polymers)

RN 67881-99-6 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
 tetramethyl-9-oxo-, inner salt, 4-oxide, homopolymer (9CI) (CA INDEX  
 NAME)

CM 1


CRN 67881-98-5  
 CMF C11 H22 N O6 P



RN 125275-25-4 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-  
 tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl  
 2-methyl-2-propenoate (9CI) (CA INDEX NAME)


CM 1

CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2


CRN 97-88-1  
 CMF C8 H14 O2



RN 134483-35-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI)  
 (CA INDEX NAME)


CM 1

CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2

CRN 100-42-5  
 CMF C8 H8



L82 ANSWER 23 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN  
 ACCESSION NUMBER: 1995:606836 HCAPLUS  
 DOCUMENT NUMBER: 123:5146  
 TITLE: Protein adsorption-preventing polymers or copolymers  
 INVENTOR(S): Sakaki, Hidejiro; Nakada, Shinji; Matsumoto, Takeo;  
 Koinuma, Yasuyoshi; Nakabayashi, Norio; Ishihara,  
 Kazuhiko  
 PATENT ASSIGNEE(S): Nippon Oils & Fats Co Ltd, Japan  
 SOURCE: Jpn. Kokai Tokkyo Koho, 6 pp.  
 CODEN: JKXXAF  
 DOCUMENT TYPE: Patent  
 LANGUAGE: Japanese  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 07083923	A2	19950331	JP 1993-228973	19930914
JP 3443891	B2	20030908		

PRIORITY APPLN. INFO.: JP 1993-228973 19930914  
 ED Entered STN: 14 Jun 1995  
 AB 2-Methacryloyloxyethyl phosphorylcholine polymer and copolymer containing 2-methacryloyloxyethyl phosphorylcholine are used for preventing protein adsorption. The (co)polymers are useful for increasing the reproductivity and accuracy of two-site method, e.g. antigen or antibody sandwich

immunoassay, for biochem. or clin. diagnosis. In example, poly-2-methacryloyloxyethyl phosphorylcholine, and 2-methacryloyloxyethyl phosphorylcholine copolymd. with Bu methacrylate, Me methacrylate, 2-hydroxyethyl methacrylate, or styrene were prepared. The prepared polymer or copolymers were used for preventing adsorption of FITC-labeled mouse anti-human carcinoembryonic antigen IgG during immunoassay.

IC ICM G01N033-531

ICS G01N033-543

CC 9-15 (Biochemical Methods)

IT **Immunoassay**

(methacryloyloxyethyl phosphorylcholine polymer or copolymer for preventing protein adsorption in two-site anal. method or immunoassay)

IT **67881-98-5DP**, 2-Methacryloyloxyethyl phosphorylcholine, polymers or copolymers **67881-99-6P 67882-00-2P**

**125275-25-4P 134483-35-5P 148569-41-9P**

RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(methacryloyloxyethyl phosphorylcholine polymer or copolymer for preventing protein adsorption in two-site anal. method or **immunoassay**)

IT **67881-98-5DP**, 2-Methacryloyloxyethyl phosphorylcholine, polymers or copolymers **67881-99-6P 67882-00-2P**


**125275-25-4P 134483-35-5P 148569-41-9P**

RL: ARU (Analytical role, unclassified); MOA (Modifier or additive use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(methacryloyloxyethyl phosphorylcholine polymer or copolymer for preventing protein adsorption in two-site anal. method or **immunoassay**)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)




RN 67881-99-6 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, homopolymer (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P




RN 67882-00-2 HCAPLUS  
 CN Ethanaminium, 2-[[hydroxy[2-[(2-methyl-1-oxo-2-propenyl)oxy]ethoxy]phosphinyl]oxy]-N,N,N-trimethyl-, inner salt, polymer with methyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 80-62-6

CMF C5 H8 O2




RN 125275-25-4 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 97-88-1


CMF C8 H14 O2



RN 134483-35-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI)  
 (CA INDEX NAME)


CM 1

CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2


CRN 100-42-5  
 CMF C8 H8



RN 148569-41-9 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with 2-hydroxyethyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)


CM 1

CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2

CRN 868-77-9  
 CMF C6 H10 O3



=> d ibib ab hitstr 24-26

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

L82 ANSWER 24 OF 70 USPATFULL on STN

ACCESSION NUMBER: 2004:203385 USPATFULL  
 TITLE: Agglutination accelerator for immunological measurement  
 INVENTOR(S): Sumida, Kyoichi, Amagasaki-shi, JAPAN  
 Wada, Koji, Amagasaki-shi, JAPAN  
 Ishihara, Kazuhiko, Bunkyo-ku, JAPAN  
 PATENT ASSIGNEE(S): WAKO PURE CHEMICALS INDUSTRIES, LTD., Chuo-ku, JAPAN  
 (non-U.S. corporation)

	NUMBER	KIND	DATE
PATENT INFORMATION:	US 2004157276	A1	20040812
APPLICATION INFO.:	US 2004-626502	A1	20040304 (10)

	NUMBER	DATE
PRIORITY INFORMATION:	JP 2001-169051	20010605
DOCUMENT TYPE:	Utility	
FILE SEGMENT:	APPLICATION	
LEGAL REPRESENTATIVE:	MERCHANT & GOULD PC, P.O. BOX 2903, MINNEAPOLIS, MN, 55402-0903	
NUMBER OF CLAIMS:	16	
EXEMPLARY CLAIM:	1	
NUMBER OF DRAWINGS:	1 Drawing Page(s)	
LINE COUNT:	1035	

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

AB An object of the present invention is to provide an immunoassay of PSA using an agglutination accelerator, which has an agglutination accelerating effect equal to or stronger than the known agglutination accelerator; hardly generates non-specific turbidity; and hardly generates salting out even in a solution with a high salt concentration. The present invention relates to an immunoassay of a prostate-specific antigen comprising performing an antigen-antibody reaction in the presence of a polymer having a monomer unit derived from a monomer represented by the following general formula [2]: ##STR1##

(wherein R.sup.1-R.sup.3 are each independently a hydrogen atom or an alkyl group optionally having a hydroxyl group; R.sup.4 is an alkylene group; R.sup.5 is an alkylene group optionally having a substituent and optionally having an oxygen atom in a chain; R.sup.6 is a hydrogen atom or a methyl group, and X is an oxygen atom or a --NH-- group), and a kit of reagent for an immunoassay comprising a reagent containing an agglutination accelerator for the immunoassay.

IT 67881-99-6P, Poly-2-methacryloyloxyethylphosphorylcholine 125275-25-4P, n-Butyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer 144514-08-9P, Stearyl methacrylate-2-Methacryloyloxyethylphosphorylcholine copolymer 313216-64-7P 478015-82-6P (agglutination-promoting agent for antigen or antibody immunoassay)


RN 67881-99-6 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, homopolymer (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P




RN 125275-25-4 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with butyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 97-88-1

CMF C8 H14 O2




RN 144514-08-9 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with octadecyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5

CMF C11 H22 N O6 P






CM 2

CRN 32360-05-7

CMF C22 H42 O2




RN 313216-64-7 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with phenylmethyl 2-methyl-2-propenoate (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 2495-37-6

CMF C11 H12 O2




RN 478015-82-6 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with 2-hydroxy-N,N,N-trimethyl-3-[(2-methyl-1-oxo-2-propenyl)oxy]-1-propanaminium chloride (9CI) (CA INDEX NAME)

CM 1

CRN 67881-98-5


CMF C11 H22 N O6 P



CM 2

CRN 13052-11-4

CMF C10 H20 N O3 . Cl

● Cl<sup>-</sup>

L82 ANSWER 25 OF 70 USPATFULL on STN

ACCESSION NUMBER: 2003:238145 USPATFULL

TITLE: Highly reproducible agglutination immunoassay method and reagents

INVENTOR(S): Shigenobu, Kayoko, Sunto-gun, JAPAN  
Shuto, Kenshiro, Tsukuba-shi, JAPAN  
Sakaki, Shujiro, Tsukuba-shi, JAPAN

	NUMBER	KIND	DATE
PATENT INFORMATION:	US 2003166302	A1	20030904
APPLICATION INFO.:	US 2003-363038	A1	20030228 (10)
	WO 2001-JP7385		20010828


	NUMBER	DATE
PRIORITY INFORMATION:	JP 2000-259964	20000829
DOCUMENT TYPE:	Utility	
FILE SEGMENT:	APPLICATION	
LEGAL REPRESENTATIVE:	VENABLE, BAETJER, HOWARD AND CIVILETTI, LLP, P.O. BOX 34385, WASHINGTON, DC, 20043-9998	
NUMBER OF CLAIMS:	22	
EXEMPLARY CLAIM:	1	
LINE COUNT:	1030	

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

AB The present invention provides an agglutination immunoassay, wherein the agglutination of insoluble carrier particles such as latex are stabilized and uniformized to give good reproducibility, and a reagent therefor. In the agglutination immunoassay which comprises allowing an antigenic substance in a sample to bind to insoluble carrier particles carrying substantially neither antigens nor antibodies thereon, and allowing an antibody or an antibody complex which reacts specifically to the antigenic substance to bind to the antigenic substance to give a selective agglutination of the insoluble carrier particles, a

homopolymer prepared by polymerization of a monomer such as 2-methacryloyloxyethyl phosphorylcholine having a phosphorylcholine group and a vinyl group, or a copolymer prepared by polymerization of a monomer having a phosphorylcholine group and a vinyl group, with a monomer having a vinyl group such as n-butyl methacrylate is used.

IT 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine  
(highly reproducible agglutination immunoassay method and reagent)  
RN 67881-98-5 USPATFULL  
CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



L82 ANSWER 26 OF 70 USPATFULL on STN

ACCESSION NUMBER: 2001:170879 USPATFULL

TITLE: Measuring method and measuring reagent of C-reactive protein

INVENTOR(S): Yokohama, Hiroaki, Tokyo, Japan  
Umehara, Harumi, Shizuoka, Japan  
Matsumori, Shigeru, Shizuoka, Japan  
Yamada, Satoshi, Ibaraki, Japan  
Shuto, Kenshiro, Ibaraki, Japan  
Sakaki, Shujiro, Ibaraki, Japan  
Suzuki, Ken, Ibaraki, Japan

	NUMBER	KIND	DATE
PATENT INFORMATION:	US 2001026927	A1	20011004
APPLICATION INFO.:	US 2001-794323	A1	20010228 (9)

	NUMBER	DATE
PRIORITY INFORMATION:	JP 2000-54096	20000229
	JP 2000-54102	20000229
DOCUMENT TYPE:	Utility	
FILE SEGMENT:	APPLICATION	
LEGAL REPRESENTATIVE:	VENABLE, P.O. Box 34385, Washington, DC, 20043-9998	
NUMBER OF CLAIMS:	24	
EXEMPLARY CLAIM:	1	
NUMBER OF DRAWINGS:	5 Drawing Page(s)	
LINE COUNT:	1407	

CAS INDEXING IS AVAILABLE FOR THIS PATENT.

AB An object of the present invention is to provide a method and a reagent for measuring the subject substances containing high concentration of C-reactive protein without dilution while avoiding prozone phenomenon.


C-reactive protein is measured with a compound having a phosphrylcholine group and a cationic group shown by the general formula (I) [in the formula (I), where R<sup>sup.1</sup>, R<sup>sup.2</sup> and R<sup>sup.3</sup> stand for a hydrogen atom, substituted or non-substituted alkyl, or substituted or non-substituted alkenyl, and X<sup>sup.-</sup> stands for an inorganic anion or an organic anion) and an antibody to C-reactive protein. Or, C-reactive protein is measured with a surface active agent having a

phosphorylcholine group, a surface active agent having a cationic group shown by the formula (II) [Y.sub.1 stands for a hydrophobic group, and R.sub.1, R.sub.2 and R.sub.3 stand for a hydrogen atom, substituted or non-substituted alkyl, or substituted or non-substituted alkenyl], and an antibody to C-reactive protein. As an antibody to C-reactive protein, an antibody carried by a water-insoluble carrier such as latex made from polystyrene is preferable. ##STR1##

IT 67881-98-5, 2-Methacryloyloxyethylphosphorylcholine  
(antibody and surface active agent having phosphorylcholine group and surface active agent having cationic group for measuring C-reactive protein)

RN 67881-98-5 USPATFULL

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



=> d iall abeq tech abex 27

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

L82 ANSWER 27 OF 70 WPIX COPYRIGHT 2005 THE THOMSON CORP on STN  
ACCESSION NUMBER: 2003-175043 [17] WPIX  
DOC. NO. CPI: C2003-045673  
TITLE: Non-specific hybridization inhibitors based on polymer containing carboxyl, sulfone and phosphorylcholine groups, applicable in e.g. gene analysis and disease diagnosis in clinical examination.  
DERWENT CLASS: A89 B04 D16  
INVENTOR(S): ASHIHARA, Y; ISOMURA, M; SAKAKI, S; SHUTO, K; TSUCHIDA, M  
PATENT ASSIGNEE(S): (NIOF) NOF CORP; (NIOF) NIPPON OILS & FATS CO LTD; (ASHI-I) ASHIHARA Y; (ISOM-I) ISOMURA M; (SAKA-I) SAKAKI S; (SHUT-I) SHUTO K; (TSUC-I) TSUCHIDA M  
COUNTRY COUNT: 23  
PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG	MAIN	IPC
WO 2002088389	A1	20021107	(200317)*	JA	26	C12Q001-68	
RW: AT BE CH CY DE DK ES FI FR GB GR IE IT LU MC NL PT SE TR							
W: CN US							
JP 2003014767	A	20030115	(200317)		11	G01N033-566	
EP 1391519	A1	20040225	(200415)	EN		C12Q001-68	
R: AT BE CH CY DE DK ES FI FR GB GR IE IT LI LU MC NL PT SE TR							
US 2004219546	A1	20041104	(200473)			C12Q001-68	
CN 1547615	A	20041117	(200516)			C12Q001-68	

#### APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
-----------	------	-------------	------

WO 2002088389	A1	WO 2002-JP4128	20020425
JP 2003014767	A	JP 2002-123774	20020425
EP 1391519	A1	EP 2002-722775	20020425
		WO 2002-JP4128	20020425
US 2004219546	A1	WO 2002-JP4128	20020425
		US 2004-476069	20040211
CN 1547615	A	CN 2002-812948	20020425

## FILING DETAILS:

PATENT NO	KIND	PATENT NO
EP 1391519	A1 Based on	WO 2002088389

PRIORITY APPLN. INFO: JP 2001-128699      20010426  
 INT. PATENT CLASSIF.:

MAIN: C12Q001-68; G01N033-566  
 SECONDARY: C08F230-02; G01N033-53  
 ADDITIONAL: C08F030-02  
 INDEX: C08F030:02

## BASIC ABSTRACT:

WO 200288389 A UPAB: 20030312

NOVELTY - A non-specific inhibitor containing a polymer (H) having at least 1 carboxyl or sulfone group and a phosphorylcholine-like group in its molecule, which has a weight-average molecular weight of 1000-5000 and exhibits an effect of inhibiting non-specific hybridization, is new.

DETAILED DESCRIPTION - INDEPENDENT CLAIMS are also included for the following:

(1) reagents for chemical examination containing the inhibitors and examination reagents; and

(2) clinical examination by contacting a specimen with an examination reagent hybridizable with a definite nucleic acid-relating substance in the presence of the inhibitor, and detection of the reaction product.

USE - The inhibitors are applicable in e.g. gene analysis and disease diagnosis in clinical examination.

ADVANTAGE - With the inhibitors for clinical examination, non-specific hybridization can be inhibited and a nucleic acids-relating substance to be assayed can be easily and highly accurately detected.

Dwg.0/0

FILE SEGMENT: CPI

FIELD AVAILABILITY: AB; GI; DCN

MANUAL CODES: CPI: A04-A; A12-L04; A12-V03C2; B04-C03; B04-E01;  
 B11-C08F2; B12-K04A; B12-K04F; D05-H09; D05-H12D1

TECH UPTX: 20030312

TECHNOLOGY FOCUS - POLYMERS - Preferred Inhibitors: The phosphorylcholine-like group is a component of a phosphorylcholine-like group-containing monomer of formula (I).

X = a divalent group;

Y = 1-6C alkyleneoxy group;

Z1 = H or R5O(C=O);

R1 = H or methyl;

R2, R3, R4 = independently H, 1-6C alkyl or hydroxyalkyl;

m = 0 or 1;

n = 1-4; and

R5 = 1-10C alkyl or 1-10C hydroxalkyl.

The polymer (H) is particularly made by polymerizing 5-95 mol.% a phosphorylcholine-like group-containing monomer, carboxyl group-containing monomer or/and sulfone group-containing monomer, and optionally 5-60 mol.% a hydrophobic monomer. The hydrophobic monomer is of formula (II)

R6 = H or methyl;  
 L1 = -C6H4-, -C6H10-, -(C=O)O-, -O-, -(C=O)NH-, -O(C=O)- or -O(C=O)O-;  
 L2 = H, -(CH2)g-L3 or ((CH2)p-O)h-L3;  
 g, h = 1-24;  
 p = 3-5; and  
 L3 = H, methyl, -C6H5- or -OC6H5.

TECHNOLOGY FOCUS - INSTRUMENTATION AND TESTING - Preferred Examination  
 Method: The amount of inhibitor applied in the hybridization system is  
 equivalent to the presence of 0.0001-20 %, by weight polymer (H) in the  
 inhibitor.

ABEX

UPTX: 20030312

EXAMPLE - 2-(Methacryloyloxy)ethylphosphorylcholine (MPC; 19.4 g) and  
 methacrylic acid (0.6 g) in water (40 g) were polymerized with succinic  
 peroxide (1.6 g) at 60 degrees C for 8 hours, under nitrogen, to give a  
 polymer powder (14.6 g; weight-average molecular weight = 153000) after  
 precipitation from acetone; characterization by GPC and nuclear magnetic  
 resonance (NMR). The polymer was tested as hybridization inhibitor, e.g.  
 at 3.2 %, by weight, in the detection of a nucleic acid after  
 hybridization.

=&gt; d ibib ed ab hitind 28

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE,  
 BIOSIS' - CONTINUE? (Y)/N:y

L82 ANSWER 28 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 1  
 ACCESSION NUMBER: 2004:220235 TOXCENTER  
 COPYRIGHT: Copyright 2005 ACS  
 DOCUMENT NUMBER: CA14117282686G  
 TITLE: Biological evaluation and drug delivery application of  
 cationically modified phospholipid polymers  
 AUTHOR(S): Palmer, Rosemary R.; Lewis, Andrew L.; Kirkwood, Laura C.;  
 Rose, Susanna F.; Lloyd, Andrew W.; Vick, Terry A.;  
 Stratford, Peter W.  
 CORPORATE SOURCE: Drug Delivery, Biocompatibles UK Ltd., Surrey, GU9 8QL,  
 UK.  
 SOURCE: Biomaterials, (2004) Vol. 25, No. 19, pp. 4785-4796.  
 CODEN: BIMADU. ISSN: 0142-9612.  
 COUNTRY: UNITED KINGDOM  
 DOCUMENT TYPE: Journal  
 FILE SEGMENT: CAPLUS  
 OTHER SOURCE: CAPLUS 2004:357432  
 LANGUAGE: English  
 ENTRY DATE: Entered STN: 20041005  
 Last Updated on STN: 20041229

ED Entered STN: 20041005

Last Updated on STN: 20041229

AB Phospholipid-like polymers based on 2-methacryloyloxyethyl  
 phosphorylcholine containing varying amts. of the cationically charged monomer  
 choline methacrylate (CMA) from 0 to 30% have been prepared Substrates  
 coated with these materials were shown to bind significantly lower amts.  
 of specific proteins compared to the uncoated control. **ELISA**  
**assays** demonstrated that fibrinogen did not bind appreciably to  
 coatings containing 0-30% CMA, whereas albumin binding was seen to increase  
 significantly as the CMA content of the coating increased. Platelet  
 activation **assays**, measurement of plasma **coagulation**  
 time and whole blood contact scanning electron microg. demonstrated that

the haemocompatibility of the coatings was shown to be unaffected by the CMA component. The CMA polymer coatings have been shown to absorb/adsorb many different drug compds. covering a wide range of mol. wts. and release these in a controlled fashion. The range of cationic polymers assessed can interact with the net neg. charge found in many large therapeutic biomols., such as DNA fragments used in gene therapy, that may be of interest in the preventative treatment of conditions such as restenosis. Coronary stents coated with 6% or 20% CMA-containing polymers have been shown to load and release this type of genetic material irresp. of mol. weight of the biomol. Ex vivo and in vivo studies have shown that these compds. can be delivered to the stented section of the vessel with very low quantities delivered outside the vessel target area.

CC 63-7  
ST Miscellaneous Descriptors  
phospholipid choline methacrylate coating protein adsorption stent drug  
delivery  
RN 208035-88-5

=> d hitstr 28

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

'HITSTR' IS NOT A VALID FORMAT

In a multiframe environment, a format can only be used if it is valid in at least one of the files. Refer to file specific help messages or the STNGUIDE file for information on formats available in individual files.

REENTER DISPLAY FORMAT FOR ALL FILES (FILEDEFAULT):end

=> d ibib ed ab hitind 29-

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, USPATFULL, WPIX, TOXCENTER, MEDLINE, BIOSIS' - CONTINUE? (Y)/N:y

YOU HAVE REQUESTED DATA FROM 42 ANSWERS - CONTINUE? Y/(N):y

L82 ANSWER 29 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 5  
ACCESSION NUMBER: 2004:228146 TOXCENTER  
COPYRIGHT: Copyright 2005 ACS  
DOCUMENT NUMBER: CA14219360387X  
TITLE: Copolymers of 2-methacryloyloxyethyl phosphorylcholine (MPC) as biomaterials  
AUTHOR(S): Haris, Parvez I.; Nakabayashi, N.; Iwasaki, Y.  
CORPORATE SOURCE: Inst. Biomater. Bioeng., Tokyo Med. Dent. Univ., Kanda, Tokyo, 101-0062, Japan.  
SOURCE: Bio-Medical Materials and Engineering, (2004) Vol, 14, No. 4, pp. 345-354.  
CODEN: BMENEO. ISSN: 0959-2989.  
COUNTRY: JAPAN  
DOCUMENT TYPE: Journal  
FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 2004:810369  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20041014  
Last Updated on STN: 20050503  
ED Entered STN: 20041014  
Last Updated on STN: 20050503  
AB A review and discussion. Copolymers of 2-methacryloyloxyethyl

phosphorylcholine (MPC) showed good hemocompatibility as hypothesized.  
The hypothesis was surfaces having phosphorylcholine groups by polymerization  
of

MPC could accumulate phospholipids from blood stream and show good blood compatibility. The authors designed and prepared a methacrylate having a phosphorylcholine group. While it was possible to introduce them by polymer reactions, polymer reaction is not always good method to prepare the desired pure surface. This must be very important point to consider for biomaterials. The hypothesis was confirmed by changing copolymer composition. The adsorption amount of phospholipids on the surfaces increased with increasing the MPC units in the copolymers. On the other hand, increasing MPC units in MPC copolymers decreased adsorption amount of peptides. There is limitation in blood compatibility tests in vitro due to unstable characteristics of blood itself. The authors evaluated them with series of blood compatibility tests, in vitro, ex vivo and in vivo, on coated PMMA beads, modified hollow fibers for hemodialysis and 2 mm small diameter blood vessels, resp. These data suggested MPC is a promising methacrylate to develop good blood contacting devices, which may not require systemic **anticoagulation**. Conventional blood compatible biomaterials were not suitable to make permeable membranes. But MPC is soluble in water and we could prepare permeable membranes to various solutes by the copolymer. Introduction of MPC copolymers on cellulose and polysulfone hollow fiber membranes gave them nonthrombogenicity but it did not give adverse effect on their permeability. These data suggested that it is possible to apply them to hemodialyzers, oxygenators and percutaneous glucose sensors to keep diabetic patients easier. MPC surfaces are good hydrogel to minimize damage on tissues by lubricating between organs and the coated devices. They do not induce denaturation of peptides, which is beneficial to keep activities of enzymes longer. And poly-MPC dissolved is applicable to stabilize several bioactive peptides in aqueous phase. So MPC polymers are useful to minimize fouling by inhibiting the adsorption of bioactive proteins. MPC has high potential to develop many varieties of new biomaterials useful in so-called biotechnol. MPC and their copolymers are com. available from NOF (Tokyo, Japan) and Biocompatibles (UK, as PC technol.).

CC 63-0

ST Miscellaneous Descriptors

review methacryloyloxyethyl phosphorylcholine copolymer biomaterial  
RN 67881-98-5 (2-Methacryloyloxyethyl phosphorylcholine)

L82 ANSWER 30 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 7

ACCESSION NUMBER: 2003:125101 TOXCENTER

DOCUMENT NUMBER: PubMed ID: 12580778

TITLE: In vivo evaluation of a MPC polymer coated continuous flow left ventricular assist system

AUTHOR(S): Kihara Shin'ichiro; Yamazaki Kenji; Litwak Kenneth N; Litwak Philip; Kameneva Marina V; Ushiyama Hiroyuki; Tokuno Toshimasa; Borzelleca David C; Umezu Mitsuo; Tomioka Jun; Tagusari Osamu; Akimoto Takehide; Koyanagi Hitoshi; Kurosawa Hiromi; Kormos Robert L; Griffith Bartley P

CORPORATE SOURCE: Department of Surgery, McGowan Institute for Regenerative Medicine, University of Pittsburgh, 300 Technology Drive, Pittsburgh, PA 15219, U.S.A. kiharas@msx.upmc.edu

SOURCE: Artificial organs, (2003 Feb) 27 (2) 188-92.

Journal Code: 7802778. ISSN: 0160-564X.

COUNTRY: United States

DOCUMENT TYPE: (EVALUATION STUDIES)

Journal; Article; (JOURNAL ARTICLE)

FILE SEGMENT: MEDLINE



OTHER SOURCE: MEDLINE 2003106561  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20030527  
Last Updated on STN: 20030527

ED Entered STN: 20030527

Last Updated on STN: 20030527

AB The aim of this study was the evaluation of the thrombogenicity and the biocompatibility of the SunMedical EVAHEART left ventricular assist system (LVAS) coated with 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer compared to a diamond-like carbon (DLC) coating. Four calves were implanted with the MPC polymer-coated LVAS. Eight calves were implanted with DLC coated LVAS. The thrombogenicity and biocompatibility of the pumps were evaluated. At explant, 60.0 +/- 37.2% (5-85%) of the pump surface area was still coated with MPC polymer after the duration of 45.0 +/- 32.0 days. In 1 out of 4 MPC and 2 out of 8 DLC coated pumps, there was a very small amount of thrombus around the seal ring; otherwise the blood contacting surfaces were free of thrombus. Major organs were normal except for a few lesions in kidneys from both groups. The MPC polymer coated EVAHEART LVAS seems to have low thrombogenicity and high biocompatibility similar to the DLC coated system. The current study demonstrated that the MPC polymer coating shows great promise for being used as an antithrombogenic substrate for the LVAS due to its ease of application, significant cost benefit, and reduction in **anticoagulation** therapy in acute postoperative period.

CT Animals

Blood Physiology

Carbon

Cattle

\*Coated Materials, Biocompatible

Coated Materials, Biocompatible: AE, adverse effects

\*Heart-Assist Devices

Heart-Assist Devices: AE, adverse effects

\*Methacrylates

Methacrylates: AE, adverse effects

\*Phosphorylcholine

Phosphorylcholine: AE, adverse effects

\*Phosphorylcholine: AA, analogs & derivatives

Research Support, Non-U.S. Gov't

Thrombosis: ET, etiology

RN 107-73-3 (Phosphorylcholine)

67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)

7440-44-0 (Carbon)

CN 0 (Coated Materials, Biocompatible); 0 (Methacrylates)

L82 ANSWER 31 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 9

ACCESSION NUMBER: 2002:147879 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13713181798G

TITLE: Improving the blood compatibility of ion-selective electrodes by employing poly(mpc-co-bma), a copolymer containing phosphorylcholine, as a membrane coating

AUTHOR(S): Berrocal, Maria J.; Johnson, R. Daniel; Badr, Ibrahim H. A.; Liu, Mingdong; Gao, Dayong; Bachas, Leonidas G.

CORPORATE SOURCE: Department of Chemistry and Center of Membrane Sciences, University of Kentucky, Lexington, KY, 40506-0055, USA.

SOURCE: Analytical Chemistry, (2002) Vol. 74, No. 15, pp. 3644-3648.

CODEN: ANCHAM. ISSN: 0003-2700.

COUNTRY: UNITED STATES

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 2002:473221  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20020702  
Last Updated on STN: 20020924

ED Entered STN: 20020702

Last Updated on STN: 20020924

AB The hydrogel poly(2-methacryloyloxyethylphosphorylcholine-co-Bu methacrylate), or poly(MPC-co-BMA), was used as a coating for polyurethane- and poly(vinyl chloride)-based membranes to develop ion-selective electrodes (ISEs) with enhanced blood compatibility. Adverse interactions of poly(MPC-co-BMA) with blood were diminished due to the phosphorylcholine functionalities of the hydrogel, which mimic the phospholipid polar groups present on the surface of many cell membranes. As demonstrated by **immunostaining**, hydrogel-coated PVC membranes soaked in platelet-rich plasma showed less adhesion and activation of platelets than uncoated PVC membranes, indicating an improvement in biocompatibility owing to the hydrogel. Furthermore, little differences in the potentiometric response characteristics, e.g., slope, detection limit, and selectivity, of ISEs employing uncoated and coated membranes were observed

CC 9-7

ST Miscellaneous Descriptors  
blood compatibility ion selective electrode

RN 7439-93-2 (Lithium)  
7439-95-4 (Magnesium)  
7440-09-7 (Potassium)  
7440-23-5 (Sodium)  
7440-70-2 (Calcium)  
14798-03-9 (Ammonium)  
107-73-3 (Phosphoryl choline)  
9002-86-2 (Polyvinyl chloride)

RN 125275-25-4

L82 ANSWER 32 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 11

ACCESSION NUMBER: 2003:84457 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13818276154D

TITLE: Importance of a biofouling-resistant phospholipid polymer to create a heparinized blood-compatible surface

AUTHOR(S): Iwasaki, Yasuhiko; Shibata, Naoya; Ninomiya, Madoka; Kurita, Kimio; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo Medical and Dental University, Tokyo, 101-0062, Japan.

SOURCE: Journal of Biomaterials Science, Polymer Edition, (2002) Vol. 13, No. 3, pp. 323-335.  
CODEN: JBSEEA. ISSN: 0920-5063.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2002:522306

LANGUAGE: English

ENTRY DATE: Entered STN: 20030415  
Last Updated on STN: 20030429

ED Entered STN: 20030415

Last Updated on STN: 20030429

AB Heparinization is believed to be one of the methods to suppress thrombus formation on blood-contacting surfaces. However, this study hypothesizes that heparinization alone might not be sufficient to provide a blood-compatible surface; i.e., a surface property that resists biofouling

is necessary to obtain an effective heparin-modified surface. 2-Methacryloyloxyethyl phosphorylcholine (MPC) polymers with 2-aminoethyl methacrylate (AEMA) were synthesized to immobilize heparin through ionic bonding. The primary amino groups of AEMA were considered to be the polymer surface because the  $\xi$ -potential of the surface was pos. when the mole fraction of the AEMA units was above 0.2. The antithrombogenic character of the polymer surface modified with heparin was evaluated by both Lee-White and microsphere column methods. The **coagulation** period of human whole blood in the absence of **anticoagulant** in glass tubing coated with the MPC polymer was longer than that in the original glass tube. Cell adhesion was completely inhibited on the MPC polymer surface after contact with human whole blood without **anticoagulant**. However, many adherent blood cells were observed on poly(2-ethylhexyl methacrylate-co-AEMA) (no MPC unit) even after heparinization. These results strongly indicate that the MPC polymer is a useful substrate where the heparin works well and that the heparin-immobilized MPC polymer has superior blood compatibility to the simple MPC polymer.

CC

63-7

ST Miscellaneous Descriptors

phospholipid polymer heparin coated antithrombogenic

RN 25719-51-1 (2-Ethylhexyl methacrylate polymer)

9005-49-6 (Heparin)

RN 182816-96-2; 503182-41-0; 503182-42-1

L82 ANSWER 33 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 12

ACCESSION NUMBER: 2000:225730 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13414198008K

TITLE: Crosslinkable coatings from phosphorylcholine-based polymers

AUTHOR(S): Lewis, A. L.; Cumming, Z. L.; Goreish, H. H.; Kirkwood, L. C.; Tolhurst, L. A.; Stratford, P. W.

CORPORATE SOURCE: Research and Development Group, Farnham Business Park, Biocompatibles Ltd., Farnham, Surrey, GU9 8QL, UK.

SOURCE: Biomaterials, (2001) Vol. 22, No. 2, pp. 99-111.

CODEN: BIMADU. ISSN: 0142-9612.

COUNTRY: UNITED KINGDOM

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2000:895492

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020313

ED Entered STN: 20011116

Last Updated on STN: 20020313

AB 2-Methacryloyloxyethyl phosphorylcholine (MPC) was synthesized and then used in the preparation of crosslinked polymer membranes with lauryl methacrylate, hydroxypropyl methacrylate and trimethoxysilylpropyl methacrylate (crosslinker) comonomers. Some phys. aspects of the membrane properties were evaluated in order to establish the basis for the synthesis of a series of post-crosslinkable polymers. These materials were made by copolymer of the constituent monomers via a free radical method, and characterized using NMR, FT-IR, viscometry and elemental anal. The optimum crosslink d. and conditions required for curing coatings of these polymers were investigated using atomic force microscopy (AFM) and showed the inclusion of 5 mol% silyl crosslinking agent to be ideal. A nanoindentation technique was employed to determine if the coating developed elasticity upon crosslinking. The biol. properties of the coatings were evaluated using a variety of protein adsorption **assays** and blood

contacting expts., and an enzyme **immunoassay** was developed to detect *E. coli* in order to assess the level of bacterial adhesion to these biomaterials. Polymers of this type were shown to be very useful as coating materials for improving the biocompatibility of, or reducing the levels of adherent bacteria to medical devices.

CC 63-7

ST Miscellaneous Descriptors

phosphorylcholine methacrylate copolymer crosslinked coating;  
biocompatibility phosphorylcholine methacrylate copolymer

RN **144514-07-8** (3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium,  
4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with  
dodecyl 2-methyl-2-propenoate)

**67881-98-5** (2-Methacryloyloxyethyl phosphorylcholine)

RN **210570-82-4**

L82 ANSWER 34 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 13

ACCESSION NUMBER: 1999:204966 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13205054803B

TITLE: The effect of the chemical structure of the phospholipid  
polymer on fibronectin adsorption and fibroblast adhesion  
on the gradient phospholipid surface

AUTHOR(S): Iwasaki, Yasuhiko; Sawada, Shin-Ichi; Nakabayashi, Nobuo;  
Khang, Gilson; Lee, Hai Bang; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo  
Medical and Dental University, Tokyo, 101-0062, Japan.

SOURCE: Biomaterials, (1999) Vol. 20, No. 22, pp. 2185-2191.

CODEN: BIMADU. ISSN: 0142-9612.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1999:714631

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020403

ED Entered STN: 20011116

Last Updated on STN: 20020403

AB The interaction between biocomponents and the polyethylene (PE) surface modified with poly[ $\omega$ -methacryloyloxyalkyl phosphorylcholine (MAPC)] was considered taking into account the surface characteristics, i.e., d., mobility, and orientation of the poly(MAPC). The PE surface, grafted gradually with the poly(MAPC) was prepared by corona irradiation method. The amount of peroxide produced on the PE surface which was determined with 1,1-diphenyl-2-picrylhydrazyl, increased with an increase in the energy of the corona. The surface d. of the poly(MAPC) was increased with an increase in the amount of the peroxides produced by the corona irradiation. The orientation and mobility of the poly(MAPC) grafted on the PE surface was evaluated with 1,6-diphenyl-1,3,5-hexatriene. The orientation of the poly[6-methacryloyloxyhexyl phosphorylcholine (MHPC)] which has six methylene chains between the phospholipid polar group and the backbone was higher than that of other poly(MAPC)s. The mobility of the poly(MAPC) decreased with an increase in the methylene chain length in the MAPC unit. The fibronectin adsorption on the gradient PE sheet grafted with poly(MAPC) was determined with enzyme-labeled **immunoassay**. The amount of adsorbed fibronectin on the PE grafted with poly[2-methacryloyloxyethyl phosphorylcholine(MPC)] and poly(MHPC) decreased with an increase in their surface d. Especially, the PE sheet grafted with the poly(MHPC) was

effectively

reduced compared with other poly(MAPC)s. On the poly[10-methacryloyloxydecyl (MDPC)], there is a min. amount of adsorbed

fibronectin. The fibronectin adsorption pattern on the PE sheet grafted with poly(MAPC) was quite different from the chemical structure of the MAPC unit. The human normal diploid fibroblasts (WI-38 cells) were cultured on the gradient PE sheet grafted with poly(MAPC) changing the concentration of seeded WI-38 cells. The adhesion behavior of the WI-38 cells was different depending on the concentration of the seeded WI-38 cells. When the concentration was low, the number of the adherent WI-38 cells had the same

tendency

as fibronectin adsorption. The gradient PE sheet grafted with the poly(MHPC) effectively reduced WI-38 cells adhesion even when the concentration of the WI-38 cells was high. The biocompatibility of polymer surfaces can be improved by highly oriented phosphorylcholine group.

CC 63-7

ST Miscellaneous Descriptors

polyethylene graft phosphorylcholine methacrylate; fibroblast adhesion  
polyethylene graft phosphorylcholine methacrylate; fibronectin  
absorption polyethylene graft phosphorylcholine methacrylate

RN 252877-49-9; 252877-50-2; 252877-51-3

L82 ANSWER 35 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 14

ACCESSION NUMBER: 1999:214849 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13208098101A

TITLE: Inhibition of fibroblast cell adhesion on substrate by  
coating with 2-methacryloyloxyethylphosphorylcholine  
polymers

AUTHOR(S): Ishihara, Kazuhiko; Ishikawa, Eri; Iwasaki, Yasuhiko;  
Nakabayashi, Nobuo

CORPORATE SOURCE: Department of Materials Science, Graduate School of  
Engineering, The University of Tokyo, Tokyo, 113-8656,  
Japan.

SOURCE: Journal of Biomaterials Science, Polymer Edition, (1999)  
Vol. 10, No. 10, pp. 1047-1061.  
CODEN: JBSEEA. ISSN: 0920-5063.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1999:784811

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020403

ED Entered STN: 20011116

Last Updated on STN: 20020403

AB Fibroblast adhesion and growth behavior were examined on various polymers coated on a poly(ethylene terephthalate) (PET) substrate. The polymers used were poly[2-methacryloyloxyethylphosphorylcholine (MPC)-co-Bu methacrylate] copolymer (PMB)s with different MPC unit compns., and poly(2-hydroxyethyl methacrylate). Surface anal. by dynamic contact angle measurement revealed that the mobility of the polymer chain on the PET substrate depended on the MPC unit composition, but there was no significant difference between the PMBs with 3-10 mol% MPC units and poly(HEMA). Fibronectin adsorption on the polymer surface from a cell culture medium was determined by **immunoassay**. The adsorbed fibronectin was evenly distributed in every polymer, however, the amount was reduced with an increase in the MPC unit composition in the PMB. This result suggested that the MPC unit could weaken the interaction between the polymer surface and proteins. When fibroblast L-929 cells, were cultured on the polymers, the cells adhered and the number of cells increased on not only the hydrophobic poly(BMA) but also on the hydrophilic poly(HEMA). However, the number of cells that adhered on the PMB surface decreased with an increase in the

MPC unit composition This was a result of the fibronectin adsorption behavior. Thus, it could be concluded that since the PMB could suppress cell adhesion proteins e.g. fibronectin, the PMB showed excellent cell adhesive resistance properties.

CC 63-7

ST Miscellaneous Descriptors

fibroblast adhesion inhibition methacryloyloxyethylphosphorylcholine  
polymer coating

RN 125275-25-4 (Butyl methacrylate-2-methacryloyloxyethylphosphoryl  
choline copolymer)  
25038-59-9 (PET)  
25249-16-5 (Poly(2-hydroxyethyl methacrylate))

L82 ANSWER 36 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 16

ACCESSION NUMBER: 1998:191577 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13006071466V

TITLE: Reduced protein adsorption on novel phospholipid polymers

AUTHOR(S): Ishihara, Kazuhiko; Iwasaki, Yasuhiko

CORPORATE SOURCE: Department of Materials Science, Graduate School of  
Engineering, The University of Tokyo, Tokyo, 113-8656,  
Japan.

SOURCE: Journal of Biomaterials Applications, (1998) Vol. 13, No.  
2, pp. 111-127.

CODEN: JBAPEL. ISSN: 0885-3282.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1998:638100

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020509

ED Entered STN: 20011116

Last Updated on STN: 20020509

AB We have synthesized phospholipid polymers containing 2-methacryloyloxyethyl phosphorylcholine (MPC) units as novel blood compatible polymers and have evaluated their interactions with blood components. It was found that in the absence of **anticoagulants**, blood clotting was delayed and blood cell adhesion and activation were effectively prevented on the MPC copolymer surface. A little amount of protein adsorbed on the MPC copolymer from human plasma was compared with conventional polymers, and the amount was reduced with increasing MPC unit fraction. To clarify the reason for the little protein adsorption on the MPC copolymer, the water structure in the hydrated polymer was examined with attention to the free water fraction. Hydration of the polymers occurred when they were immersed in water. The thermal anal. of these hydrated polymers revealed that the free water fractions in the poly(MPC-co-Bu methacrylate(BMA)) and poly(MPC-co-n-dodecyl methacrylate) were significantly larger than those in the poly(2-hydroxyethyl methacrylate) (HEMA). The conformation of proteins adsorbed on poly(HEMA) changed considerably but that on poly(MPC-co-BMA) was almost the same as the native state. We concluded from these results that the proteins are hardly adsorbed and do not change their original conformation on the polymer surfaces which possess a high free water fraction such as phospholipid polymers.

CC 63-7

ST Miscellaneous Descriptors

protein adsorption phospholipid polymer methacrylate

RN 125275-25-4 (3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium,  
4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with  
butyl 2-methyl-2-propenoate)

144514-07-8 (3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium,  
4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with  
dodecyl 2-methyl-2-propenoate)  
9003-63-8 (Poly(butyl methacrylate))  
25249-16-5 (Poly(2-hydroxyethyl methacrylate))  
67881-98-5 (2-Methacryloyloxyethyl phosphorylcholine)

L82 ANSWER 37 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN DUPLICATE 17  
ACCESSION NUMBER: 1995:55014 TOXCENTER  
DOCUMENT NUMBER: PubMed ID: 7622528  
TITLE: Adhesion and cytokine production by monocytes on  
poly(2-methacryloyloxyethyl phosphorylcholine-co-alkyl  
methacrylate)-coated polymers  
AUTHOR(S): DeFife K M; Yun J K; Azeez A; Stack S; Ishihara K;  
Nakabayashi N; Colton E; Anderson J M  
CORPORATE SOURCE: Department of Pathology, Case Western Reserve University,  
Cleveland, Ohio 44106, USA  
CONTRACT NUMBER: HL 33849 (NHLBI)  
HL 48771 (NHLBI)  
SOURCE: Journal of biomedical materials research, (1995 Apr) 29  
(4) 431-9.  
Journal Code: 0112726. ISSN: 0021-9304.  
COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
FILE SEGMENT: MEDLINE  
OTHER SOURCE: MEDLINE 95348144  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20011116  
Last Updated on STN: 20011116  
ED Entered STN: 20011116  
Last Updated on STN: 20011116  
AB Human monocytes isolated from peripheral venous blood were **assayed**  
for their ability to adhere to various polymers. The culture supernatants  
were also **assayed** for the cytokines, interleukin-1 beta  
(IL-beta), interleukin-6 (IL-6), and tumor necrosis factor-alpha  
(TNF-alpha). The polymers evaluated for adherence and cytokine production  
included Pellethane, polyethylene and poly[n-butyl methacrylate (BMA)]  
coated with poly[2-methacryloyloxyethyl phosphorylcholine (MPC)-co-alkyl  
methacrylate] copolymers. In some experiments the test polymers were  
adsorbed with fibrinogen or IgG prior to the addition of monocytes. MPC  
copolymer-coated materials inhibited monocyte and macrophage adhesion  
after 1 and 8 days of culture relative to corresponding uncoated polymers  
and tissue culture polystyrene (TCPS). The degree of inhibition by coated  
Pellethane compared to uncoated Pellethane was the greatest, while  
inhibition of adhesion by coated poly(BMA) was the least compared to  
uncoated poly(BMA). However, adhesion was significantly decreased on both  
coated and uncoated poly(BMA) by day 8. While IL-1 beta, IL-6, and  
TNF-alpha release was variably influenced by polymer coating, release was  
consistently inhibited relative to TCPS on day 1. However, cytokine  
production was not inhibited compared to corresponding uncoated polymers  
on day 1. With or without protein preadsorption, IL-1 beta release was  
not detectable in the supernatants of any polymer on day 8, IL-6  
production was diminished on day 8, and TNF-alpha production was sustained  
on day 8. Overall, MPC copolymer-coated and uncoated poly(BMA) were the  
least stimulating, while TCPS was the most stimulating. (ABSTRACT TRUNCATED  
AT 250 WORDS)  
CT Biocompatible Materials  
Cell Adhesion  
Cells, Cultured  
Humans

\*Interleukin-1: SE, secretion  
\*Interleukin-6: SE, secretion  
\*Methacrylates  
\*Monocytes: CY, cytology  
\*Phosphorylcholine: AA, analogs & derivatives  
  Polymethacrylic Acids  
  Polystyrenes  
  Polyurethanes  
  Research Support, U.S. Gov't, Non-P.H.S.  
  Research Support, U.S. Gov't, P.H.S.  
  Surface Properties  
\*Tumor Necrosis Factor-alpha: SE, secretion  
RN 107-73-3 (Phosphorylcholine)  
  125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate))  
  9003-63-8 (polybutyl methacrylate)  
CN 0 (Biocompatible Materials); 0 (Interleukin-1); 0 (Interleukin-6); 0 (Methacrylates); 0 (Polymethacrylic Acids); 0 (Polystyrenes); 0 (Polyurethanes); 0 (Tumor Necrosis Factor-alpha)

L82 ANSWER 38 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN  
ACCESSION NUMBER: 2005:77416 TOXCENTER  
COPYRIGHT: Copyright 2005 ACS  
DOCUMENT NUMBER: CA14223435647T  
TITLE: Structural Study of DNA Condensation Induced by Novel Phosphorylcholine-Based Copolymers for Gene Delivery and Relevance to DNA Protection  
AUTHOR(S): Chim, Y. T. A.; Lam, J. K. W.; Ma, Y.; Armes, S. P.; Lewis, A. L.; Roberts, C. J.; Stolnik, S.; Tendler, S. J. B.; Davies, M. C.  
CORPORATE SOURCE: Laboratory of Biophysics and Surface Analysis, School of Pharmacy, The University of Nottingham, Nottingham, NG7 2RD, UK.  
SOURCE: Langmuir, (2005) Vol. 21, No. 8, pp. 3591-3598.  
CODEN: LANGD5. ISSN: 0743-7463.  
COUNTRY: UNITED KINGDOM  
DOCUMENT TYPE: Journal  
FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 2005:180198  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20050308  
Last Updated on STN: 20050531  
ED Entered STN: 20050308  
Last Updated on STN: 20050531  
AB Poly[2-(dimethylamino)ethyl methacrylate-b-2-methacryloyloxyethyl phosphorylcholine] (DMA-MPC) is currently under investigation as a new vector candidate for gene therapy. The DMA block has been previously demonstrated to condense DNA effectively. The MPC block contains a phosphorylcholine (PC) headgroup, which can be found naturally in the outside of the cell membrane. This PC-based polymer is extremely hydrophilic and acts as a biocompatible steric stabilizer. In this study, we assess in detail the morphologies of DNA complexes obtained using the diblock copolymer series DMAxMPC30 (where the mean d.p. of the MPC block was fixed at 30 and the DMA block length was systematically varied) using TEM and liquid atomic force microscopy (AFM). Both techniques indicate more compact complex morphologies (more efficient condensation) as the length of the cationic DMA block increases. However, the detailed morphologies of the DMAxMPC30-DNA complexes observed by TEM in vacuo and by AFM in aqueous medium are different. This phenomena is believed to be related to the highly hydrophilic nature of the MPC block. TEM studies revealed that the



morphol. of the complexes changes from loosely condensed structures to highly condensed rods, toroids, and oval-shaped particles as the DMA moiety increases. In contrast, morphol. changes from plectonemic loops to flowerlike and rectangular blocklike structures, with an increase in highly condensed central regions, are observed by in situ AFM studies. The relative population of each structure is clearly dependent on the polymer mol. composition. Enzymic degradation **assays** revealed that only the DMA homopolymer provided effective DNA protection against DNase I degradation, while other highly condensed copolymer complexes, as judged from TEM and gel electrophoresis, only partially protected the DNA. However, AFM images indicated that the same highly condensed complexes have less condensed regions, which we believe to be the initiation sites for enzymic attack. This indicates that the open structures observed by AFM of the DNA complexation by the DMAxMPC30 copolymer series are closer to in vivo morphol. when compared to TEM.

CC 63-6

ST Miscellaneous Descriptors

DNA phosphorylcholine copolymer gene delivery  
RN 409334-34-5 (2-(Dimethylamino)ethyl methacrylate-2-methacryloyloxyethyl phosphorylcholine block copolymer)

L82 ANSWER 39 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2005:246857 TOXCENTER

COPYRIGHT: Copyright (c) 2005 The Thomson Corporation

DOCUMENT NUMBER: 42-16504

TITLE: Novel biocompatible phosphorylcholine-based self-assembled nanoparticles for drug delivery

AUTHOR(S): Salvage, JP; Rose, SF; Phillips, GJ; Lloyd, AW; Lewis, AL; et al

CORPORATE SOURCE: Univ Brighton, Sch Pharm &amp; Biomol Sci, Moulsecoomb, Brighton BN2 4GJ, E Sussex, England

SOURCE: Journal of Controlled Release (Netherlands), (2005) Vol. 104, pp. 259-270. 35 Refs.

CODEN: JCREEC. ISSN: 0168-3659.

DOCUMENT TYPE: Journal

FILE SEGMENT: IPA

OTHER SOURCE: IPA 2005:16478

LANGUAGE: English

ENTRY DATE: Entered STN: 20050920

Last Updated on STN: 20050920

ED Entered STN: 20050920

Last Updated on STN: 20050920

AB Major challenges associated with nano-sized drug delivery systems include removal from systemic circulation by phagocytic cells and controlling appropriate drug release at target sites. 2-methacryloyloxyethyl phosphorylcholine (MPC) has been copolymerised in turn with two pH responsive comonomers (2-(diethylamino)ethyl methacrylate (DEA) and 2-(diisopropylamino)ethyl methacrylate (DPA), to develop novel biocompatible drug delivery vehicles. Micelles were prepared from a series of copolymers with varying block compositions and their colloidal stability and dimensions were assessed over a range of solution pH using photon correlation spectroscopy. The drug loading capacities of these micelles were evaluated using Orange OT dye as a model compound. The cytotoxicity of the micelles was assessed using an in vitro **assay**. The MPC-DEA diblock copolymers formed micelles at around pH 8 and longer DEA block lengths allowed higher drug loadings. However, these micelles were not stable at physiological pH. In contrast, MPC-DPA diblock copolymers formed micelles of circa 30 nm diameter at physiological pH. In vitro **assays** indicated that these MPC-DPA diblock copolymers had negligible cytotoxicities. Thus novel non-toxic

biocompatible micelles of appropriate size and good colloidal stability with pH-modulated drug uptake and release can be readily produced using MPC-DPA diblock copolymers. (C) 2005 Elsevier B.V All rights reserved.

SC 9 Pharmaceutics

ST Miscellaneous Descriptors

2-Methacryloyloxyethyl phosphorylcholine; micelles  
2-(Diisopropylamino)ethyl methacrylate; micelles  
N,N-Diethylaminoethyl methacrylate; micelles  
Nanoparticles; 2-methacryloyloxyethyl phosphorylcholine  
Copolymers; 2-methacryloyloxyethyl phosphorylcholine  
Stability; 2-methacryloyloxyethyl phosphorylcholine  
Micelles; 2-methacryloyloxyethyl phosphorylcholine  
Particle size; 2-methacryloyloxyethyl phosphorylcholine  
Toxicity; 2-methacryloyloxyethyl phosphorylcholine  
Nanoparticles; 2-(diisopropylamino)ethyl methacrylate  
Copolymers; 2-(diisopropylamino)ethyl methacrylate  
Stability; 2-(diisopropylamino)ethyl methacrylate  
Micelles; 2-(diisopropylamino)ethyl methacrylate  
Particle size; 2-(diisopropylamino)ethyl methacrylate  
Toxicity; 2-(diisopropylamino)ethyl methacrylate  
Nanoparticles; n,n-diethylaminoethyl methacrylate  
Copolymers; n,n-diethylaminoethyl methacrylate  
Stability; n,n-diethylaminoethyl methacrylate  
Micelles; n,n-diethylaminoethyl methacrylate  
Particle size; n,n-diethylaminoethyl methacrylate  
Toxicity; n,n-diethylaminoethyl methacrylate

RN 67881-98-5 (2-Methacryloyloxyethyl phosphorylcholine)  
16715-83-6 (2-(Diisopropylamino)ethyl methacrylate)  
105-16-8 (N,N-Diethylaminoethyl methacrylate)

L82 ANSWER 40 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2005:44317 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA14218341817W

TITLE: In situ modification on cellulose acetate hollow fiber  
membrane modified with phospholipid polymer for biomedical  
application

AUTHOR(S): Ye, Sang Ho; Watanabe, Junji; Iwasaki, Yasuhiko; Ishihara,  
Kazuhiko

CORPORATE SOURCE: Department of Materials Engineering, School of  
Engineering, The University of Tokyo, 7-3-1 Hongo,  
Bunkyo-ku, Tokyo, 113-8656, Japan.

SOURCE: Journal of Membrane Science, (2005) Vol. 249, No. 1-2, pp.  
133-141.

CODEN: JMESDO. ISSN: 0376-7388.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2005:124799

LANGUAGE: English

ENTRY DATE: Entered STN: 20050215

Last Updated on STN: 20050426

ED Entered STN: 20050215

Last Updated on STN: 20050426

AB The hollow fiber membrane (HFM) made from synthetic polymers need  
improvement in terms of hemocompatibility or biocompatibility, for use in  
the medical field. In this study, cellulose acetate (CA) HFM modified  
with the water-soluble amphiphilic 2-methacryloyloxyethyl phosphorylcholine  
(MPC) copolymer (poly (MPC-co-Bu methacrylate) (PMB80, MPC:BMA = 80:20  
(mol%)) was prepared by a dry-jet wet spinning process. The PMB80 was

coated on the CA HFM surface in situ during the phase inversion of the dope solution by using a PMB80 solution as an inner **coagulant**. The CA/PMB80 coating HFM showed no phys. structure changes in comparison with the CA HFM prepared using the same preparative conditions. The structure and permeability of the CA/PMB80 coating HFM was controllable by changing the preparative conditions. From the results of the X-ray photoelectron spectroscopic (XPS) observations, the amount of modification was changed with the concentration of PMB80 in the **coagulant**. The XPS signal attributed to the phosphorus atom of the PMB80 remained even after 1 mo of rinsing with distilled water. Also, the CA/PMB80 coated HFM showed good permeability and a low membrane fouling property in comparison with the non-modified CA HFM, due to the low protein adsorption property of the PMB80.

CC 63-8

ST Miscellaneous Descriptors

cellulose acetate hollow fiber membrane phospholipid polymer biomedical

RN 125275-25-4 (Butyl methacrylate-2-methacryloyloxyethyl phosphorylcholine copolymer)

RN 9004-35-7

L82 ANSWER 41 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:108846 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA14118296374R

TITLE: Amphiphilic block copolymers based on poly(2-acryloyloxyethyl phosphorylcholine) prepared via RAFT polymerisation as biocompatible nanocontainers

AUTHOR(S): Stenzel, Martina H.; Barner-Kowollik, Christopher; Davis, Thomas P.; Dalton, Helen M.

CORPORATE SOURCE: Centre for Advanced Macromolecular Design, School of Chemical Engineering and Industrial Chemistry, The University of New South Wales, Sydney, NSW 2031, Australia.

SOURCE: Macromolecular Bioscience, (2004) Vol. 4, No. 4, pp. 445-453.

CODEN: MBAIBU. ISSN: 1616-5187.

COUNTRY: AUSTRALIA

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2004:371833

LANGUAGE: English

ENTRY DATE: Entered STN: 20040511

Last Updated on STN: 20041229

ED Entered STN: 20040511

Last Updated on STN: 20041229

AB Amphiphilic block copolymers composed of poly(Bu acrylate) and poly(2-acryloyloxyethyl phosphorylcholine) have been prepared using reversible addition fragmentation transfer (RAFT) polymerization The conversion of

the polymerization was determined using online FT NIR spectroscopy. NMR spectroscopy

was used not only to support the results obtained from FT NIR spectroscopy but also prove the formation of micelles. Due to the strong aggregation tendency of these block copolymers and the resulting difficulties concerning the mol. weight anal. test expts. were carried out replacing poly(2-acryloyloxyethyl phosphorylcholine) with poly(2-hydroxyethyl acrylate). Micelle size and the aggregation behavior were investigated using dynamic light scattering. The sizes of the nanocontainers obtained were found to be influenced by the block length as well as the solvent leading to micelles in the range between 40 and 160 nm. The toxicity of

the RAFT agent used was then analyzed by cell growth inhibition tests.

CC 35-7

ST Miscellaneous Descriptors

RAFT polymn block amphiphilic polyacrylate synthesis; aggregation  
micelle hydrodynamic radius acryloyloxyethyl phosphorylcholine block  
copolymer; cytotoxicity **assay** phenylmethylthiothioxomethylthi  
o propanoic acid RAFT initiator

RN 765206-20-0 (Butyl acrylate-2-acryloyloxyethyl phosphorylcholine  
diblock copolymer)

765276-02-6 (Butyl acrylate-2-hydroxyethyl acrylate diblock copolymer)

RN 497931-76-7

L82 ANSWER 42 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:84930 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA14115248491E

TITLE: Phosphorylcholine-containing polymers for use in cell  
encapsulation

AUTHOR(S): Yang, Ying; Zhang, Sifu; Jones, Graham; Morgan, Noel; El  
Haj, Alicia J.

CORPORATE SOURCE: Centre for Science and Technology in Medicine, School of  
Medicine, Keele University/University Hospital of North  
Staffordshire, Stoke-on-Trent, Staffs, UK.

SOURCE: Artificial Cells, Blood Substitutes, and Biotechnology,  
(2004) Vol. 32, No. 1, pp. 91-104.

CODEN: ACBSDA.

COUNTRY: UNITED KINGDOM

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2004:278630

LANGUAGE: English

ENTRY DATE: Entered STN: 20040413

Last Updated on STN: 20041229

ED Entered STN: 20040413

Last Updated on STN: 20041229

AB A model system for encapsulation of pancreatic islets which has potential  
properties for improving biocompatibility and **immunosuppression**  
was investigated. In vitro and in vivo studies have shown that  
phosphorylcholine-containing polymers have high biocompatibility due to low  
adsorption of proteins and reduced thrombus formation. Encapsulation of  
islets isolated from rats with a compound membrane composed of  
phosphorylcholine-containing polymers and cellulose acetate led to rapid  
insulin production and diffusion across the membrane in response to glucose  
challenge. The phosphorylcholine-containing polymer had a mol. weight of about  
1.3 + 10<sup>4</sup> Da. The polymer-coated membrane excluded larger mols.  
such as IgG (mol. weight 150 kDa), thereby acting as a phys. **immuno**  
-barrier, but allowed smaller mols. such as glucose and insulin to pass  
through.

CC 63-5

ST Miscellaneous Descriptors

phosphorylcholine encapsulation pancreatic islet

RN 9004-10-8 (Insulin)

9004-35-7 (Cellulose acetate)

125275-25-4 (2-MethacryloyloxyethylPhosphorylcholine-butyl  
methacrylate copolymer)

L82 ANSWER 43 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2004:149311 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA14103042760H

TITLE: In vitro and ex vivo blood compatibility study of  
2-methacryloyloxyethyl phosphorylcholine (MPC)  
copolymer-coated hemodialysis hollow fibers

AUTHOR(S): Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo  
Medical and Dental University, Chiyoda-ku, Tokyo,  
101-0062, Japan.

SOURCE: Journal of Artificial Organs, (2003) Vol. 6, No. 4, pp.  
260-266.  
CODEN: JAORFN. ISSN: 1434-7229.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2003:1000315

LANGUAGE: English

ENTRY DATE: Entered STN: 20040707  
Last Updated on STN: 20041214

ED Entered STN: 20040707  
Last Updated on STN: 20041214

AB To identify the advantages of 2-methacryloyloxyethyl phosphorylcholine  
(MPC) copolymer-coated polysulfone (PSf) hollow fibers for hemodialyzer  
and hemofilter minimodules with hollow fibers were made and blood  
compatibility was evaluated in vitro and ex vivo. Three types of hollow  
fibers, i.e., pure PSf (no additives), PSf alloyed with  
poly(1-vinyl-2-pyrrolidone) (PVPy), and PSf coated with the MPC copolymer,  
were processed in wet conditions. Com. available hollow fibers (APS) were  
used as a control sample. The PSf hollow fibers have a condensed  
structure. A porous structure was observed when the PVPy was alloyed before  
wet processing, and no effect of the innercoated MPC copolymer on the  
porous structure was observed. One-tenth-sized minimodules of the  
conventional hemodialyzer were fabricated with 200 fibers each. The  
solute permeability of the hollow fibers was evaluated using 10% bovine  
serum in a buffer solution containing cytochrome C, which is a model protein of  
B2-microglobulin. After circulation for 2.5 h, the solute  
permeability of APS and PVPy-alloyed PSf hollow fibers decreased to 50%  
compared with their initial values. In contrast, the value for the hollow  
fibers innercoated with the MPC copolymer maintained its initial level.  
The inner surface of the dialysis membranes was observed with a transmission  
electron microscope and a layer of adsorbed protein on the PSf, APS, and  
PVPy-alloyed PSf hollow fibers was observed, but not on the MPC  
copolymer-coated fibers. Blood cell adhesion was then evaluated by  
circulation of whole rabbit blood without any **anticoagulant** ex  
vivo. Many adherent cells were observed on the PVPy-alloyed PSf hollow  
fibers; however, blood cells did not adhere or aggregate on the MPC  
copolymer-coated hollow fibers. From these results, we concluded that the  
in-situ coating of MPC copolymer on PSf hollow fibers is effective in  
preventing blood **coagulation** and maintaining the solute  
permeability of the fibers.

CC 63-7

ST Miscellaneous Descriptors  
blood compatibility methacrylic phosphorylcholine copolymer polysulfone  
hemodialysis hollow fiber

RN 9007-43-6 (Cytochrome C)  
9003-39-8 (Poly(1-vinyl-2-pyrrolidone))

RN 393587-07-0

L82 ANSWER 44 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2003:306221 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA14009133678H

TITLE: Nonthrombogenic hemodialyzer with MPC copolymer  
AUTHOR(S): Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Ishihara, Kazuhiko  
CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo  
Medical and Dental University, Kanda-surugadai,  
Chiyoda-ku, Tokyo, 101-0062, Japan.  
SOURCE: Advances in Science and Technology (Faenza, Italy), (2003)  
Vol. 41, No. Materials in Clinical Applications VI, pp.  
161-170.  
CODEN: ASETE5.  
COUNTRY: JAPAN  
DOCUMENT TYPE: Journal  
FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 2003:445247  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20031230  
Last Updated on STN: 20050628

ED Entered STN: 20031230

Last Updated on STN: 20050628

AB Development of non-thrombogenic dialysis membrane was aimed. Copolymers of 2-methacryloyloxyethyl phosphorylcholine (MPC) showed good hemocompatibility. They were introduced on the blood-contacting surface of hemodialysis cellulose and polysulfone membranes without adverse effect on mech. properties and permeability. Ex-vivo single pass of rabbit blood through the mini-modules without an **anticoagulant** was performed for 0.5 h and blood cells did not attach on the surface. It was concluded that the modified hollow fibers with MPC copolymers are promising to develop a hemodialyzer, which does not require **anticoagulants**. MPC polymer surface could form self-assembled biomimetic membrane by accumulating phospholipids from blood stream and show non-thrombogenicity. MPC is available from NOF Co., in Tokyo.

CC 63-7

ST Miscellaneous Descriptors

methacryloyloxyethyl phosphorylcholine copolymer hemodialyzer membrane  
cellulose polysulfone blood biocompatibility; platelet antithrombogenic  
hemodialysis membrane implant methacryloyloxyethyl phosphorylcholine  
copolymer

RN 60-27-5 (Creatinine)

9007-43-6 (Cytochrome C)

**67881-99-6** (2-Methacryloyloxyethyl-phosphorylcholine polymer)

9004-67-5 (Methyl cellulose)

RN **142146-61-0**; **182816-96-2**

L82 ANSWER 45 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2003:196361 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13908122556V

TITLE: Preparation and blood compatibility of  
phosphorylcholine-bonded O-butyrylchitosan

AUTHOR(S): Zhu, Aiping; Shan, Bing; Yuan, Youling; Shen, Jian

CORPORATE SOURCE: Department of Polymer Science and Engineering, Nanjing  
University, Nanjing, 210093, Peop. Rep. China.

SOURCE: Polymer International, (2003) Vol. 52, No. 1, pp. 81-85.  
CODEN: PLYIEI. ISSN: 0959-8103.

COUNTRY: CHINA

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2003:134342

LANGUAGE: English

ENTRY DATE: Entered STN: 20030812

Last Updated on STN: 20030819

ED Entered STN: 20030812  
Last Updated on STN: 20030819

AB 2-Methacryloyloxyethyl phosphorylcholine (MPCE) was synthesized using phosphorus trichloride, ethylene glycol, 2-hydroxyethyl methacrylate and triethylamine, and then used in the preparation of O-butyrylchitosan-bonded MPCE (MPCE-BCS) by Michael addition of MPCE to amino groups of O-butyrylchitosan. The structures of MPCE and MPCE-BCS were characterized by FTIR and <sup>1</sup>H NMR. The blood-compatibility of MPCE-BCS was evaluated by means of blood clotting and platelet adhesion **assays**. The blood-clotting **assay** indicated that O-butyrylchitosan was haemocompatible. Both the blood-clotting **assay** and platelet adhesion **assay** confirmed that MPCE-BCS had excellent antithrombogenicity.

CC 63-5

ST Miscellaneous Descriptors  
chitosan phosphorylcholine deriv blood compatibility  
**anticoagulant**

RN 9012-76-4Q (Chitosan, reaction product with 2-Methacryloyloxyethyl phosphorylcholine)  
106-31-0 (Butyric anhydride)  
107-21-1 (Ethylene glycol)  
121-44-8 (Triethylamine)  
868-77-9 (2-Hydroxyethyl methacrylate)  
7719-12-2 (Phosphorus trichloride)  
9012-76-4 (Chitosan)  
822-39-9 (2-Chloro-1,3,2-dioxaphospholane)  
6609-64-9 (2-Chloro-2-oxo-1,3,2-dioxaphospholane)

RN **312490-87-2**; 124384-94-7; 82793-19-9

L82 ANSWER 46 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2000:171883 TOXCENTER  
COPYRIGHT: Copyright 2005 ACS  
DOCUMENT NUMBER: CA13322313528Y  
TITLE: Photoinduced graft polymerization of 2-methacryloyloxyethyl phosphorylcholine on polyethylene membrane surface for obtaining blood cell adhesion resistance

AUTHOR(S): Ishihara, K.; Iwasaki, Y.; Ebihara, S.; Shindo, Y.; Nakabayashi, N.

CORPORATE SOURCE: Graduate School of Engineering, Department of Materials Science, The University of Tokyo, Tokyo, 113-8656, Japan.

SOURCE: Colloids and Surfaces, B: Biointerfaces, (2000) Vol. 18, No. 3,4, pp. 325-335.  
CODEN: CSBBEQ. ISSN: 0927-7765.

COUNTRY: JAPAN  
DOCUMENT TYPE: Journal  
FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 2000:507676  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20011116  
Last Updated on STN: 20020409

ED Entered STN: 20011116  
Last Updated on STN: 20020409

AB Phospholipid polymer, poly[2-methacryloyloxyethyl phosphorylcholine (MPC)], was grafted with polyethylene (PE) membrane using photoinduced polymerization technique to make the membrane resistant to cell adhesion. The water contact angle on the PE membrane grafted with poly(MPC) decreased with an increase in the photopolymerization time. This decrease corresponded to the increase in the amount of poly(MPC) grafted on the PE surface. The same graft polymerization procedure was applied using other hydrophilic monomers, such

as acrylamide (AAm), N-vinylpyrrolidone (VPy) and methacryloyl poly(ethylene glycol) (MPEG). These monomers were also polymerized to form grafted chains on the PE membrane, and the grafting was confirmed with XPS. Anal. of amount and distribution of plasma proteins at the plasma-contacting surface of the original and the modified PE membranes were analyzed using **immunogold assay**. The grafting of poly(MPC) and poly(VPy) on PE membrane reduced the plasma protein adsorption significantly compared with that on the original PE membrane. However, the PE membranes grafted with poly(AAm) or poly(MPEG) did not show any effects on protein adsorption. Platelet adhesion on the original and modified PE membranes from platelet-rich plasma was also examined. A large number of platelets adhered and activated on the original PE membrane. Grafting with poly(AAm) did not suppress platelet adhesion, but grafting with poly(MPC) or poly(VPy) on the PE membrane was effective in preventing platelet adhesion. It is concluded that the introduction of the phosphorylcholine group on the surface could decrease the cell adhesion to substrate polymer.

CC 63-7

ST Miscellaneous Descriptors

photopolymer methacryloxyloxyethyl phosphorylcholine polyethylene  
biocompatibility; blood adhesion photopolymer  
methacryloxyloxyethylphosphorylcholine polyethylene

RN 108144-73-6 (3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium,  
4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with  
2-propenamide, graft)

**252877-49-9** (3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium,  
4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with  
ethene, graft)

RN 176587-89-6; 220830-40-0

L82 ANSWER 47 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2002:129717 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13625390798A

TITLE: Bioinspired phospholipid polymer biomaterials for making  
high performance artificial organs

AUTHOR(S): Ishihara, K.

CORPORATE SOURCE: Department of Materials Science, Graduate School of  
Engineering, The University of Tokyo, Tokyo, 113-8656,  
Japan.SOURCE: Science and Technology of Advanced Materials, (2000) Vol.  
1, No. 3, pp. 131-138.

CODEN: STAMCV. ISSN: 1468-6996.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2001:678942

LANGUAGE: English

ENTRY DATE: Entered STN: 20020612

Last Updated on STN: 20020618

ED Entered STN: 20020612

Last Updated on STN: 20020618

AB A review. Novel polymer biomaterials, which can be used in contact with  
blood, are prepared with strong inspiration from the surface structure of  
biomembrane. That is, the polymers with a phospholipid polar group in the  
side chain, 2-methacryloxyloxyethyl phosphorylcholine (MPC) polymers were  
synthesized. The MPC polymers can inhibit surface-induced clot formation  
effectively, when they are in contact with blood even in the absence of an  
**anticoagulant**. This phenomenon was due to the reduction of plasma  
protein and suppression of denaturation of adsorbed proteins, that is the



MPC polymers interact with blood components very mildly. As the mol. structure of the MPC polymer was easily designed by changing the monomer units and their composition, it could be applied to surface modification of artificial organs and biomedical devices for improving blood and tissue compatibility. Thus, the MPC polymers are useful polymer biomaterials for manufacturing high performance artificial organs and biomedical devices to provide safe medical treatments.

CC 63-0

ST Miscellaneous Descriptors

review phospholipid polymer biomaterial; methacryloyloxyethyl  
phosphocholine polymer biomaterial review

RN 67881-99-6 (Poly(2-methacryloyloxyethylphosphocholine))

L82 ANSWER 48 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2000:156977 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13324340029T

TITLE: Control of cell adhesion and proliferation on MPC-BMA  
copolymer surface

AUTHOR(S): Watanabe, Akihiko; Iwasaki, Yasuhiko; Nakabayashi, Nobuo;  
Ishihara, Kazuhiko

CORPORATE SOURCE: Dep. of Org. Mater., Div. of Biomater. Inst. of Biomater.  
and Bioeng., Tokyo Med. and Dent. Univ., Japan.

SOURCE: Seitai Zairyo Kogaku Kenkyusho Hokoku (Tokyo Ika Shika  
Daigaku), (2000) Vol. 33, pp. 38-43.  
CODEN: SZKHF9.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 2000:393744

LANGUAGE: Japanese

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020403

ED Entered STN: 20011116

Last Updated on STN: 20020403

AB A review with 12 refs. For many medical applications of biomaterials, reduction or elimination of cell adhesion is desirable to enhance biocompatibility (Drumheller and Hubell, 1995). Adhesion of cells to an implant surface results from the protein adsorbed there. If protein adsorption can be prevented, cellular attachment will be suppressed. Copolymers of methacryloyloxyethyl phosphorylcholine (MPC) have affinity for phospholipids due to the phosphorylcholine polar groups on the MPC copolymer surface (Ishihara et al., 1990). When MPC copolymers are placed in contact with plasma, phospholipids are adsorbed and accumulated. They rearrange to create an organized layer which interacts mildly with proteins, thus preventing the adsorption of proteins on the material surface (Ishihara et al., 1992). The organized lipid layer on the material mimics a biomembrane surface. MPC copolymers have been investigated for blood-contacting applications and were found to possess excellent resistance against protein and platelet adhesion when exposed to platelet-rich plasma and whole blood even in the absence of an **anticoagulant** (Ishihara et al., 1991 and 1992). These copolymers may be suitable to for improving implant surfaces in applications where inhibition of cellular attachment is desired. In the present study, the adhesion of fibroblast cells in vitro to surfaces coated with an MPC copolymer was examined in comparison with a non-coated surface.

CC 63-0

ST Miscellaneous Descriptors

review biocompatible implant methacryloyloxyethylphosphorylcholine  
copolymer

RN 125275-25-4

L82 ANSWER 49 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN  
ACCESSION NUMBER: 1999:184688 TOXCENTER  
COPYRIGHT: Copyright 2005 ACS  
DOCUMENT NUMBER: CA13124327452Q  
TITLE: Modification of polysulfone with phospholipid polymer for improvement of the blood compatibility. Part 2. Protein adsorption and platelet adhesion  
AUTHOR(S): Ishihara, Kazuhiko; Fukumoto, Kikuko; Iwasaki, Yasuhiko; Nakabayashi, Nobuo  
CORPORATE SOURCE: Department of Materials Science, Graduate School of Engineering, The University of Tokyo, Tokyo, 113-8656, Japan.  
SOURCE: Biomaterials, (1999) Vol. 20, No. 17, pp. 1553-1559.  
CODEN: BIMADU. ISSN: 0142-9612.  
COUNTRY: JAPAN  
DOCUMENT TYPE: Journal  
FILE SEGMENT: CAPLUS  
OTHER SOURCE: CAPLUS 1999:569969  
LANGUAGE: English  
ENTRY DATE: Entered STN: 20011116  
Last Updated on STN: 20020509  
ED Entered STN: 20011116  
Last Updated on STN: 20020509  
AB Protein adsorption and platelet adhesion from human plasma on polysulfone (PSf) membranes modified with 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer were studied. The modification was carried out by blending of the MPC polymer in the PSf. The amount of protein adsorbed on the PSf/MPC polymer blend membrane was significantly decreased with an increase in the composition of the blended MPC polymer. The distribution of the specific proteins adsorbed on the membrane surface was also determined by a gold-colloid immunoassay. Albumin,  $\gamma$ -globulin and fibrinogen were observed on every membrane surface after contact with plasma. However, in the case of the blended membrane, the d. of the adsorbed proteins decreased compared with that of original PSf membrane. That is, the MPC polymer blended in the membrane could function as a protein-adsorption-resistant additive. The number of platelets adhered on the PSf membrane was reduced, and change in the morphol. of adherent platelets was also suppressed by the modification with the MPC polymer. Therefore, the PSf/MPC polymer blend membrane had improved blood compatibility compared with the PSf membrane.  
CC 63-7  
ST Miscellaneous Descriptors  
polysulfone membrane phospholipid polymer biocompatibility; protein adsorption polysulfone membrane phospholipid polymer; platelet adhesion polysulfone membrane phospholipid polymer  
RN 125275-25-4; 144514-07-8; 28776-63-8

L82 ANSWER 50 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN  
ACCESSION NUMBER: 1999:209043 TOXCENTER  
COPYRIGHT: Copyright 2005 ACS  
DOCUMENT NUMBER: CA13206069285H  
TITLE: The development of in vitro biocompatibility tests for the evaluation of intraocular biomaterials  
AUTHOR(S): Lloyd, A. W.; Dropcova, S.; Faragher, R. G. A.; Gard, P. R.; Hanlon, G. W.; Mikhalovsky, S. V.; Olliff, C. J.; Denyer, S. P.  
CORPORATE SOURCE: Drug Delivery & Biomaterials Research Group, School of Pharmacy and Biomolecular Sciences, University of

SOURCE: Brighton, Brighton, BN2 4GJ, UK.  
Journal of Materials Science: Materials in Medicine,  
(1999) Vol. 10, No. 10/11, pp. 621-627.  
CODEN: JSMMEI. ISSN: 0957-4530.

COUNTRY: UNITED KINGDOM

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1999:742634

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116  
Last Updated on STN: 20020403

ED Entered STN: 20011116  
Last Updated on STN: 20020403

AB Recent developments in ocular implant technol. require the in vitro evaluation of ocular compatibility in early stage development programs. This requires an understanding and appreciation of the biol. interactions which occur in the ocular environment and their relevance with respect to the clin. complications associated with surgical implantation of devices. This paper describes the development of a series of clin. reflective in vitro **assays** for assessing the potential ocular compatibility of novel intraocular lens materials. Staphylococcus epidermidis attachment, fibrinogen adsorption, mouse embryo fibroblast 3T3 adhesion and proliferation, primary rabbit lens cell adhesion, human peripheral blood macrophage adhesion and granulocyte activation tests were employed to evaluate two widely used intraocular biomaterials poly(Me methacrylate) (PMMA) and silicone, and a novel biomimetic phosphorylcholine-based coating (PC). The performance of these materials in the in vitro **assays** was compared to their ability to reduce postoperative inflammation in vivo in a rabbit model. The results demonstrated that the in vitro **assays** described here are predictive of in vivo ocular compatibility. These **assays** offer a more relevant means of assessing the ocular compatibility of biomaterials than those presently required by the authorities for regulatory approval of medical devices and implants.

CC 63-7

ST Miscellaneous Descriptors  
intraocular biomaterial biocompatibility polymer; cell adhesion  
intraocular biomaterial biocompatibility

RN 9011-14-7 (PMMA)  
67881-98-5Q (2-Methacryloyloxyethylphosphorylcholine, copolymers)

L82 ANSWER 51 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1999:148015 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13110134558V

TITLE: Competitive adsorption between phospholipid and plasma protein on a phospholipid polymer surface

AUTHOR(S): Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Nakatani, Masako; Mihara, Takashi; Kurita, Kimio; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Tokyo, 101-0062, Japan.

SOURCE: Journal of Biomaterials Science, Polymer Edition, (1999) Vol. 10, No. 5, pp. 513-529.  
CODEN: JBSEEA. ISSN: 0920-5063.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1999:310610

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020423

ED Entered STN: 20011116

Last Updated on STN: 20020423

AB The competitive adsorption of proteins and phospholipids on  $\omega$ -methacryloyloxyalkyl phosphorylcholine (MAPC) polymer was evaluated in this study. Albumin, fibrinogen, and dimyristoylphosphatidylcholine (DMPC) were used as model components. The amount of DMPC adsorbed on the MAPC polymers increased with an increase in the MAPC unit composition of the polymer. The methylene chain length of the MAPC unit was another factor influencing the DMPC adsorption when the MAPC unit composition of the MAPC polymer was low. The state of albumin and DMPC liposome adsorbed on the 2-methacryloyloxyethyl phosphorylcholine (MPC) polymer was determined by dynamic contact angle (DCA) measurement. The adsorption strength of albumin on the MPC polymer was weaker than that on the poly[n-Bu methacrylate (BMA)], i.e., the albumin was detached from the MPC polymer during the rinsing process. On the poly(BMA) surface, no difference in the shape of the DCA loops before and after contact with the DMPC liposomal suspension was observed. Fibrinogen adsorption on the MAPC polymer was detected by gold-colloid labeled **immunoassay**. The amount of fibrinogen adsorbed on every MAPC polymer surface was reduced by addition of the DMPC liposome in the fibrinogen solution. The number of platelets

adhered on the MAPC polymer was also decreased when the DMPC liposome was present in the fibrinogen solution during pretreatment. We concluded that phospholipids were preferentially adsorbed on the MAPC polymer surface compared with plasma protein and that the adsorbed phospholipids played an important role in showing an excellent blood compatibility on the MAPC polymer.

CC 63-7

ST Miscellaneous Descriptors

methacryloyloxyalkyl phosphorylcholine polymer adsorption phospholipid protein

RN 125275-25-4 (Butyl methacrylate-2-methacryloyloxyethyl phosphorylcholine copolymer)

18656-38-7 (Dimyristoylphosphatidylcholine)

RN 158760-97-5; 197913-20-5

L82 ANSWER 52 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1998:102775 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA12809102541F

TITLE: Zwitterionic and reactive group-containing vinyl polymers for blood-compatible surface coatings

AUTHOR(S): Bowers, Roderick W. J.; Jones, Stephen A.; Stratford, Peter W.; Charles, Stephen A.

CORPORATE SOURCE: ASSIGNEE: Biocompatibles Ltd.

PATENT INFORMATION: US 5705583 A 6 Jan 1998

SOURCE: (1998) U.S., 29 pp., Cont.-in-part of U.S. Ser. No. 175,348.

CODEN: USXXAM.

COUNTRY: UNITED KINGDOM

DOCUMENT TYPE: Patent

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1998:31161

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020605

ED Entered STN: 20011116

Last Updated on STN: 20020605

AB Polymers of  $\geq 1$  radical polymerizable monomers have pendant groups

bearing a center of permanent pos. charge and other pendant groups capable of stably binding the polymer to a surface, addnl. reactive groups in the polymer may serve as points for attachment of ligands to the polymer. The polymers may contain pendant groups which bind the polymer to a surface by physisorption, covalent bonding or ionic interactions. Thus, a coating of poly(2(methacryloyloxyethyl)-2'(trimethylammonium)ethyl phosphate inner salt-n-dodecyl methacrylate) (1:2) applied to PVC tubing showed no effect on blood pumping with or without **anticoagulant**, compared to an untreated PVC tubing which required **anticoagulant**.

CC 35-4

ST Miscellaneous Descriptors

zwitterionic acrylic polymer coating; ammonium phosphate zwitterionic acrylic polymer; biocompatible coating acrylic inner salt polymer; dodecyl methacrylate copolymer coating; crosslinkable monomer zwitterionic acrylic polymer; blood compatible coating zwitterionic acrylic polymer

RN 7719-12-2 (Phosphorous trichloride)

7429-90-5 (Aluminum)

9002-86-2 (PVC)

9002-88-4 (Polyethylene)

12597-69-2 (Steel)

107-21-1 (1,2-Ethanediol)

868-77-9 (2-Hydroxyethyl methacrylate)

920-46-7 (Methacryloyl chloride)

1120-71-4 (Propane sultone)

25265-75-2 (Butanediol)

41862-94-6 (Dodec-7-yn-1-ol)

RN 144514-07-8; 144514-08-9; 146109-91-3;

150120-09-5; 150120-10-8; 150120-11-9;

150120-12-0; 150120-14-2; 150120-16-4;

150120-17-5; 150120-18-6; 150120-19-7;

166195-20-6; 201359-42-4; 201359-43-5;

201359-44-6; 201425-85-6; 67881-98-5;

150120-13-1; 150120-15-3; 150205-71-3; 822-39-9; 6609-64-9;

82793-19-9; 92035-97-7; 144026-20-0; 144026-21-1; 150205-72-4; 818-61-1;

2867-47-2; 30030-25-2

L82 ANSWER 53 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1999:114772 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA13024329064T

TITLE: Preparation of phospholipid-accumulated surface for creation of a new biocompatible material

AUTHOR(S): Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Inst. Med. Dent. Eng., Tokyo Med. Dent. Univ., Tokyo, 101-0062, Japan.

SOURCE: Iyo Kizai Kenkyusho Hokoku (Tokyo Ika Shika Daigaku), (1998) Vol. 32, pp. 14-22.

CODEN: IKKHBS. ISSN: 0082-4739.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1999:86719

LANGUAGE: Japanese

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020521

ED Entered STN: 20011116

Last Updated on STN: 20020521

AB A review with 28 refs. Chemical compds. having phospholipid polar groups are interesting from the viewpoint not only of polymer chemical, but also of

biol., medical, and life sciences. Phospholipids, one of the main components of biomembrane, form a bilayer structure used to construct a biomembrane. Recently, many researchers have been attempting to prepare a phospholipid-accumulated surfaces to develop a new biocompatible material using phospholipid liposome, self-assembling monolayer (SAM), or polymeric phospholipids. In this review, we report the preparation of phospholipid-accumulated surfaces and their properties. We also prepared a methacrylate having a phosphorylcholine group, 2-methacryloyloxyethyl phosphorylcholine (MPC), to make a new biocompatible polymeric material. The competitive adsorption of proteins and phospholipids on the MPC polymer was evaluated using albumin, fibrinogen, and dimyristoyl phosphatidylcholine (DMPC) as model components, resp. The amount of DMPC adsorbed on the MPC polymers increased with an increase in the MPC composition in the polymer. The state of albumin and DMPC liposome adsorbed on the MPC polymer was determined by dynamic contact angle (DCA) measurement. The albumin adsorbed on the MPC polymer was much easily detached compared with that on the poly(Bu methacrylate). On the poly(Bu methacrylate) surface, no difference in the shape of the DCA loops before and after contact with the DMPC liposomal suspension was observed. Fibrinogen adsorption on the MPC polymer was detected by gold-colloid labeled **immunoassay**. The amount of fibrinogen adsorbed on the MPC polymer surface with 10 mol% MPC was reduced by addition of the DMPC liposome in the fibrinogen solution. The application of MPC polymer for improving the tribol. property of elastomer, new material for orthopedic bearing, is also reported in this review. We concluded that the creation of a phospholipid-accumulated surface is very useful method for obtaining biocompatibility on biomedical polymers.

CC 63-0

ST Miscellaneous Descriptors

review phospholipid surface biocompatible medical goods

RN 67881-99-6 (2-Methacryloyloxyethyl phosphorylcholine polymer)

L82 ANSWER 54 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1997:191266 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA12719267962M

TITLE: Reduction of surface-induced platelet activation on phospholipid polymer

AUTHOR(S): Iwasaki, Yasuhiko; Mikami, Asako; Kurita, Kimio; Yui, Nobuhiko; Ishihara, Kazuhiko; Nakabayashi, Nobuo

CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Tokyo, 101, Japan.

SOURCE: Journal of Biomedical Materials Research, (1997) Vol. 36, No. 4, pp. 508-515.

CODEN: JBMRBG. ISSN: 0021-9304.

COUNTRY: JAPAN

DOCUMENT TYPE: Journal

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1997:588729

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020618

ED Entered STN: 20011116

Last Updated on STN: 20020618

AB  $\omega$ -Methacryloyloxyalkyl phosphorylcholine (MAPC) polymers which were synthesized with attention to the surface structure of a biomembrane show excellent blood compatibility, i.e., resistance to protein adsorption and blood cell adhesion. To clarify the stability of platelets in contact with the MAPC polymer surfaces, cytoplasmic free calcium concentration ([Ca<sup>2+</sup>]<sub>i</sub>)

in the platelets was measured. A platelet suspension was passed through a column packed with various polymer beads after treatment with plasma, and the  $[Ca^{2+}]_i$  in the platelets eluted from the column was measured. The  $[Ca^{2+}]_i$  in contact with the MAPC polymers, i.e., poly[2-methacryloyloxyethyl phosphorylcholine-co-Bu methacrylate (BMA)] (PMEB) and poly(6-methacryloyloxyhexyl phosphorylcholine-co-BMA) (PMHB), was less than that in contact with poly(BMA). However, poly(10-methacryloyloxydecyl phosphorylcholine-co-BMA) (PMDB) was not effective in suppressing the increase in  $[Ca^{2+}]_i$ , and thus was at the same level as in the poly(BMA). Platelets in contact with PMEB or PMHB were less activated compared with those in contact with PMDB and poly(BMA). The state of the platelets adhered to these polymer surfaces, both morphol. and immunol., was examined SEM observation of the polymer surface after contact with a platelet suspension revealed that many platelets adhered and changed their shape on the poly(BMA). The nos. of adherent platelets were reduced on all MAPC polymer surface. The relative amount of  $\alpha$ -granule membrane glycoprotein (GMP-140) which appears on the cell membrane by activation of platelets on the PMEB surfaces was less than that on poly(BMA) and poly(2-hydroxyethyl methacrylate). Thus, PMEB and PMHB suppressed not only platelet adhesion but also activation of the platelets in contact with these surfaces.

CC 63-7  
ST Miscellaneous Descriptors  
phospholipid polymer surface platelet activation  
RN 125275-25-4; 158760-97-5

L82 ANSWER 55 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1996:183031 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA12512150849X

TITLE: Preparation of self-assembled biomimetic membranes and their functions

AUTHOR(S): Nakabayashi, Nobuo

CORPORATE SOURCE: Institute Medical and Dental Engineering, Tokyo Medical and Dental University, Tokyo, 101, Japan.

SOURCE: Advanced Biomaterials in Biomedical Engineering and Drug Delivery Systems, [Iketani Conference on Biomedical Polymers], 5th, Kagoshima, Japan, Apr. 18-22, 1995, (1996) pp. 193-197.  
CODEN: 63CXA6.

COUNTRY: JAPAN

DOCUMENT TYPE: Conference

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1996:470886

LANGUAGE: English

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020730

ED Entered STN: 20011116

Last Updated on STN: 20020730

AB A review, with 18 refs. New hypothesis to prepare nonthrombogenic materials has been proposed and the supporting evidence was discussed. That was a surface which could adsorb phospholipids, would be biocompatible and a methacrylate having affinity with phospholipids was designed. 2-Methacryloyloxyethyl phosphorylcholine (MPC) was prepared, copolymd. with several monomers and their evaluation was carried out. It was found that polymers having phosphorylcholine groups, phospholipid polymers, have good affinity with phospholipids and could adsorb them on the surface. Liposomal structure was kept when the phospholipid polymers were soaked in liposomal solution. The structure was identified by XPS, comparison of the gel-liquid crystalline transition temperature of phospholipid liposome with

differential scanning calorimetry, and desorption of phospholipids. Data suggested that self-assembled biomimetic membrane was prepared on the MPC copolymers. Their unusual but interesting property was inhibition of protein adsorption even in plasma solution and blood. They did not adsorb and activate platelets. Preparation of dialysis membranes which do not require **anticoagulants** is also possible. So it was concluded that MPC copolymers are promising basic biocompatible biomaterials.

CC 63-0

ST Miscellaneous Descriptors

RN 67881-98-5Q (2-Methacryloyloxyethyl phosphorylcholine membrane biomaterial review methacryloyloxyethyl phosphorylcholine membrane biomaterial)

L82 ANSWER 56 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1994:121497 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA12008079736X

TITLE: Biocompatible water-soluble cellulose derivatives, their manufacture and uses

AUTHOR(S): Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: ASSIGNEE: NOF Corp.

PATENT INFORMATION: WO 9316117 A1 19 Aug 1993

SOURCE: (1993) PCT Int. Appl., 19 pp.

CODEN: PIXXD2.

COUNTRY: JAPAN

DOCUMENT TYPE: Patent

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1994:79736

LANGUAGE: Japanese

ENTRY DATE: Entered STN: 20011116

Last Updated on STN: 20020917

ED Entered STN: 20011116

Last Updated on STN: 20020917

AB The title cellulose derivs. are manufactured by grafting a water-soluble cellulose

substrate with 2-methacryloyloxyethylphosphorylcholine (I), and are useful for dialysis membranes for purification of blood, etc. Thus, acetylating a cellulose powder with H<sub>2</sub>SO<sub>4</sub>-catalyzed Ac<sub>2</sub>O in AcOH, and deacetylating the derivative with aqueous Na<sub>2</sub>CO<sub>3</sub> and NaOH gave a water-soluble cellulose (II) which

was then purified using membrane to remove low mol. weight products. Mixing 10 mL 0.5% aqueous solution of the purified II with 0.17 g ammonium Ce nitrate, 3

mL 0.1N HNO<sub>3</sub>, then with 1.2 g I, and stirring under Ar for 1 h at 40° gave a grafted product (III) bearing 11.3% groups derived from I, and having mol. weight (polyethylene glycol-conversion, GPC-method-based) 1.4 × 10<sup>5</sup>. Passing an aqueous solution of the III through viscose-derived hollow-fibers gave treated fibers bearing the III 8.6 µg/cm<sup>2</sup>; a dialysis module formed from the fibers showed low adhesion of platelets during purification of blood specimens.

CC 43-3

ST Miscellaneous Descriptors

**anticoagulant** cellulose methacryloyloxyethylphosphorylcholine graft manuf; dialysis membrane **anticoagulant** treatment cellulose graft

RN 142146-61-0

L82 ANSWER 57 OF 70 TOXCENTER COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1994:121480 TOXCENTER

COPYRIGHT: Copyright 2005 ACS

DOCUMENT NUMBER: CA12008078993S



TITLE: Biocompatible supported membranes useful for dialyzers, etc.

AUTHOR(S): Kamo, Jun; Nakabayashi, Norio; Ishihara, Kazuhiko

CORPORATE SOURCE: ASSIGNEE: Nakabayashi Norio

PATENT INFORMATION: JP 93177119 A2 20 Jul 1993

SOURCE: (1993) Jpn. Kokai Tokkyo Koho, 4 pp.  
CODEN: JKXXAF.

COUNTRY: JAPAN

DOCUMENT TYPE: Patent

FILE SEGMENT: CAPLUS

OTHER SOURCE: CAPLUS 1994:78993

LANGUAGE: Japanese

ENTRY DATE: Entered STN: 20011116  
Last Updated on STN: 20020917

ED Entered STN: 20011116  
Last Updated on STN: 20020917

AB The title microporous membranes are derived from hydrophobic copolymers of (aliphatic hydrocarboxy)alkyl methacrylates and 2-methacryloyloxyethyl phosphorylcholine (I). Thus, impregnation of a polyethylene hollow-fiber membrane (having pore volume 71 vol%; as support) with a 0.7% MeOH-THF 1:1 solution of a I (31.0 mol%)-Bu methacrylate copolymer (II; mol. weight 38,000) for 5 min gave a biocompatible, functionalized membrane with II content 6.3% which showed adhesion of blood platelets 8.0% in a dialysis assay, vs. 39.5% for a supported Bu methacrylate homopolymer in place of the II.

CC 38-3

ST Miscellaneous Descriptors  
membrane antithrombotic methacrylate phosphorylcholine ester copolymer;  
biocompatible membrane methacrylate phosphorylcholine ester copolymer;  
dialyzer membrane methacryloyloxyethyl phosphorylcholine copolymer  
biocompatible; prosthetic methacryloyloxyethyl phosphorylcholine  
copolymer biocompatible; ammonium phosphate inner salt polymer  
biocompatible

RN 125275-25-4 (n-Butyl methacrylate-2-methacryloyloxyethyl  
phosphorylcholine copolymer)

RN 134980-16-8; 134980-17-9

L82 ANSWER 58 OF 70 MEDLINE on STN DUPLICATE 10

ACCESSION NUMBER: 2002098928 MEDLINE

DOCUMENT NUMBER: PubMed ID: 11829441

TITLE: The vascular prosthesis without pseudointima prepared by antithrombogenic phospholipid polymer.

AUTHOR: Yoneyama Toshikazu; Sugihara Ken-ichi; Ishihara Kazuhiko; Iwasaki Yasuhiko; Nakabayashi Nobuo

CORPORATE SOURCE: The Second Department of Surgery, School of Medicine, Tokyo Medical and Dental University, Japan.

SOURCE: Biomaterials, (2002 Mar) 23 (6) 1455-9.  
Journal code: 8100316. ISSN: 0142-9612.

PUB. COUNTRY: England; United Kingdom

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 200207

ENTRY DATE: Entered STN: 20020207  
Last Updated on STN: 20020730  
Entered Medline: 20020729

ED Entered STN: 20020207  
Last Updated on STN: 20020730  
Entered Medline: 20020729

AB On the luminal surface of the common synthetic vascular prostheses, blood

**coagulation** can occur and a thrombus membrane is formed when blood flow passes through it. The thrombus membrane should be organized according to the wound healing process and it becomes a pseudointima which could serve as a blood conduit. However, the small-diameter vascular prosthesis may be quickly occluded by the initial thrombus. Therefore, no clinically applicable small-diameter prostheses have been developed to date. 2-Methacryloyloxyethyl phosphorylethanolamine (MPC) polymers resemble the structure of an outer cell membrane similar to the fluid mosaic model and demonstrate excellent antithrombogenicity. The purpose of this study is to develop a clinically applicable small-diameter prosthesis based on the new concept of the MPC polymer. We prepared vascular prostheses (2mm ID) from polymer blend composed of segmented polyurethane and the MPC polymer. The prostheses were placed in rabbit carotid arteries. The luminal surface retrieved at eight weeks after implantation was clear without thrombus and pseudointima. We now realize that the vascular prosthesis having the MPC polymer can be applied as a small-diameter prosthesis because it functions without thrombus and pseudointima formation.

CT Animals

Arteries: PA, pathology

\*Biocompatible Materials

\*Methacrylates: CH, chemistry

Microscopy, Electron, Scanning

\*Phosphorylcholine: AA, analogs & derivatives

\*Phosphorylcholine: CH, chemistry

\*Polymers: CH, chemistry

Rabbits

Time Factors

\*Tunica Intima: CH, chemistry

Tunica Intima: UL, ultrastructure

RN 107-73-3 (Phosphorylcholine); 67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)

CN 0 (Biocompatible Materials); 0 (Methacrylates); 0 (Polymers)

L82 ANSWER 59 OF 70

MEDLINE on STN

ACCESSION NUMBER: 2005301982 MEDLINE

DOCUMENT NUMBER: PubMed ID: 15948416

TITLE: Effects of surface modification of intraocular lenses on foreign body reaction.

AUTHOR: Okajima Yasuhiko; Saika Shizuya; Sawa Mitsuru

CORPORATE SOURCE: Department of Ophthalmology, Nihon University School of Medicine, Japan.

SOURCE: Nippon Ganka Gakkai zasshi, (2005 May) 109 (5) 267-73.

Journal code: 7505716. ISSN: 0029-0203.

PUB. COUNTRY: Japan

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: Japanese

FILE SEGMENT: Priority Journals

ENTRY MONTH: 200507

ENTRY DATE: Entered STN: 20050614

Last Updated on STN: 20050731

Entered Medline: 20050729

ED Entered STN: 20050614

Last Updated on STN: 20050731

Entered Medline: 20050729

AB PURPOSE: In order to improve biocompatibility, we investigated the effects of surface modification by 2-methacryloyloxyethyl phosphorylcholine (MPC) on the foreign body reaction of intraocular lens (IOLs). MATERIALS AND METHODS: Materials of the IOLs were polymethylmethacrylate, hydrophobic acryl, and MPC surface-modified hydrophobic IOLs (MPC modified acryl). In an in vitro study, cultured macrophages sampled from mouse intra-abdominal

exudate were cultured on a plate for each IOL material. The cell density and morphology of attached cells on the IOL materials were investigated. In an in vivo study, each IOL material was implanted in the peritoneal space of mice and foreign body reaction was investigated with a light microscope and a scanning electron microscope. RESULTS: In the in vitro study, the cells on the MPC modified acryl IOL material were remarkably fewer than those on the plates of the other two IOL materials. Regarding the implanted IOL materials, MPC modified acryl IOL material showed more polynuclear giant foreign body cells in the early period than the other two IOL materials. CONCLUSION: MPC surface modification can reduce the foreign body reaction of IOLs and has the potential to improve biocompatibility of IOL materials.

CT Check Tags: Male  
 Animals  
 Cell Adhesion  
 Cells, Cultured  
 \*Coated Materials, Biocompatible  
 English Abstract  
 \*Foreign-Body Reaction  
 \*Lenses, Intraocular  
 \*Macrophages, Peritoneal: IM, immunology  
 Macrophages, Peritoneal: PA, pathology  
 \*Methacrylates  
 Mice  
 Mice, Inbred C57BL  
 \*Phosphorylcholine  
 \*Phosphorylcholine: AA, analogs & derivatives  
 RN 107-73-3 (Phosphorylcholine); 67881-98-5 (2-methacryloyloxyethyl  
 phosphorylcholine)  
 CN 0 (Coated Materials, Biocompatible); 0 (Methacrylates)

L82 ANSWER 60 OF 70 MEDLINE on STN  
 ACCESSION NUMBER: 2002392253 MEDLINE  
 DOCUMENT NUMBER: PubMed ID: 12072015  
 TITLE: Beneficial effects of synthetic phospholipid polymer,  
 poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl  
 methacrylate), on stratum corneum function.  
 AUTHOR: Kanekura T; Nagata Y; Miyoshi H; Ishihara K; Nakabayashi N;  
 Kanzaki T  
 CORPORATE SOURCE: Department of Dermatology, Kagoshima University Faculty of  
 Medicine, Japan.. takurok@m2.kufm.kagoshima-u.ac.jp  
 SOURCE: Clinical and experimental dermatology, (2002 May), 27 (3)  
 230-4.  
 Journal code: 7606847. ISSN: 0307-6938.  
 PUB. COUNTRY: England: United Kingdom  
 DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
 LANGUAGE: English  
 FILE SEGMENT: Priority Journals  
 ENTRY MONTH: 200209  
 ENTRY DATE: Entered STN: 20020727  
 Last Updated on STN: 20020911  
 Entered Medline: 20020910  
 ED Entered STN: 20020727  
 Last Updated on STN: 20020911  
 Entered Medline: 20020910  
 AB The effects of a newly synthesized phospholipid polymer,  
 poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate)  
 [poly(MPC-co-BMA)], on the water barrier function and water-holding  
 capacity of the stratum corneum were examined by measuring transepidermal  
 water loss (TEWL) and electrical conductance of the skin surface. On the

backs of four NC mice, the epidermal permeability barrier was abrogated by cellophane tape stripping 30 times. The skin was then treated with 0.1% poly(MPC-co-BMA) or distilled water twice daily for the following 3 days. Poly(MPC-co-BMA) reduced TEWL significantly compared with the control after the first treatment ( $P = 0.044$ ) and this effect was observed for 3 days. In human skin, water-holding capacity was measured at 5, 10, 15, 30 min and 1, 2, and 4 h after the application of poly(MPC-co-BMA) or distilled water to both volar forearms of 21 healthy volunteers. Skin treated with poly(MPC-co-BMA) showed significantly greater ability to retain water at all time points. Poly(MPC-co-BMA) is the first synthetic material that can enhance both the water barrier function and water-holding capacity of the stratum corneum. Our results indicate that this substance may be useful clinically in the treatment of dry skin.

CT Check Tags: Female; Male

Adult

Animals

Dose-Response Relationship, Drug

\*Epidermis: DE, drug effects

Epidermis: ME, metabolism

Galvanic Skin Response: DE, drug effects

Humans

\*Methacrylates: PD, pharmacology

Mice

Middle Aged

Patch Tests

\*Phosphorylcholine: AA, analogs & derivatives

**Phosphorylcholine: IM, immunology**

\*Phosphorylcholine: PD, pharmacology

Water: ME, metabolism

\*Water Loss, Insensible: DE, drug effects

RN 107-73-3 (Phosphorylcholine); **125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate))**;  
7732-18-5 (Water)

CN 0 (Methacrylates)

L82 ANSWER 61 OF 70 MEDLINE on STN

ACCESSION NUMBER: 2001018842 MEDLINE

DOCUMENT NUMBER: PubMed ID: 10898238

TITLE: New polymeric biomaterials-phospholipid polymers with a biocompatible surface.

AUTHOR: Ishihara K

CORPORATE SOURCE: Department of Materials Science, Graduate School of Engineering, The University of Tokyo, Japan.

SOURCE: Frontiers of medical and biological engineering : international journal of the Japan Society of Medical Electronics and Biological Engineering, (2000) 10 (2) 83-95.

Journal code: 9011464. ISSN: 0921-3775.

PUB. COUNTRY: Netherlands

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 200011

ENTRY DATE: Entered STN: 20010322

Last Updated on STN: 20010322

Entered Medline: 20001109

ED Entered STN: 20010322

Last Updated on STN: 20010322

Entered Medline: 20001109

AB New biomedical polymers were designed with attention to the surface of

biological membranes, i.e. the surface was completely covered with phospholipid polar groups. The polymers with a phosphorylcholine group, 2-methacryloyloxyethyl phosphorylcholine (MPC) co-polymerized with hydrophobic alkyl group, could interact with phospholipids in plasma selectively and strongly. The adsorbed phospholipids on the polymer surface were concentrated, organized each other and then formed a self-assembled biomimetic membrane surface. The surface showed excellent resistance for both protein adsorption and blood cell adhesion, i.e. the MPC polymer showed good blood compatibility. Based on these characteristics of the MPC polymer, it was applied to improve the biocompatibility and biostability of an implantable glucose sensor. The relative output current of the sensor covered with the MPC polymer membrane was maintained as the initial level even after 14 days of subcutaneous implantation in a rat. Therefore, it is concluded that the MPC polymer membrane is an excellent material for implantable biomedical devices.

CT Adsorption

Animals

\*Biocompatible Materials: CH, chemistry

Biosensing Techniques

\*Blood Coagulation

Cell Adhesion

Glucose: AN, analysis

Humans

Membranes, Artificial

Methacrylates: CS, chemical synthesis

\*Methacrylates: CH, chemistry

Microspheres

\*Phosphorylcholine: AA, analogs &amp; derivatives

Phosphorylcholine: CS, chemical synthesis

Phosphorylcholine: CH, chemistry

Proteins: CH, chemistry

Rats

Research Support, Non-U.S. Gov't

Surface Properties

RN 107-73-3 (Phosphorylcholine); 125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate)); 50-99-7 (Glucose)

CN 0 (Biocompatible Materials); 0 (Methacrylates); 0 (Proteins)

L82 ANSWER 62 OF 70 MEDLINE on STN

ACCESSION NUMBER: 95130614 MEDLINE

DOCUMENT NUMBER: PubMed ID: 7829565

TITLE: Selective adhesion of platelets on a polyion complex composed of phospholipid polymers containing sulfonate groups and quarternary ammonium groups.

AUTHOR: Ishihara K; Inoue H; Kurita K; Nakabayashi N

CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Japan.

SOURCE: Journal of biomedical materials research, (1994 Nov) 28 (11) 1347-55.

Journal code: 0112726. ISSN: 0021-9304.

PUB. COUNTRY: United States

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 199502

ENTRY DATE: Entered STN: 19950307

Last Updated on STN: 19950307

Entered Medline: 19950223

ED Entered STN: 19950307  
Last Updated on STN: 19950307  
Entered Medline: 19950223

AB We investigated the effects of electrical charges on cell-polymer interactions of poly[2-methacryloyloxyethyl phosphorylcholine(MPC)-co-n-butyl methacrylate (BMA)] (PMB) having excellent blood compatibility, by copolymerizing anionic or cationic methacrylates with MPC and BMA. A polyion complex (PIC) composed of anionic and cationic MPC copolymers was also prepared. When the cell adhesion on these polymer surfaces from rabbit whole blood was evaluated, we observed a considerable reduction in cell adhesion on the MPC copolymers compared with that on poly(BMA), even when the MPC copolymer was electrically charged. On the other hand, many platelets selectively adhered to the PIC surface from whole blood, but the adherent platelets maintained a discoid shape. The amount of adenosine triphosphate (ATP) in platelets adherent on the PMB or the PIC from a platelet-rich plasma (PRP) was more than 75% of that in the original PRP, which indicated that the activity of these platelets remained high. However, in the platelets adherent to poly(BMA), only a small amount of ATP remained. Protein adsorption on the polymer surface from human plasma was investigated using a gold-colloid-labeled **immunoassay** against albumin gamma-globulin, and fibrinogen. Many of these proteins adsorbed on poly(BMA), whereas a small amount of protein was observed on the MPC copolymers that had an electrical charge. Albumin adsorption and suppression of gamma-globulin and fibrinogen adsorption were found on the PIC. Therefore, the introduction of electrical charges in the PMB did not have an adverse effect on cell adhesion and protein adsorption.(ABSTRACT TRUNCATED AT 250 WORDS)

CT Adenosine Triphosphate: ME, metabolism  
Adsorption  
Animals  
Blood Platelets: ME, metabolism  
Blood Platelets: PH, physiology  
\*Materials Testing  
\*Methacrylates  
Microscopy, Electron, Scanning  
\*Phosphorylcholine: AA, analogs & derivatives  
\*Platelet Adhesiveness  
\*Polymers: CS, chemical synthesis  
Polymers: CH, chemistry  
\*Proteins: PK, pharmacokinetics  
Rabbits  
Research Support, Non-U.S. Gov't

RN 107-73-3 (Phosphorylcholine); 56-65-5 (Adenosine Triphosphate);  
**67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)**

CN 0 (Methacrylates); 0 (Polymers); 0 (Proteins); 0 (butyl methacrylate)

L82 ANSWER 63 OF 70 MEDLINE on STN  
ACCESSION NUMBER: 95085403 MEDLINE  
DOCUMENT NUMBER: PubMed ID: 7993191  
TITLE: Improvement of hemocompatibility on a cellulose dialysis membrane with a novel biomedical polymer having a phospholipid polar group.  
AUTHOR: Ishihara K; Fukumoto K; Miyazaki H; Nakabayashi N  
CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Japan.  
SOURCE: Artificial organs, (1994 Aug) 18 (8) 559-64.  
Journal code: 7802778. ISSN: 0160-564X.  
PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English

FILE SEGMENT: Priority Journals  
 ENTRY MONTH: 199501  
 ENTRY DATE: Entered STN: 19950124  
 Last Updated on STN: 19980206  
 Entered Medline: 19950106

ED Entered STN: 19950124  
 Last Updated on STN: 19980206  
 Entered Medline: 19950106

AB To improve surface hemocompatibility on cellulose hollow fibers for hemodialysis, newly designed hemocompatible polymers with a phospholipid polar group, 2-methacryloyloxyethyl phosphorylcholine (MPC) polymers, were introduced on the surface through two different methods: direct grafting of MPC on the surface, or coating of a water-soluble cellulose grafted with MPC. The MPC was polymerized using cerium ion as an initiator in the cellulose hollow fibers, and the poly(MPC) chains were grafted directly on the surface. Another modification of the cellulose hollow fibers was attempted by coating them with a water-soluble graft copolymer composed of a poly(MPC) side chain and a cellulose backbone. The coating process from an aqueous solution of the graft copolymer was very convenient, and the graft copolymer on the surface was not detached even after water circulated into the hollow fibers. These cellulose hollow fibers modified with MPC polymers displayed excellent hemocompatibility such as prevention of blood cell adhesion and aggregation after contact with blood without an **anticoagulant**. The permeability of the hollow fibers did not decrease as a result of these modifications. From these results, it is clearly suggested that introduction of the MPC units was effective for improving the hemocompatibility of the hollow fibers for hemodialysis.

CT Animals  
 \*Biocompatible Materials  
 \*Blood Physiology  
 \*Cellulose  
 Humans  
 Materials Testing  
 Membranes, Artificial  
 \*Methacrylates  
 \*Phosphorylcholine: AA, analogs & derivatives  
 Platelet Adhesiveness  
 Platelet Aggregation  
 Polymers  
 Rabbits  
 \*Renal Dialysis: IS, instrumentation  
 Renal Dialysis: MT, methods  
 Research Support, Non-U.S. Gov't  
 Solubility  
 Surface Properties

RN 107-73-3 (Phosphorylcholine); 67881-98-5 (2-methacryloyloxyethyl phosphorylcholine); 9004-34-6 (Cellulose)

CN 0 (Biocompatible Materials); 0 (Methacrylates); 0 (Polymers)

L82 ANSWER 64 OF 70 MEDLINE on STN

ACCESSION NUMBER: 95197310 MEDLINE

DOCUMENT NUMBER: PubMed ID: 7890439

TITLE: Polymeric biomaterials: influence of phosphorylcholine polar groups on protein adsorption and complement activation.

AUTHOR: Yu J; Lamba N M; Courtney J M; Whateley T L; Gaylor J D; Lowe G D; Ishihara K; Nakabayashi N

CORPORATE SOURCE: Bioengineering Unit, University of Strathclyde, Glasgow, UK.

SOURCE: International journal of artificial organs, (1994 Sep) 17

(9) 499-504.

Journal code: 7802649. ISSN: 0391-3988.

PUB. COUNTRY: Italy  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 199504  
ENTRY DATE: Entered STN: 19950427  
Last Updated on STN: 19950427  
Entered Medline: 19950418

ED Entered STN: 19950427

Last Updated on STN: 19950427

Entered Medline: 19950418

AB The introduction to polymeric biomaterials of phosphorylcholine polar groups represents an approach towards the development of materials with improved blood compatibility. In this respect, two biomaterials, one a copolymer of butyl methacrylate and 2-methacryloyloxyethylphosphorylcholine (MPC), (poly(BMA-co-MPC) and the other, MPC-grafted Cuprophane, were examined with respect to their influence on protein adsorption and complement activation. Protein adsorption was studied by measurement of the adsorption of radiolabelled single proteins (albumin and fibrinogen), while complement activation was measured using **radioimmunoassay** for C3a des Arg. The investigation demonstrated that the polymers containing phosphorylcholine polar groups can achieve a marked reduction in protein adsorption and complement activation and supports the utilization of phosphorylcholine polar groups as a means of improving the compatibility of biomaterials for blood-contacting applications.

CT Adsorption: DE, drug effects

Albumins: DE, drug effects

\*Albumins: ME, metabolism

\*Biocompatible Materials: CH, chemistry

Biocompatible Materials: PD, pharmacology

Cellulose: AA, analogs &amp; derivatives

Cellulose: CH, chemistry

Complement 3a: ME, metabolism

\*Complement Activation: DE, drug effects

Fibrinogen: DE, drug effects

\*Fibrinogen: ME, metabolism

Humans

Membranes, Artificial

Methacrylates: CH, chemistry

Phosphorylcholine: AA, analogs &amp; derivatives

Phosphorylcholine: CH, chemistry

\*Phosphorylcholine: PD, pharmacology

Polymers

RN 107-73-3 (Phosphorylcholine); **67881-98-5 (2-methacryloyloxyethyl****phosphorylcholine)**; 80295-42-7 (Complement 3a); 9001-32-5

(Fibrinogen); 9004-34-6 (Cellulose); 9050-09-3 (cuprammonium cellulose)

CN 0 (Albumins); 0 (Biocompatible Materials); 0 (Methacrylates); 0

(Polymers); 0 (butyl methacrylate)

L82 ANSWER 65 OF 70 MEDLINE on STN

ACCESSION NUMBER: 94266935 MEDLINE

DOCUMENT NUMBER: PubMed ID: 8207035

TITLE: Hemocompatibility on graft copolymers composed of  
poly(2-methacryloyloxyethyl phosphorylcholine) side chain  
and poly(n-butyl methacrylate) backbone.

AUTHOR: Ishihara K; Tsuji T; Kurosaki T; Nakabayashi N

CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical  
and Dental University, Japan.



SOURCE: Journal of biomedical materials research, (1994 Feb) 28 (2)  
225-32.  
Journal code: 0112726. ISSN: 0021-9304.

PUB. COUNTRY: United States

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 199407

ENTRY DATE: Entered STN: 19940721  
Last Updated on STN: 19940721  
Entered Medline: 19940712

ED Entered STN: 19940721  
Last Updated on STN: 19940721  
Entered Medline: 19940712

AB To improve the hemocompatibility on hydrophobic biomedical materials by a simple coating technique, graft copolymers composed of a hydrophilic side chain with phospholipid polar groups and a hydrophobic backbone were synthesized. The hydrophilic chain had phospholipid polar groups, poly[2-methacryloyloxyethyl phosphorylcholine (MPC)], and the hydrophobic backbone was poly[n-butyl methacrylate (BMA)]. Because the graft copolymers obtained could dissolve in ethanol, they could be used as a coating material. When the poly(MPC-graft-BMA) was coated onto a poly(BMA) membrane, the composition of the MPC units on the surface was maintained in the bulk graft copolymer even after immersion in water. Protein adsorption on the membrane coated with the graft copolymer from human plasma detected by a gold-colloid labeled **immunoassay** was drastically decreased compared with that on glass and the original membrane. Moreover, blood cell adhesion, activation, and aggregation on the membrane after contact with human citrated whole blood were suppressed by the coating of the graft copolymer. These results clearly show that the poly(MPC-graft-BMA) is a suitable material for improving hemocompatibility on the biomedical devices because of its protein adsorption and cell adhesion resistant properties.

CT Adsorption  
Blood Proteins: PK, pharmacokinetics  
Cell Adhesion  
Erythrocytes: PH, physiology  
Humans  
\*Materials Testing  
Methacrylates: CS, chemical synthesis  
\*Methacrylates: ST, standards  
Microscopy, Electron, Scanning  
\*Phosphorylcholine: AA, analogs & derivatives  
Phosphorylcholine: CS, chemical synthesis  
Phosphorylcholine: ST, standards  
Research Support, Non-U.S. Gov't

RN 107-73-3 (Phosphorylcholine); 125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate))

CN 0 (Blood Proteins); 0 (Methacrylates)

L82 ANSWER 66 OF 70 MEDLINE on STN

ACCESSION NUMBER: 94064695 MEDLINE

DOCUMENT NUMBER: PubMed ID: 8245045

TITLE: Effects of phospholipid adsorption on nonthrombogenicity of polymer with phospholipid polar group.

AUTHOR: Ishihara K; Oshida H; Endo Y; Watanabe A; Ueda T; Nakabayashi N

CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Japan.

SOURCE: Journal of biomedical materials research, (1993 Oct) 27

(10) 1309-14.  
 Journal code: 0112726. ISSN: 0021-9304.

PUB. COUNTRY: United States  
 DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
 LANGUAGE: English  
 FILE SEGMENT: Priority Journals  
 ENTRY MONTH: 199312  
 ENTRY DATE: Entered STN: 19940201  
 Last Updated on STN: 19950206  
 Entered Medline: 19931223

ED Entered STN: 19940201  
 Last Updated on STN: 19950206  
 Entered Medline: 19931223

AB Polymers with phospholipid polar groups, 2-methacryloyloxyethyl phosphorylcholine (MPC) polymers, have excellent nonthrombogenic properties. The effects of adsorption of phospholipids on platelet adhesion and activation on the MPC copolymer with n-butyl methacrylate (BMA) were investigated with particular attention to the structure of the phospholipids adsorbed onto the polymer surface. The electrical nature of the phospholipids adsorbed on the polymer surface affected the thrombogenicity of the polymer. On the MPC polymer surface treated with an aqueous liposomal solution of acidic phospholipids, phosphatidylserine, platelet adhesion and activation occurred to a greater extent when compared to a poly(MPC-co-BMA) surface. However, on the MPC polymer surface treated with electrically neutral phosphatidylcholines, reduced thrombogenicity could be observed. Therefore, the adsorption of the phosphatidylcholines was an important factor in reducing the thrombogenicity on the polymers. Moreover, by comparison of the poly(MPC-co-BMA) to a poly(BMA), platelet adhesion and activation on these polymer surfaces depended on the adsorption state of the phosphatidylcholines. The amount of phosphatidylcholine adsorbed on the poly(MPC-co-BMA) increased with an increase in the MPC mole fraction of the copolymer. This indicates that the MPC moieties have affinity for the phosphatidylcholines. We conclude that the poly(MPC-co-BMA) can adsorb large amounts of phosphatidylcholines and that these phospholipids organize themselves. The organized adsorption layer of the phosphatidylcholines on the surface, which construct biomembrane-like surfaces, can reduce platelet adhesion and activation effectively.

CT 1,2-Dipalmitoylphosphatidylcholine  
 Adsorption  
 Animals  
 \*Biocompatible Materials  
   **Blood Coagulation**  
   Liposomes  
 \*Methacrylates  
   Microscopy, Electron, Scanning  
 \*Phospholipids  
 \*Phosphorylcholine: AA, analogs & derivatives  
 \*Platelet Activation  
 \*Platelet Adhesiveness  
   Rabbits  
   Research Support, Non-U.S. Gov't

RN 107-73-3 (Phosphorylcholine); 125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate)); 2644-64-6 (1,2-Dipalmitoylphosphatidylcholine)

CN 0 (Biocompatible Materials); 0 (Liposomes); 0 (Methacrylates); 0 (Phospholipids)

L82 ANSWER 67 OF 70 MEDLINE on STN  
 ACCESSION NUMBER: 93131994 MEDLINE

DOCUMENT NUMBER: PubMed ID: 1484061  
TITLE: Hemocompatibility of human whole blood on polymers with a phospholipid polar group and its mechanism.  
AUTHOR: Ishihara K; Oshida H; Endo Y; Ueda T; Watanabe A; Nakabayashi N  
CORPORATE SOURCE: Institute for Medical and Dental Engineering, Tokyo Medical and Dental University, Japan.  
SOURCE: Journal of biomedical materials research, (1992 Dec) 26 (12) 1543-52.  
Journal code: 0112726. ISSN: 0021-9304.  
PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 199302  
ENTRY DATE: Entered STN: 19930226  
Last Updated on STN: 19980206  
Entered Medline: 19930218

ED Entered STN: 19930226

Last Updated on STN: 19980206

Entered Medline: 19930218

AB The hemocompatibility of a polymer containing a phospholipid polar group, poly(2-methacryloyloxyethyl phosphorylcholine (MPC)-co-n-butyl methacrylate(BMA)), with human whole blood was evaluated. When human whole blood without an **anticoagulant** was contacted with polymers, the blood cell adhesion and aggregation on the polymer without the MPC moiety was extensive, and considerable fibrin deposition was observed. This phenomenon was suppressed with an increase in the polymer MPC composition. Thus, the MPC moiety in the copolymer plays an important role in the nonthrombogenic behavior of the copolymer. These results were also confirmed by the whole blood **coagulation** time on the polymer surface which was determined by Lee-White method. The adsorption of phospholipids and proteins from human plasma on poly(MPC-co-BMA) was investigated to clarify the mechanism of the nonthrombogenicity observed with the polymer. The amount of phospholipids was increased; whereas, adsorbed proteins were decreased with an increase in the MPC composition. From these results, we concluded that the phospholipids adsorbed on poly(MPC-co-BMA) play the most important role in the nonthrombogenicity of the MPC copolymer.

CT Check Tags: In Vitro

Adsorption

\*Biocompatible Materials: CH, chemistry

**Blood Coagulation**

\*Blood Physiology

Humans

\*Methacrylates: CH, chemistry

Methacrylates: ME, metabolism

Microscopy, Electron, Scanning

Microspheres

\*Phospholipids: CH, chemistry

\*Phosphorylcholine: AA, analogs & derivatives

Phosphorylcholine: CH, chemistry

Phosphorylcholine: ME, metabolism

Polyhydroxyethyl Methacrylate: AN, analysis

\*Polymers: CH, chemistry

Research Support, Non-U.S. Gov't

Thrombosis: BL, blood

RN 107-73-3 (Phosphorylcholine); 125275-25-4 (poly(2-methacryloyloxyethyl phosphorylcholine-co-n-butyl methacrylate)); 25249-16-5 (Polyhydroxyethyl Methacrylate)

CN 0 (Biocompatible Materials); 0 (Methacrylates); 0 (Phospholipids); 0 (Polymers)

L82 ANSWER 68 OF 70 BIOSIS COPYRIGHT (c) 2005 The Thomson Corporation on STN

ACCESSION NUMBER: 2005:317261 BIOSIS

DOCUMENT NUMBER: PREV200510103200

TITLE: Novel biocompatible phosphorylcholine-based self-assembled nanoparticles for drug delivery.

AUTHOR(S): Salvage, Jonathan P.; Rose, Susanna F.; Phillips, Gary J.; Hanlon, Geoffrey W.; Lloyd, Andrew W. [Reprint Author]; Ma, Iris Y.; Armes, Stephen P.; Billingham, Norman C.; Lewis, Andrew L.

CORPORATE SOURCE: Univ Brighton, Sch Pharm and Biomol Sci, Biomed Mat Res Grp, Moulsecoomb, Brighton BN2 4GJ, E Sussex, UK  
a.w.lloyd@brighton.ac.uk

SOURCE: Journal of Controlled Release, (MAY 18 2005) Vol. 104, No. 2, pp. 259-270.

CODEN: JCREEC. ISSN: 0168-3659.

DOCUMENT TYPE: Article

LANGUAGE: English

ENTRY DATE: Entered STN: 17 Aug 2005

Last Updated on STN: 17 Aug 2005

ED Entered STN: 17 Aug 2005

Last Updated on STN: 17 Aug 2005

AB Major challenges associated with nano-sized drug delivery systems include removal from systemic circulation by phagocytic cells and controlling appropriate drug release at target sites. 2-methacryloyloxyethyl phosphorylcholine (MPC) has been copolymerised in turn with two pH responsive comonomers (2-(diethylamino)ethyl methacrylate (DEA) and 2-(diisopropylamino)ethyl methacrylate (DPA), to develop novel biocompatible drug delivery vehicles. Micelles were prepared from a series of copolymers with varying block compositions and their colloidal stability and dimensions were assessed over a range of solution pH using photon correlation spectroscopy. The drug loading capacities of these micelles were evaluated using Orange OT dye as a model compound. The cytotoxicity of the micelles was assessed using an in vitro assay. The MPC-DEA diblock copolymers formed micelles at around pH 8 and longer DEA block lengths allowed higher drug loadings. However, these micelles were not stable at physiological pH. In contrast, MPC-DPA diblock copolymers formed micelles of circa 30 nm diameter at physiological pH. In vitro assays indicated that these MPC-DPA diblock copolymers had negligible cytotoxicities. Thus novel non-toxic biocompatible micelles of appropriate size and good colloidal stability with pH-modulated drug uptake and release can be readily produced using MPC-DPA diblock copolymers. (c) 2005 Elsevier B.V All rights reserved.

CC Biophysics - Bioengineering 10511

IT Major Concepts

Methods and Techniques; Biomaterials

IT Chemicals & Biochemicals

2-methacryloyloxyethyl phosphorylcholine; 2-(diethylamino)ethyl methacrylate; 2-(diisopropylamino)ethyl methacrylate

IT Methods & Equipment

photon correlation spectroscopy: laboratory techniques, spectrum analysis techniques; phosphorylcholine-based self-assembled nanoparticle: drug delivery device

IT Miscellaneous Descriptors

drug delivery; pH-responsive

RN 67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)

105-16-8 (2-(diethylamino)ethyl methacrylate)

L82 ANSWER 69 OF 70 BIOSIS COPYRIGHT (c) 2005 The Thomson Corporation on  
STN

ACCESSION NUMBER: 2005:70886 BIOSIS  
DOCUMENT NUMBER: PREV200500069303  
TITLE: A novel methodology for pre-screening anti-thrombogenicity  
of artificial organs under physiologically identical  
pulsatile environments.  
AUTHOR(S): Iwasaki, K.; Takeuchi, Y.; Saeki, W.; Umezu, M.; Ishihara,  
K.; Imachi, K.  
SOURCE: International Journal of Artificial Organs, (July 2004)  
Vol. 27, No. 7, pp. 567. print.  
Meeting Info.: 31st Annual Congress of the European Society  
for Artificial Organs (ESAO). Warsaw, Poland. September  
08-11, 2004. European Society for Artificial Organs.  
CODEN: IJAODS. ISSN: 0391-3988.  
DOCUMENT TYPE: Conference; (Meeting)  
Conference; Abstract; (Meeting Abstract)  
LANGUAGE: English  
ENTRY DATE: Entered STN: 16 Feb 2005  
Last Updated on STN: 16 Feb 2005  
ED Entered STN: 16 Feb 2005  
Last Updated on STN: 16 Feb 2005  
CC General biology - Symposia, transactions and proceedings 00520  
Biochemistry studies - Proteins, peptides and amino acids 10064  
Biophysics - Bioengineering 10511  
Cardiovascular system - Physiology and biochemistry 14504  
Blood - Blood and lymph studies 15002  
Blood - Blood cell studies 15004  
Immunology - General and methods 34502  
IT Major Concepts  
Biomaterials; Cardiovascular System (Transport and Circulation)  
IT Parts, Structures, & Systems of Organisms  
aorta: circulatory system; atrium: circulatory system; platelet: blood  
and lymphatics  
IT Chemicals & Biochemicals  
2-methacryloyloxyethyl phosphorylcholine; IgG [immunoglobulin  
G]; albumin; fibrinogen; polyurethane  
IT Methods & Equipment  
gold colloid labeling immunoassay: laboratory techniques;  
ventricular assist device: prosthetic  
RN 67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)

L82 ANSWER 70 OF 70 BIOSIS COPYRIGHT (c) 2005 The Thomson Corporation on  
STN

ACCESSION NUMBER: 2000:113997 BIOSIS  
DOCUMENT NUMBER: PREV200000113997  
TITLE: Chemical modification of silk fibroin with  
2-methacryloyloxyethyl phosphorylcholine. II.  
Graft-polymerization onto fabric through  
2-methacryloyloxyethyl isocyanate and interaction between  
fabric and platelets.  
AUTHOR(S): Furuzono, T.; Ishihara, K.; Nakabayashi, N.; Tamada, Y.  
[Reprint author]  
CORPORATE SOURCE: National Institute of Sericultural and Entomological  
Science, 2-1 Owashi, Tsukuba, Ibaraki, 305-8634, Japan  
SOURCE: Biomaterials, (Feb., 2000) Vol. 21, No. 4, pp. 327-333.  
print.  
CODEN: BIMADU. ISSN: 0142-9612.  
DOCUMENT TYPE: Article

LANGUAGE: English  
ENTRY DATE: Entered STN: 29 Mar 2000  
Last Updated on STN: 3 Jan 2002

ED Entered STN: 29 Mar 2000  
Last Updated on STN: 3 Jan 2002

AB 2-Methacryloyloxyethyl phosphorylcholine (MPC) was grafted onto silk fabric in a two-step heterogeneous system through the vinyl bonds of 2-methacryloyloxyethyl isocyanate (MOI) modified on the fabric. First, habutae silk fabric was modified with the MOI monomer in anhydrous dimethyl sulfoxide using di-n-butyltin (IV) dilaurate and hydroquinone at 35degreeC. The saturated weight gain of modified MOI monomer on the fabric was 7.3 wt% versus the original silk. Second, graft polymerization with MPC onto the MOI modified silk was conducted using 2,2'-azo bis(2-(2-imidazolin-2-yl)propane dihydrochloride) (VA-044) as an azo polymerization initiator. The weight of the grafted MPC eventually gained was about 26.0 wt%. The MOI-modified and MPC-grafted silk fabrics were analyzed by Fourier transform infrared (FT-IR) spectroscopy. To confirm the improved biocompatibility of the silk fabric, platelet adhesion was preliminarily tested measuring lactate dehydrogenase. The number of platelets adhering to polyMPC-grafted silk fabric decreased by about one tenth compared to original and MOI-modified silk after 60 min of contact with human platelet-rich plasma (1.0 X 10<sup>6</sup> platelets cm<sup>-2</sup>).

CC Blood - General and methods 15001  
Biochemistry methods - General 10050  
Biochemistry studies - General 10060  
Biophysics - General 10502  
Pathology - Therapy 12512

IT Major Concepts  
Biomaterials; Chemistry; Blood and Lymphatics (Transport and Circulation)

IT Chemicals & Biochemicals  
2-methacryloyloxyethyl isocyanate; 2-methacryloyloxyethyl phosphorylcholine; silk fibroin

IT Methods & Equipment  
platelet adhesion **assay**: analytical method

IT Miscellaneous Descriptors  
chemical modification; fabric-platelet interaction; graft-polymerization; silk fabric

ORGN Classifier  
Hominidae 86215  
Super Taxa  
Primates; Mammalia; Vertebrata; Chordata; Animalia  
Organism Name  
human  
Taxa Notes  
Animals, Chordates, Humans, Mammals, Primates, Vertebrates

RN 30674-80-7 (2-methacryloyloxyethyl isocyanate)  
67881-98-5 (2-methacryloyloxyethyl phosphorylcholine)

=> d his l81

(FILE 'HCAPLUS, TOXCENTER, WPIX, MEDLINE, BIOSIS, CANCERLIT, EMBASE,  
PASCAL, JICST-EPLUS, DRUGU, BIOTECHNO, BIOTECHDS, SCISEARCH, CONF,  
CONFSCI, DISSABS' ENTERED AT 15:25:01 ON 20 SEP 2005)

L81 6 DUP REM L80 (6 DUPLICATES REMOVED)

=> d que l81

L63 QUE ABB=ON PLU=ON ?ASSAY? OR ?IMMUNO? OR ?AGGLUT? OR E  
LISA OR RIA OR ?COAG?

L76 1325 SEA SUMIDA, K?/AU

L77 17872 SEA WADA, K?/AU

L78 14175 SEA ISHIHARA, K?/AU

L79 3656 SEA (L76 OR L77 OR L78) AND L63

L80 12 SEA L79 AND WAKO/CS,SO,PA

L81 6 DUP REM L80 (6 DUPLICATES REMOVED)

=> d ibib ed ab l81 1-6

YOU HAVE REQUESTED DATA FROM FILE 'HCAPLUS, BIOSIS, PASCAL' - CONTINUE? (Y)/N:y

L81 ANSWER 1 OF 6 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 1

ACCESSION NUMBER: 2005:284034 HCAPLUS

DOCUMENT NUMBER: 142:332458

TITLE: Latex **agglutination immunoassay**  
reagent for the determination of prostate specific  
antigen

INVENTOR(S): Sumida, Kyoichi; Fujita, Minoru; Adachi,  
Hiromichi

PATENT ASSIGNEE(S): Wako Pure Chemical Industries, Ltd., Japan

SOURCE: U.S. Pat. Appl. Publ., 9 pp.

CODEN: USXXCO

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 2005069967	A1	20050331	US 2004-867912	20040615
JP 2005106609	A2	20050421	JP 2003-340028	20030930

PRIORITY APPLN. INFO.: JP 2003-340028 A 20030930

ED Entered STN: 03 Apr 2005

AB The present invention relates to (1) a reagent for an **immunoassay** of a target substance existing in a free form and a bound form in a specimen, comprising a latex 1 which is immobilized with a monoclonal antibody 1 for the target substance, and a latex 2 which has a different mean particle size from the latex 1 and is immobilized with a monoclonal antibody 2 having a different recognition site for the target substance from the antibody 1; (2) an **immunoassay** method comprising reacting the target substance with the reagent of (1) and determining an amount of the substance based on the result of an **agglutination** reaction among the target substance, the latex 1 and the latex 2; (3) a reagent kit comprising a reagent of (1) and a reagent containing an **agglutination** accelerator for an antigen-antibody reaction.

L81 ANSWER 2 OF 6 HCAPLUS COPYRIGHT 2005 ACS on STN DUPLICATE 2

ACCESSION NUMBER: 2002:958988 HCAPLUS  
DOCUMENT NUMBER: 138:21783  
TITLE: **Agglutination**-promoting agent for antigen or antibody **immunoassay**  
INVENTOR(S): Kakuta, Kyoichi; Wada, Hiroshi; **Ishihara, Kazuhiko**  
PATENT ASSIGNEE(S): **Wako** Pure Chemical Industries, Ltd., Japan  
SOURCE: Jpn. Kokai Tokkyo Koho, 13 pp.  
CODEN: JKXXAF  
DOCUMENT TYPE: Patent  
LANGUAGE: Japanese  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2002365296	A2	20021218	JP 2001-169051	20010605
US 2004157276	A1	20040812	US 2004-626502	20040304
PRIORITY APPLN. INFO.:			JP 2001-169051	A 20010605

ED Entered STN: 18 Dec 2002

AB Provided are aggregation-promoting compds. for use in **agglutination immunoassay**. These compds are branched polymers or copolymers having basic (monomer) structure of OPO2-O-R4-N(R1R2R3) where R1-3 are independently H, OH or alkyl group; and R4 is an alkyl group or alkylene group. The **agglutination immunoassay** reagent comprises carrier- or latex-immobilized antibody or antigen. The **agglutination immunoassay** is useful for determination of antigen or antibody, e.g. C-reactive protein, rheumatic factor and prostate-specific antigen.

L81 ANSWER 3 OF 6 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2003:173633 HCAPLUS  
DOCUMENT NUMBER: 138:217873  
TITLE: Method of stabilizing substance altering in aqueous medium with heavy water  
INVENTOR(S): **Wada, Koji**; Hanada, Toshiro  
PATENT ASSIGNEE(S): **Wako** Pure Chemical Industries, Ltd., Japan  
SOURCE: PCT Int. Appl., 149 pp.  
CODEN: PIXXD2  
DOCUMENT TYPE: Patent  
LANGUAGE: Japanese  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2003018614	A1	20030306	WO 2002-JP8458	20020822
W:	AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM			
RW:	GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, SK, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG			
PRIORITY APPLN. INFO.:			JP 2001-253729	A 20010824

ED Entered STN: 07 Mar 2003



AB Disclosed are a stabilizer for a substance altering in aqueous media which comprises heavy water; a method of stabilizing a substance altering in aqueous media, characterized by causing the substance to coexist with heavy water; and a reagent for determining or detecting a substance contained in a sample comprising a substance altering in aqueous media and heavy water. The stabilizer stabilizes a substance altering in various aqueous media. Also provided is a composition containing the substance thus stabilized. The reagent is for determining or detecting a substance contained in the sample stabilized. A reagent for measuring GPT containing L-alanine 200, NADH 0.3, tris buffer (pH = 9) 25 mM, LDH 2400 IU, and sodium adipate 0.1 % in D2O was prepared

REFERENCE COUNT: 7 THERE ARE 7 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L81 ANSWER 4 OF 6 BIOSIS COPYRIGHT (c) 2005 The Thomson Corporation on STN  
DUPLICATE 4

ACCESSION NUMBER: 2002:405026 BIOSIS  
DOCUMENT NUMBER: PREV200200405026  
TITLE: Decreased marble burying behavior in female mice lacking neuromedin-B receptor (NMB-R) implies the involvement of NMB/NMB-R in 5-HT neuron function.  
AUTHOR(S): Yamada, Kazuyuki [Reprint author]; Wada, Etsuko; Yamano, Mariko; Sun, Ying-Jie; Ohara-Imazumi, Mica; Nagamatsu, Shinya; Wada, Keiji  
CORPORATE SOURCE: Division of Animal Experiment, Advanced Technology Development Center, Brain Science Institute, Riken, 2-1 Hirosawa, Wako-City, Saitama, 351-0198, Japan yamada@ncnp.go.jp  
SOURCE: Brain Research, (28 June, 2002) Vol. 942, No. 1-2, pp. 71-78. print.  
CODEN: BRREAP. ISSN: 0006-8993.  
DOCUMENT TYPE: Article  
LANGUAGE: English  
ENTRY DATE: Entered STN: 24 Jul 2002  
Last Updated on STN: 24 Jul 2002

ED Entered STN: 24 Jul 2002

Last Updated on STN: 24 Jul 2002

AB Neuromedin B (NMB) is a mammalian bombesin-like peptide distributed widely in the central nervous system. This peptide exerts its function via the NMB receptor (NMB-R). Female NMB-R-deficient mice were used to study the role that NMB/NMB-R may play in 5-HT neuron function since this relationship was suggested in previous in vitro studies. As 5-HT neurons are thought to modulate marble burying behavior, a role for NMB-R in this behavior was assessed. Relative to wild-type mice, NMB-R-deficient mice showed decreased marble burying behavior. However, depletion of 5-HT by treatment with p-chlorophenylalanine (p-CPA) increased burying behavior in NMB-R-deficient mice suggesting that increased levels of 5-HT in the brain cause a decrease in burying behavior in NMB-R-deficient mice. While HPLC analysis showed that 5-HT content in the whole brain does not differ between NMB-R-deficient and wild-type mice, an immunohistochemical analysis of brain sections showed that 5-HT expression in the dorsal raphe (DR) nucleus is elevated in NMB-R-deficient mice. Furthermore, a quantitative RT-PCR analysis revealed that 5-HT1A-receptor gene expression is downregulated in NMB-R-deficient mice at the whole brain level. These behavioral and biological results suggest that NMB/NMB-R may modulate 5-HT neuronal activity by affecting DR function.

L81 ANSWER 5 OF 6 BIOSIS COPYRIGHT (c) 2005 The Thomson Corporation on STN  
ACCESSION NUMBER: 2004:203901 BIOSIS  
DOCUMENT NUMBER: PREV200400204444

TITLE: Analysis of lipid raft domains enriched in BACE1, possible interaction domains between amyloid precursor protein and BACE1.

AUTHOR(S): Sakurai, T. [Reprint Author]; Okuno, M.; Kaneko, K.; **Wada, K.**; Nukina, N.

CORPORATE SOURCE: Lab. for Structural Neuropathology, Brain Sci. Inst., RIKEN, **Wako**, Japan

SOURCE: Society for Neuroscience Abstract Viewer and Itinerary Planner, (2003) Vol. 2003, pp. Abstract No. 730.10. <http://sfn.scholarone.com>. e-file.

Meeting Info.: 33rd Annual Meeting of the Society of Neuroscience. New Orleans, LA, USA. November 08-12, 2003. Society of Neuroscience.

DOCUMENT TYPE: Conference; (Meeting)  
Conference; Abstract; (Meeting Abstract)

LANGUAGE: English

ENTRY DATE: Entered STN: 14 Apr 2004  
Last Updated on STN: 14 Apr 2004

ED Entered STN: 14 Apr 2004

Last Updated on STN: 14 Apr 2004

AB Lipid rafts are membrane microdomains enriched in cholesterol and sphingolipids, and act as platforms for conducting a variety of cellular functions. In vitro and in vivo studies pointed out high cellular cholesterol levels as a promoting factor for the processing of amyloid precursor protein (APP) to amyloid-beta peptide. Since APP and amyloidogenic APP-cleaving enzymes, BACE1 and gamma-secretase, show cholesterol-dependent association with lipid rafts, the domain will be a critical site for amyloidogenic processing. We hypothesized that raft-dependent interaction between APP and BACE1 is regulated not only by lipid environment but also by proteins in the same domains. To search for candidate raft proteins, we isolated lipid raft fractions from mouse brains or primary cultured neurons by a modified method, and performed **immuno**-isolation with anti-APP or anti-BACE1 antibody. By mass-spectrometric and **immunological** analyses, only a limited number of proteins were detected both in APP-containing rafts and in BACE1-containing ones, suggesting that APP and BACE1 are located in different raft domains and interact transiently. Systematic identification and characterization of enriched proteins in each fraction will provide new insights into the regulatory mechanisms of interaction between APP and BACE1 in lipid rafts.

L81 ANSWER 6 OF 6 PASCAL COPYRIGHT 2005 INIST-CNRS. ALL RIGHTS RESERVED. on  
STN DUPLICATE 3

ACCESSION NUMBER: 2002-0458932 PASCAL

COPYRIGHT NOTICE: Copyright .COPYRGT. 2002 INIST-CNRS. All rights reserved.

TITLE (IN ENGLISH): Variability in cholesteryl ester transfer protein in healthy Japanese hyper-HDL-cholesterolemic subjects

AUTHOR: YOSHIDA Akihiro; KODAMA Michiteru; NOMURA Hideki; KOBAYASHI Norifumi; **SUMIDA Kyoichi**; NAITO Michitaka

CORPORATE SOURCE: Department of Clinical Laboratory, Nakatsugawa Municipal Hospital, Nakatsugawa, Japan; Department of Geriatrics, Nagoya University Graduate School of Medicine, Nagoya, Japan; Osaka Research Laboratories, **Wako** Pure Chemical Industries, Osaka, Japan; Division of Nutrition and Health, Graduate School of Life Studies, Sugiyama Jogakuen University, Nagoya, Japan

SOURCE: Internal medicine : (Tokyo. 1992), (2002), 41(5),

357-359, 13 refs.  
ISSN: 0918-2918  
DOCUMENT TYPE: Journal  
BIBLIOGRAPHIC LEVEL: Analytic  
COUNTRY: Japan  
LANGUAGE: English  
AVAILABILITY: INIST-11214, 354000108274270100

UP 20021001

AB Objective Hyper-high density lipoprotein (HDL)-cholesterolemia has been considered to be anti-atherogenic and is referred to as longevity syndrome. However, hyper-HDL-cholesterolemia induced by a cholesteryl ester transfer protein (CETP) deficiency may not be athero-protective, rather being atherogenic in nature. In a rural area in central Japan, the incidence of hyper-HDL-cholesterolemia has been found to be rather high (3.1% of healthy people). We studied healthy Japanese people in this area with hyper-HDL-cholesterolemia, particularly in relation to CETP. Methods Serum lipids were analyzed, and CETP mass was determined with an enzyme immunoassay method. Materials Blood was drawn after an overnight fast from 17 Japanese (5 males and 12 females) with serum HDL-cholesterol (C)  $\geq 100$  mg/dl. Results Serum CETP mass in hyper-HDL-cholesterolemic subjects was distributed in a wide range. Serum CETP mass was positively correlated with low-density lipoprotein (LDL)-C, apolipoprotein (Apo) B, and LDL-C/HDL-C, with statistical significance. CETP was also positively correlated with LDL-C/Apo B. Conclusion These results suggest that hyper-HDL-cholesterolemia may not be a single clinical entity, but a mixture of various pathophysiological conditions, and that the ratio of LDL-C to HDL-C and the size of LDL may be important factors in classifying these conditions.

=> file stnguide

FILE 'STNGUIDE' ENTERED AT 15:39:41 ON 20 SEP 2005

USE IS SUBJECT TO THE TERMS OF YOUR CUSTOMER AGREEMENT

COPYRIGHT (C) 2005 AMERICAN CHEMICAL SOCIETY, JAPAN SCIENCE


AND TECHNOLOGY CORPORATION, AND FACHINFORMATIONSZENTRUM KARLSRUHE

FILE CONTAINS CURRENT INFORMATION.

LAST RELOADED: Sep 16, 2005 (20050916/UP).

=>

**THIS PAGE BLANK (USPTO)**



## RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Nippon Shokubai Kagaku	1992			JP 419561 A	
Nof Corporation	1998			JP 10114800 A	HCAPLUS
Shionogi & Co Ltd				EP 141627 A	HCAPLUS
Shionogi & Co Ltd	1985			JP 6091983 A	

L82 ANSWER 19 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2000:215700 HCAPLUS

DOCUMENT NUMBER: 132:262394

TITLE: Polymer/enzyme-conjugate and polymer/enzyme/antibody-conjugate for enzyme immunoassay

INVENTOR(S): Sakaki, Shujiro; Yamada, Satoru; Shudo, Kenshiro; Nakabayashi, Nobuo; Ishihara, Kazuhiko

PATENT ASSIGNEE(S): Nippon Oil and Fats Co., Ltd., Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 15 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 2000093169	A2	20000404	JP 1998-274782	19980929

PRIORITY APPLN. INFO.: JP 1998-274782, 19980929

ED Entered STN: 04 Apr 2000

AB Polymer/enzyme-conjugate and polymer/enzyme/substance with biol. specific binding ability-conjugate are provided for the use in a highly sensitive enzyme immunoassay. This polymer/enzyme-conjugate is prepared by chemical binding an enzyme for immunol. measurement (e.g., peroxidase) with a polymer synthesized by polymerizing the monomer constituent containing a hydrophilic monomer possessing a phosphorylcholin-analog group (e.g., 2-methacryloyloxyethylphosphorylcholine (MPC) (I)) and a monomer possessing a chemical reactive group (e.g., methacrylate, 2-aminoethyl(meth)acrylate). The substance with biol. specific binding ability used for the conjugate is either antibody, biotin, avidin, or antigen. Various samples of polymer/horse radish peroxidase/biotin or IgG-conjugate prepared by this method exhibited an excellent solubility and 1.8-36 times higher sensitivity than the cases where no polymer was used to make conjugates.

IC ICM C12N011-08

ICS G01N033-532; C08F008-00; C08F220-06; C08F220-34; C08F230-02

CC 9-10 (Biochemical Methods)

Section cross-reference(s): 7

IT Immunoassay

(enzyme; polymer/enzyme-conjugate and polymer/enzyme/antibody-conjugate for enzyme immunoassay)

IT 7659-36-1, 2-Aminoethylmethacrylate 7659-38-3, 2-Aminoethylacrylate

18358-13-9, Methacrylate, reactions 67881-98-5,

2-Methacryloyloxyethylphosphorylcholine


RL: RCT (Reactant); RACT (Reactant or reagent)  
(polymer/enzyme-conjugate and polymer/enzyme/antibody-conjugate for  
enzyme **immunoassay**)

IT **67881-98-5**, 2-Methacryloyloxyethylphosphorylcholine

RL: RCT (Reactant); RACT (Reactant or reagent)  
(polymer/enzyme-conjugate and polymer/enzyme/antibody-conjugate for  
enzyme **immunoassay**)

RN 67881-98-5 HCAPLUS

CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



L82 ANSWER 20 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 2000:599226 HCAPLUS

DOCUMENT NUMBER: 133:293102

TITLE: Water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in immunoassay system

AUTHOR(S): Sakaki, Shujiro; Iwasaki, Yasuhiko; Nakabayashi, Nobuo; Ishihara, Kazuhiko

CORPORATE SOURCE: Institute of Biomaterials and Bioengineering, Tokyo Medical and Dental University, Tokyo, 101-0062, Japan

SOURCE: Polymer Journal (Tokyo) (2000), 32(8), 637-641

CODEN: POLJB8; ISSN: 0032-3896

PUBLISHER: Society of Polymer Science, Japan

DOCUMENT TYPE: Journal

LANGUAGE: English

ED Entered STN: 29 Aug 2000

AB The purpose of this study is the development of a novel synthetic blocking reagent for the ELISA method. The water-soluble amphiphilic phospholipid polymer, poly[2-methacryloyloxyethyl phosphorylcholine (MPC)-co-styrene (St)], was synthesized, and the function of the poly(MPC-co-St) as a blocking reagent was compared with conventional blocking reagents which are made of proteins such as bovine serum albumin (BSA) and casein. The poly(MPC-co-St) solution functioned at the same level as BSA solution and casein


solution for preventing non-specific antibody adsorption ( $p > 0.01$ ). When the 1.0% BSA solution and 1.0% casein solution were used as a blocking reagent, the remaining activity of the immobilized antibody decreased about 50% after 20 days. On the other hand, in 0.01% and 0.1% poly(MPC-co-St) solns., the activity remained 76% and 91% of the initial value, resp. The effects of poly(MPC-co-St) on the stabilization of the immobilized antibody depended on its concentration. These results indicated that the poly(MPC-co-St) had the ability to inhibit denaturation of protein, i.e., proteins in the ELISA system kept their native structure. We concluded that the water-soluble amphiphilic poly-(MPC-co-St) is an effective synthetic blocking reagent in the ELISA method.

CC 9-10 (Biochemical Methods)


IT **Immunoassay**

(enzyme-linked immunosorbent assay; water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in immunoassay system)


**134483-35-5P**  
 RL: ARU (Analytical role, unclassified); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation)  
 (water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in **immunoassay** system)  
 IT 100-42-5, reactions **67881-98-5**, 2-Methacryloyloxyethyl phosphorylcholine  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in **immunoassay** system)  
 IT **134483-35-5P**  
 RL: ARU (Analytical role, unclassified); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation)  
 (water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in **immunoassay** system)  
 RN 134483-35-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide, polymer with ethenylbenzene (9CI) (CA INDEX NAME)  
 CM 1  
 CRN 67881-98-5  
 CMF C11 H22 N O6 P



CM 2  
 CRN 100-42-5  
 CMF C8 H8



IT **67881-98-5**, 2-Methacryloyloxyethyl phosphorylcholine  
 RL: RCT (Reactant); RACT (Reactant or reagent)  
 (water-soluble 2-methacryloyloxyethyl phosphorylcholine copolymer as a novel synthetic blocking reagent in **immunoassay** system)  
 RN 67881-98-5 HCAPLUS  
 CN 3,5,8-Trioxa-4-phosphaundec-10-en-1-aminium, 4-hydroxy-N,N,N,10-tetramethyl-9-oxo-, inner salt, 4-oxide (9CI) (CA INDEX NAME)



RETABLE

Referenced Author (RAU)	Year (RPY)	VOL (RVL)	PG (RPG)	Referenced Work (RWK)	Referenced File
Anderton, B	1980	2	122	Immunol Today	
Donovan, J	1975	250	1966	J Biol Chem	HCAPLUS
Farr, A	1981	47	129	J Immunol Method	HCAPLUS
Imagawa, M	1982	4	41	J Appl Biochem	HCAPLUS
Ishihara, K	1991	25	1397	J Biomed Mater Res	HCAPLUS
Ishihara, K	1998	39	323	J Biomed Mater Res	HCAPLUS
Ishihara, K	1994	32	859	J Polym Sci Part A	HCAPLUS
Ishihara, K	1990	22	355	Polym J	HCAPLUS
Ishihara, K	1999	31	1231	Polym J	HCAPLUS
Ishikawa, E	1983	4	209	J Immunoassay	HCAPLUS
Isikawa, E	1983	18	219	Immunoenzym Tech	
Kojima, M	1991	12	121	Biomaterials	HCAPLUS
Lu, D	1991	3	127	J Biomater Sci Polym	HCAPLUS
Orci, L	1985	28	528	Diabetologia	HCAPLUS
Osborn, M	1978	77	R27	J Cell Biol	MEDLINE
Sakaki, S	1999	47	523	J Biomed Mater Res	HCAPLUS
Shaw, D	1980			Introduction to coll	
Wisdom, G	1976	22	1243	Enzyme-Immunoassay	HCAPLUS

L82 ANSWER 21 OF 70 HCAPLUS COPYRIGHT 2005 ACS on STN

ACCESSION NUMBER: 1998:274659 HCAPLUS

DOCUMENT NUMBER: 129:26995

TITLE: Polymeric solid support-immobilized antigen or antibody and its use

INVENTOR(S): Sakaki, Shujiro; Shudo, Kenjiro; Yamada, Akira; Matsuyama, Kazuo; Nakabayashi, Nobuo; Ishihara, Kazuhiko

PATENT ASSIGNEE(S): Nippon Oil and Fats Co., Ltd., Japan; Nakabayashi, Norio; Ishihara, Kazuhiko; Foundation for Scientific Technology Promotion

SOURCE: Jpn. Kokai Tokkyo Koho, 9 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 10114800	A2	19980506	JP 1996-271126	19961014
PRIORITY APPLN. INFO.:			JP 1996-271126	19961014

ED Entered STN: 13 May 1998

AB The disclosed antigens or antibodies are coupled to polymeric solid support through phosphorylcholine groups and used for immunoassay. The phosphorylcholine-containing polymer is e.g. polymer comprising 2-methacryloyloxyethyl-2'-(trimethylammonio)ethylphosphate (MCP). Copolymers of MCP and Bu methacrylate, methylmethacrylate or 2-hydroxyethyl methylmethacrylate were prepared, coated with anti-mouse antibody for immunoassay.

IC ICM C07K017-08  
ICS G01N033-543; C07K016-00

CC 15-2 (Immunochemistry)  
Section cross-reference(s): 9

IT **Immunoassay**  
(polymeric solid support-immobilized antigen or antibody and its use)

IT 67881-98-5DP, polymers and copolymers 67882-00-2P  
125275-25-4P 134483-35-5P 148569-41-9P